


3.0 - Basic Structural Components:

V. Wiring Harnesses

The wiring harnesses used by Motion Concept's are also manufactured by Perpetual Motion Enterprises (PME). A wide variety of wiring harnesses are needed in order to bring a TRx power positioning system to life. The common wiring harnesses used on our seating system would include *power cables*, *actuator harnesses*, and *limit harnesses*. With some of our more complex power positioning systems there is often a need for additional harnesses and cables to operate the systems safely and effectively.

Due to the often complex wiring configurations, Motion Concepts has created detailed wiring schematic diagrams & instructions for each of our seating systems to assist Technicians. These wiring schematics have been assembled in three separate Wiring Diagram Manuals* for each of the major powerbase manufacturers (Invacare, Quantum and Quickie). Please refer to **Figure 4.0** (p.31) for a sample of a wiring diagram from one of our multi-function relay boxes.

*please contact Motion Concepts for a copy of the individual Wiring Manuals.



Power Cable
eg. C78 (Invacare) power cable



Actuator Harness
eg. C45B-18 (4 motor) harness



Limit Harness
eg. C44B dual limit harness

Power Cables- draw power from the base (via the batteries or through the controller) in order to operate the power positioning system. The power cable is the first link in the series of wiring harnesses and typically connects into the actuator harness. Power cables vary in style and configuration depending on the power base manufacturer and may also vary for different powerbase models from the same manufacturer.

Actuator Harnesses- link the power cables to the control/relay boxes and distribute the power to the individual motors on the TRx Power Positioning Systems. There are three basic types of actuator harnesses used with our seating systems: *Single Motor*, *Dual Motor* and *Four Motor*.

Limit Harnesses- are controlled via the relay box and may be connected to the relay box via an input port, or they may be hard-wired directly to the relay box. Motion Concepts uses two styles of standard "plug-in" limit harnesses: Dual Limit and Four Limit. (**Note:** *These harnesses are connected to the actual limit switches (i.e.; mercury switches/ microswitches/ roller switches)* Please refer to **Section 4.0-part II. Programming & Adjustment** (p.43) for detailed information on limit switches.


NOTE: Detailed information on each individual relay/control box can be found in the Motion Concepts Electronics Manual (provided separately)

NOTES

3.0 - Basic Structural Components:

V. Wiring Harnesses

Figure 4.0 - Sample Wiring Diagram**



****Note:** Figure 4.0 above illustrates a valid wiring diagram for our Helix 2.0 Control Box at the time of printing the technical training manual. Wiring Diagrams are created & controlled by Motion Concepts for each of the major powerbase manufacturers (Invacare, Pride and Quickie). Up-to-date wiring diagrams are provided for reference with each specific seating system. Please contact our Technical Service Dept. regarding any wiring issues or concerns.