



# 2006 NEW PRODUCT feature sheet

## Uplift System

*Newly Revised and Improved!*

Newly revised to provide added vertical tilt, the Uplift now provides up to **80° lift** — up from 45° from the original Uplift. This increased angle assists with proper positioning, as well as self and stand assist transfers.

Other features include:

- ▶ Power module lifts back of seat pan forward 32°, 37°, 45° or **80°!**
- ▶ As seat lifts, the back and legrest are mechanically linked to open the back angle while the legrest extends to ground level
- ▶ Available with synchronized Power Recline module
- ▶ Optional Center Mount Power Articulating & Elevating Foot Platform\*
- ▶ Elevating module available\*\*

\* Available only with elevate module

\*\* Elevating module recommended for longer leg lengths

Now  
with **80° lift**



Pictured at full 80° angle



Pictured without  
elevate option



New High Back

Pictured with standard  
manual uplift, power  
elevating seat option  
and standard center  
mount platform



*Whether you're coming or going, these new systems get you there simply, safely and comfortably.*

## Horizontal Power Positioning System

### Features

- ▶ Adjustable center portion with telescoping end supports at each side to accommodate various shapes and sizes
- ▶ Standard growable, telescoping end supports are fixed or power adjustable to provide 25° movement up or down and aid with visibility, transfers and repositioning
- ▶ Sturdy, weight sustainable, manual flip down transfer side rails
- ▶ Ability to drive in lying down position on either stomach or back
- ▶ Center portion tilts 27°, 32° or 45°
- ▶ All 3 sections come in standard sizes and are fully width and depth adjustable:
  - Center body adjusts from 14" x 14" – 20" x 24" in 3" increments
  - End supports adjust from 14" x 13" – 14" x 24" in 3" increments
  - Telescoping transfer side rails adjust 8" x 18" – 8" x 24" in 3" increments
- ▶ Joystick can be positioned at desired spot on any of the three sections
- ▶ Single Function ¼ drive speed standard
- ▶ Standard 20.25" STF
- ▶ Cushion: 2" thick (1" visco, 1" HR)
- ▶ Side storage pouch for extra storage

### Options:

- ▶ Power adjustable end supports with 25° adjustability\*
- ▶ Power flip down transfer side rails
- ▶ 7" power elevation module
- ▶ Custom STF of 18.25" – 30" without elevating module
- ▶ Custom STF of 18.75" – 24" with elevating module
- ▶ Fixed foot stop
- ▶ Vent applications via RFQ
- ▶ Custom sized end supports via RFQ

Fits on a wide variety of popular bases from

