

Interface system:



Colin, 38, is a Project Manager for Public Works, a government office. He has been in a wheelchair since he was injured in a motorcycle accident 20 years ago.

Interface system *standard features & options*

STANDARD FEATURES

- ▶ Weight capacities up to 250 lbs
- ▶ 55° tilt
- ▶ Seat sizes 15" - 22"
- ▶ Sturdy CG Tilt fits most manufacturers' base/rehab seat kits
- ▶ Low seat-to-floor heights
- ▶ Easily retrofits to existing rehab seat
- ▶ Hassle-free and installed by our trained technicians at our facility
- ▶ Can be stocked for minimal turnaround time
- ▶ Lifetime structural warranty

POSITIONING OPTIONS

- ▶ CG Tilt Interface System uses some base manufacturers' rehab seats

ELECTRONICS

All electronics allow users to control their system through either a separate switch, toggle and/or their wheelchair drive controls

- ▶ Standard electronics included
- ▶ Enhanced electronics available
- ▶ Electronics options (ie. Lazarus Jump Starter, M16 Splitter, M06 Attendant Power Switch, Stealth Egg Switch)

With all the changes to funding codes these days, securing appropriate funding for tilt can be difficult. We created our CG Tilt Interface Kit especially for challenging funding climates. It utilizes the power base manufacturers' base/rehab seat, and adds our sturdy, low STF, CG Tilt mechanism.

Features include:


- ▶ Works in conjunction with manufacturers' seat frames and front rigging, including Quantum's Synergy Seat and Specialty Seat, Quickie's Rehab seat and Invacare's ASBA Seat
- ▶ Seat sizes 15" - 22"
- ▶ STF as low as 17.75"
- ▶ Sleek, attractive design

Available on:



3G Storm
TDX3

Freestyle F11/M11



1121
1122
Q600
Q610
Q6000



Options

 <p>MaTRx-V Cushion</p>	 <p>MaTRx-BP Back</p>	 <p>Legrests</p>	 <p>Headrests</p>
--	--	--	--

Available Product Feature Reference

 <p>TILT</p>	 <p>ELECTRONICS</p>
--	--