

product pictionary:

specialized positioning options



power headrest

- ▶ Headrest can retract backward up to 4" and lower out of the way up to 3"
- ▶ Fits any standard Motion Concepts Back
- ▶ Can be controlled with separate switch or as an additional switch on multi-function system
- ▶ Available via RFQ



power sliding back

- ▶ Back moves up and down to maximum 5" independent, and in synchronization with Recline
- ▶ Used instead of ESR
- ▶ Can be used with most standard Motion Concepts backs greater than 20" in height
- ▶ Requires enhanced electronics for independent operation



power swing away chin level mounting system

- ▶ Can be used for mounting of drive controls and/or power positioning controls
- ▶ Uses a switch boom, swings away to free client's face from obstruction
- ▶ Can be used on all standard Motion Concepts backs
- ▶ Left and right hand swing-in available
- ▶ Allows client to adjust height of chin control in 1.5" range



power swing away backpack kit

- ▶ Backpack swings around 270° from home position at rear to beside armrest in order to provide users with access to stored items
- ▶ Available in right and left swing down versions
- ▶ Available on Recline systems, Tilt-only via RFQ
- ▶ 50 lb pack storage capacity