

SuperLowSP system:



Luke, an engineer, was injured in a mountain bike accident at age 23.

Now there's no need to compromise! You can get the modularity, function and durability you need in a low STF Tilt system that also fits into the newest funding codes and allowables with our sturdy SuperLow SP Tilt system.

Features include:

- ▶ Single pivot 45° tilt
- ▶ STF as low as 17.75" with optional Recline and Anterior Assist
- ▶ STF as low as 15.75" with **tilt-only** SuperLow systems
- ▶ Available with 168° Recline
- ▶ ESR and/or Power Sliding Back available
- ▶ 300 lb weight capacity, depending on base

SuperLowSP system *standard features & options*

STANDARD FEATURES

- ▶ All systems come with seat assembly, seat pan, rigid back pan, height adjustable armrest, standard arm pads and Fixed 70° Center Mount Front Rigging
- ▶ Weight capacities up to 300 lbs
- ▶ Seat Sizes 16" - 21" W x 15" - 20" D, depending on base
- ▶ 3" of adjustability
- ▶ Lifetime structural warranty

POSITIONING OPTIONS

- ▶ SuperLow Tilt Only System (45° Single Pivot) - STF 15.75" - 17.75", depending on base model
- ▶ SP Tilt System (45°)
- ▶ Optional Recline module (not available with SuperLow Tilt)
- ▶ Optional Anterior Assist (not available with SuperLow Tilt)
- ▶ Optional Precline module available with Recline

SEATING OPTIONS

- ▶ Single Post Back (not available with Recline)
- ▶ MaTRx-V cushions
- ▶ MaTRx PB Backs
- ▶ Contour Backs
- ▶ Onyx Package
- ▶ Headrests

OPTIONAL ACCESSORIES

- ▶ Fixed and Swing Away Laterals
- ▶ Contour Arm Troughs and Hand Pads
- ▶ Armrest Storage Compartment - Glove Box
- ▶ Arm Pad Storage - Cup Holder

FRONT RIGGING OPTIONS

All power front rigging may increase STF by approximately 1"

- ▶ Center Mount Power Elevating & Articulating Foot Platform
- ▶ Pivot Plus Power Foot platform
- ▶ Pivot Plus Elevating Legrests - Power & Manual
- ▶ Fixed Center Mount Foot Platform
- ▶ Heavy Duty 70° Lift & Swing Footrests
- ▶ Split Power Elevating Articulating Foot Platform

ELECTRONICS

All electronics allow users to control their system through either a separate switch, toggle and/or their wheelchair drive controls


- ▶ Standard electronics included for Tilt and Tilt/Recline systems
- ▶ Enhanced electronics available
- ▶ Electronics options (ie. Lazarus Jump Starter, M16 Splitter, M06 Attendant Power Switch, Stealth Egg Switch)

OTHER OPTIONS

- ▶ O₂ Holders - Articulating and Non-Articulating
- ▶ Vent Trays available through RFQ

Available on:

3G Storm
TDX3,4,5
M91
P220
P222
P525
1121
1122
1400
1402
1420
Blast
Q600
Q6000



Options



Available Product Feature Reference

