

# F-SERIES TILT SYSTEM:

Shown on an Invacare TDX-SP with optional transfer handles  
 Seat-to-floor height starting at 16½"



**NEW** Cantilever (Independently Height & Angle Adjustable) Flip Back Armrests with Gel Armpads

## AVAILABLE ON:



Storm 3G  
 TDX 3, 4, 5  
 TDX SP  
 M71  
 M91



Rhythm  
 Groove  
 Freestyle F11/M11  
 P200/P220/P222



1103  
 1107  
 1121  
 1122  
 1400 Series  
 1113 ATS  
 Dynamo ATS  
 600 Series  
 6000 Series

Call for additional bases



Tilt Position



Optional Transfer Handles



Armrest Spacer

## PRODUCT REFERENCE

See pages 4 - 5 for icon references



**A Fabulous, Fully capable and Functional tilt system made especially for today's unique Funding situations with all the Fundamental Features and Flexibility you have Found with Motion Concepts products.**

### Highlights:

- ▶ 55° of CG Tilt
- ▶ Seat widths of 16" - 22" in 2" increments
- ▶ Infinitely adjustable seat depths from
- ▶ CG shift up to 7"

- ▶ Seat-to-floor heights from 16½" depending on the base
- ▶ Mounting brackets for base to base interchangeability
- ▶ Standard with fixed center mount foot platform, yet compatible with all front rigging including power
- ▶ Slotted seat pan and steel back pan
- ▶ Cantilever height and angle adjustable flip back armrests
- ▶ Slotted telescoping legrest receivers offer infinite seat depth adjustment within a 1" - 5" range
- ▶ Standard electronics with single push button, dual push button or low amp toggle control
- ▶ Elevating seat module available