

ARMPAD OPTIONS:

1. STANDARD ARM PADS 2. WATERFALL ARM PADS



- STANDARD ARM PADS**
- >Standard on all Motion Concepts systems
 - >Sold as a pair
 - >Available in 2 models:
 - Full:** 2¼" wide x 14" long
 - Desk:** 2¼" wide x 10" long

- WATERFALL ARM PADS**
- >Pad is shaped to 'fall' over armrest bar
 - >May be mounted inward or outward
 - >Sold as a pair
 - >Available in 2 models:
 - Full:** 2¼" wide x 14" long
 - Desk:** 2¼" wide x 10" long

GEL ARM PADS



- >Fabric covered polymer gel arm pads
- >Provides extra protection for highly sensitive arms
- >Available in 3 models:
- Narrow:** 2" wide x 12" long
- Wide:** 3½" wide x 12" long
- Long:** 3½" wide x 14" long

MULTIPOSITION FLAT ARM PAD WITH BUMPER GUARDS



- >Soft, gel-like feel provides extra comfort and protection for arms
- >Highly durable skinned polyurethane with bumper guards protect against cracking and peeling
- >Triple track extrusion allows adjustability of mounting position and attachment location for other accessories
- >Available in 2 models:
- Full:** 3" wide x 14" long
- Desk:** 3" wide x 10" long

ERGONOMIC ARM THROUGH FLAT HAND PAD



- >Molded, skinned polyurethane fits the natural arm curvature
- >Unique cut-out at elbow provides users with flexibility of various arm positions
- >Triple track extrusion allows a variety of mounting positions
- >Both Arm Through and arm pad are available in left & right
- Size:** 4½" wide x 13" long
- >Hand pad is an optional add-on; can be interfaced with other manufacturers' hand piece accessories
- Pap Size:** 4½" wide x 5" long