

INTRODUCING!

ROVIX₃

ontrax

ENHANCED DRIVING MODULE

HIGH SPEED MOTORS



ontrax

ENHANCED DRIVING MODULE

- Utilizes gyroscopic technology to deliver exceptional straight line driving performance
- Better straight line tracking requires fewer directional commands and reduces consumer fatigue
- Especially beneficial for individuals that utilize a non-proportional input device for drive control, such as a head array or sip-n-puff

HIGH SPEED MOTORS

- Increases the top speed of the ROVI X3 to 6.2mph
- Faster acceleration
- Top speed is accomplished with a minimal reduction in torque for balanced performance in a variety of settings
- OnTraxx comes standard with the 6.2 mph package. It is optional with the 5.8 mph package



Call Motion Concepts to order: 1.888.433.6818 | motionconcepts.com

