



For questions/inquiries, contact your local Invacare Territory Business Manager.

To place your orders, request a quote or speak with technical services, please contact Motion Concepts at 1.866.748.7943









## MODULAR MAXX System Vitals:

- Industry's most modular system, featuring up to 50° of CG Tilt, 170° of Recline with Extended Shear Reduction (ESR) and up to 12" of nearly silent Power Seat Elevation
- Utelizes ROVI Tracking Module that decreases the number of user driver control commands to maintain the wheelchair on its intended course when going over uneven or sloping terrain
- Available with the full range of Motion Concepts' accessories and Invacare® Matrx® Seating products, mounted via our all new Ultra Rail and Ultra Track slotted interfaces
- 300 lb. weight capacity
- Seat-to-Floor Heights as low as 17" with 12" of Power Seat Elevation
- 11" 22" of seat width and 11" 23" of seat depth adjustment

8A, Off Board

ADP CODE	WP0001287	<b>Battery Charger</b>
Suspension System	Active Ride Control (ARC)	Batteries Weight
Base Width	23.25"	Range Ø
Base Length	36.5"	Weight Capacity
Seat-to-floor Height	as low as 17"	Base Weight
Drive Wheel Size ®	14"	Ground Clearanc
Caster Wheel Size (front/rear) <sup>(n)</sup>	6"	ROVI X3 5.0 Maxi
Turning Radius	20"	ROVI X3 6.0 Maxi
Drivetrain	4-pole EAD Motors	(1) All X3 wheels come w

2 x M34

3 x M50 option for Vents

Batteries Weight (each)	42 lbs/19kg
Range Ø	17.7 miles/28.5km
Weight Capacity	300 lbs/136 kg
Base Weight	157 lbs/71kg
Ground Clearance	3"
ROVI X3 5.0 Maximum Speed	5 mph
ROVI X3 6.0 Maximum Speed	6.2 mph

. (1) All X3 wheels come with solid inserts

## **Fender Colors**

**Batteries** 

Base



## **Modular Maxx Power Positioning System**

POWER POSITIONING SYSTEMS	SYSTEM DESCRIPTION	ADP CODE
Maxx Adjustable Rehab Seat	no power options	WEPN
Modular Maxx CG Tilt Only System	up to 50°	WP0001295
Modular Maxx CG Tilt and Recline System with ESR	up to 50° Tilt, 170° Recline with (ESR) Extended Shear Reduction	WP0001297
Modular Maxx CG Tilt and Recline System with PSB	up to 50° Tilt, 170° Recline with (PSB) Power Sliding Back	WP0001297
Modular Maxx Recline System with ESR	170° Recline with (ESR) Extended Shear Reduction	WP0001296
Modular Maxx Recline System with PSB	170° Recline with (PSB) Power Sliding Back	WP0001296
Modular Maxx Elevating Lift Seat Only System	12" of Power Seat Elevation	WP0001295
600 Series Multi-Function Control Box		SECZG5005
Seat Width	11" - 22" width adjustment ranges 11"-15"(1), 16" - 20", 19" - 22"	
Seat Depth	1" - 23" depth adjustment ranges 11"-15"(1), 16" - 21", 19" - 23"	
Additional Modules	Recline Body Fit; Power Anterior Tilt up to 45°; 5° & 10° Fixed Anterior	
	Assist; 5° Pre-Tilt; 8° or 30° Precline; Latera Tilt up to 15° and Single	
	Post Angle Adjustable Fold Down Back	

(1) Pediatric Seat sizes offer up to 50° of CG Tilt, 165° of Recline and Seat-to-floor is increased by 1"







<sup>(2)</sup> Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems