

### Size really does matter ...

Although the jury's still out on the importance of size in certain respects, when it comes to the fit, comfort and stability of power positioning and seating for larger clients, it's clear that size is key.

And that's where Motion Concepts comes in. Motion Concepts presents complete bariatric solutions for all your clients needs. From our newest UltraLow Bariatric Power Positioning systems to our MaTRx Seating Series' Heavy Duty Cushions and Backs, our thoughtfully designed offerings are designed to fit all shapes and sizes.

#### UltraLow Bariatric features include:

- ▶ Standard 500 lb weight capacity (depending on base) and seat sizes from 19" - 24" (width and depth); larger sizes available thru RFQ
- ▶ CG shift up to 7" for maximum weight distribution
- ▶ Innovative slotted channels allow whole system to be easily adjustable on base (3" - 5" fore and aft)\*
- ▶ Independent, telescoping front rigging and sliding seat pan for easy seat depth adjustment
- ▶ Up to 45° tilt, 155° power recline and 30° precline
- ▶ Endomorph back canes available with tilt only (as shown)
- ▶ Seat-to-floor heights from 19.25"

#### MaTRx Seating Series includes:

- ▶ **MaTRx V Heavy Duty Cushion** accommodates sizes up to 30" W and up to 26" D; 650 lb weight capacity
- ▶ **MaTRx PB Heavy Duty Back** accommodates width sizes 20", 22", 24" in 16" or 20" heights; 500 lb weight capacity with 5"/7" Contour Depth

So go ahead and show your larger clients that when it comes to their power positioning and seating needs, you understand the importance of size. Provide them with a customized, perfect fit that combines comfort, stability and durability with an UltraLow Bariatric power positioning system and heavy duty cushions and backs from MaTRx Seating Series.

### ✓ Bariatric seating and positioning


The Bariatric UltraLow/MaTRx-V Heavy Duty Cushion/PB Heavy Duty Back are available for most power positioning systems on bases from:



UltraLow Bariatric



\*Slotted Channels



Concepts  
**Motion**

US: 700 Ensminger Rd. Suite 112  
Tonawanda, NY 14150  
Toll Free Tel: 1.888.433.6818  
Toll Free Fax: 1.888.433.6834

Canada: 84 Citation Dr.,  
Concord, ON L4K 3C1  
Toll Free Tel: 1.866.748.7943  
Tel: 905.695.0134  
Fax: 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)