

July 19, 2004

**Motion Concepts**

**ON A ROLL**  
A weekly e-bulletin from Motion Concepts

## Free as a bird ...

As a therapist or mobility provider, you're always searching for superior power positioning solutions for your clients. However, given the realities of today's precarious funding climate, you may not always feel free to utilize the best options available. And that's where Motion Concepts can help.

Motion Concepts provides superior power positioning solutions that fit seamlessly on the mobility industry's most popular bases. And Quickie's latest offering — the F-11 Freestyle — is no exception.

The F-11 Freestyle has garnered positive feedback for its smaller, lighter, more maneuverable mid-wheel-drive package with independent all-wheel suspension, 3-bar linkage system and "Overland Suspension". And while the Freestyle's standard seating does meet some clients' basic requirements, it doesn't include the versatile power positioning options that enhance the quality of their lives.

The all new CG Tilt Low-Pro from Motion Concepts picks up where most standard seating systems leave off. It provides an extraordinary rehab CG Tilt seating system alternative that features a low seat-to-floor height power positioning system based on a much smaller subframe than our traditional TRx system. This means the Low-Pro can take your clients as much as 2" lower than other CG tilt seating systems, provides 55° of tilt and 174° of recline — all at an amazing 17 $\frac{7}{8}$ " STF. Best of all, it's bundled specifically to fit into the new HCPCS codes for Power Positioning — and it fits seamlessly on the Quickie Freestyle!

Our exclusive EASE feature is also included with the Low-Pro tilt recline system. Plus, the CG tilt Low-Pro is available as a tilt-only version, which comes with two amazing, no-charge options: a pair of handy, Armrest Storage Compartments and an eye-catching, Single Post Back. Plus, an exclusive, low STF Power Elevating Seat Module is also available with seat elevation heights of 7" for a 225 lb. client or 6" for a 280 lb. client.

So, when your clients need a superior, low STF CG tilt system that's compatible with their Quickie Freestyle, feel free to set your sights higher with a new Low-Pro from Motion Concepts.

**U.S.**  
700 Ensminger Rd.,  
Suite 112  
Tonawanda, NY 14150  
TEL 1.888.433.6818  
FAX 1.888.433.6834

**Canada**  
84 Citation Dr.,  
Concord, ON  
L4K 3C1  
TEL 905.695.0134  
FAX 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)

