


## We are all about support . . .



**The all NEW Cane Mounted Height, Width and Angle Adjustable Flip Back Armrests offer a cleaner, more streamlined look as well as unsurpassed functionality and support!**

- Simple, independent height and angle adjustment to aid in individualized support and comfort positioning
- Width adjustment from a default of 0" to 0.5" or 1" per side; maximum 2" wider than the seat pan
- Available in two height profile configurations: Straight or Curved
  - Straight accommodates arm heights of 9.5" to 13"
  - Curved accommodates arm heights of 12.5" to 16"
- Flips back easily for unobstructed lateral transfers
- Designed to interface seamlessly with all of our available arm pads:
  - Standard | Modular | Waterfall | Flat | Gel | Ergonomic Arm Trough

**Straight Arm Configuration**



**Curved Arm Configuration**

