

### Put your best foot forward

Tired of trying to fit the proverbial square peg into a round hole? At Motion Concepts, we understand how awkward it can be to try to fit clients into standard leg configurations. That's why we created the Split Power Elevating Articulating Foot Platform.

Our Split Power Elevating Articulating Foot Platform is designed to create customized leg positioning along with comfort, security and stability. It also provides a natural, comfortable low pivot that helps facilitate easier transfers. So clients get all the benefits of independent elevating legrests with a smaller overall footprint comparable to that of a center mount!

But that's not all! This unique front rigging provides unparalleled flexibility and can be affixed anywhere across the front of the seat pan! This is ideal for larger clients because they get the wider seat frame they need without having to position their legs uncomfortably either fixed far apart or too close together on standard front rigging!

All in a stylish, sturdy package that's as easy to position as it is to use. How's that for clients being able to put their best foot forward?.

#### Features include:

- ▶ Can be installed at multiple positions across the front of the chair
- ▶ Ideal for bariatric applications
- ▶ Small overall footprint comparable to center mount
- ▶ Independent angle adjustable legrests provide 65° or 75° range of motion, from 90° or 70° resting to 165° extended
- ▶ Extended articulation of 3" or 6", depending on extension length
- ▶ Footplates flip back for transfers
- ▶ Footplate plantar flexes as it elevates to minimize stretching hamstrings and calf muscles

It's never been easier for your clients to 'get one foot in the door'. Now you can provide them with the comfort and security of customized leg positioning and help them put their best foot forward. Check out the Split Power Elevating Articulating Foot Platform today.

### ✓ Split Power Elevating Articulating Foot Platform

The Split Power Elevating Articulating Foot Platform is available for most power positioning systems on bases from:



**Split Power Elevating Articulating Foot Platform**

Concepts  
**Motion**

US: 700 Ensminger Rd. Suite 112  
Tonawanda, NY 14150  
Toll Free Tel: 1.888.433.6818  
Toll Free Fax: 1.888.433.6834

Canada: 84 Citation Dr.,  
Concord, ON L4K 3C1  
Toll Free Tel: 1.866.748.7943  
Tel: 905.695.0134  
Fax: 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)