

... ANOTHER SOLUTION

SINGLE AND LOVING IT!

Two single post options with multiple, functional benefits - It doesn't matter which one you choose, our single post backs are designed to deliver the function you have come to expect from Motion Concepts in a clean, appealing, non-rehab style. Both options are available with tilt or tilt/elevate systems and carry a weight capacity of 250 lbs.

SINGLED OUT:

NEW Single Post Manual Recline Back



- ✓ Convenient attendant or client releases offer infinite manual recline from 85-150 degrees
- ✓ Passive Gas Shock Assist ensures a smooth recline and return from reclined position
- ✓ The Single Post Manual Recline Back does not fold down

Removal of the dual back canes results in a smaller overall footprint and better turning radius while adding the comfort of angle adjustability and convenience of easier transport. Both styles are offered with our standard, steel Rehab Back, optional Matrix Contour or High Backs or a solid interface plate for ease of aftermarket back mounting.

Single Post Angle Adjustable Fold Down Back



- ✓ An attendant operated, manual release lever provides back angle adjustment from 91 to 127 degrees, in 6 degree increments
- ✓ Folds forward for easy transport in short ceilinged vehicles and many small commuter airliners
- ✓ Available as a no charge option with all Ultra Low CG tilt only systems and with nominal up-charge on F-Series tilt only systems

www.motionconcepts.com

usa
700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
Toll Free Tel: 1.888.433.6818
Toll Free Fax: 1.888.433.6834

canada
84 Citation Dr.,
Concord, ON L4K 3C1
Toll Free Tel: 1.866.748.7943
Tel: 905.695.0134
Fax: 905.695.0138


Motion Concepts