

Shown with seating system  
(open sideways).



## A life of EASE

**A**s a dealer or service technician, you'd probably just as soon forget the last time you tried to change a battery, access the motors or perform some type of service on a power base with a tilt system.

With most systems, this would usually involve a whole procedure, including removing a bunch of bolts, unplugging a bird's nest of wires and, after solving the problem, putting it all back together again.

Well, at Motion Concepts, we realize that life in the field is challenging enough without adding the need for lengthy service calls or tiresome installation manuals. So, we've taken care to design our EASE system with you in mind.

## EASE does it

Our exclusive EASE system (Easy Access Service Equipment) was created by our own highly regarded Senior Designer, *Son Ma* — who spends half of each week working in a dealership.

Standard on almost everything we build, the EASE is designed to provide instant access to the battery or electronics for repair or service. Simply remove two easy-to-reach bolts, which will allow the full seating system to pivot open in order to access the wiring, batteries and inner workings of the power base.

Depending on your particular base configuration, you can access the EASE by flipping your entire system forwards, backwards, or in some cases, sideways. So, you only need one service person to lift the complete seating system.

This is just one of the unique enhancements that makes Motion Concepts and our products truly unforgettable.

**U.S.**  
700 Ensminger Rd.,  
Suite 112  
Tonawanda, NY 14150  
TEL 1.888.433.6818  
FAX 1.888.433.6834

**Canada**  
84 Citation Dr.,  
Concord, ON  
L4K 3C1  
TEL 905.695.0134  
FAX 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)



Easy access to battery  
and chair electronics



Shown with seating system  
(open backwards).