

✓ **MaTRx Flo-tech cushion**

### Go with the Flo ...

*Cushions Shaped Around People*

Sometimes, even when you think your mind is made up, you just gotta 'go with the Flo'. Especially when the 'Flo' in question is the MaTRx Flo-tech cushion — the market leading, fluid pressure reduction cushion designed exclusively for users with the highest risk of skin breakdown.

The Flo-tech's unique, dual layer fluid sacs are specially designed to provide immersion to sensitive bony prominences and prevent "bottoming out".

#### Features include:

- ▶ Comfortable, molded foam base with pre-ischial contour and trochanteric shelf
- ▶ Moisture resistant "Platilon" inner liner protects foam from moisture
- ▶ Stretchy, water resistant, breathable outer cover — includes hook and loop fasteners, lifting strap and side pocket
- ▶ Optional fluid pack supplements/add-on components available to accommodate individual postural needs
- ▶ Sizes from 15" - 20" W x 16" - 21" D; 300 lb. weight capacity
- ▶ 3 year warranty

Check out the new MaTRx Flo-tech cushion from Motion Concepts today. Your clients will see in no time why going with 'the flo' just might turn out to be one of the most comfortable seating experiences they've ever had.



**Concepts  
Motion**

*US:* 700 Ensminger Rd. Suite 112  
Tonawanda, NY 14150  
Toll Free Tel: 1.888.433.6818  
Toll Free Fax: 1.888.433.6834

*Canada:* 84 Citation Dr.,  
Concord, ON L4K 3C1  
Toll Free Tel: 1.866.748.7943  
Tel: 905.695.0134  
Fax: 905.695.0138

*www.motionconcepts.com*