

## No need to lie about it ...

... or lead anyone up the garden path when the truth is so obvious. No matter how complex your clients' needs may be, they would certainly like to be out enjoying their lives as much as possible. Now you can help them get wherever their travels take them simply, safely and comfortably.

The new Horizontal Tilt Power Positioning System from Motion Concepts features a unique, adjustable center portion with standard growable telescoping end supports to accommodate various shapes and sizes.

### Features include:

- ▶ Sturdy, weight sustainable, manual or optional power flip down transfer side rails
- ▶ Ability to drive in the prone, supine or recumbent position
- ▶ Center portion tilts up to 45°
- ▶ All 3 sections come in standard sizes and are fully width and depth adjustable:
  - ▶ Center body adjusts from 14" x 14" – 20" x 24" in 3" increments
  - ▶ End supports adjust from 14" x 13" – 20" x 24" in 3" increments
  - ▶ Telescoping transfer side rails adjust 8" x 18" – 8" x 24" in 3" increments
- ▶ Joystick can be positioned at desired spot on any of the three sections
- ▶ Single Function ¼ drive speed standard
- ▶ Optional 7" power elevation module
- ▶ Standard 20.25" Seat-to-floor height (STF)
  - ▶ Custom STF of 18.75" to 30" (N/A with elevating module)
- ▶ Optional Cushion: 2" thick (1" visco, 1" HR)
- ▶ Side storage pouch for extra storage
- ▶ Vent applications and custom sized end supports via RFQ

Check out the innovative Horizontal Tilt Power Positioning system, and lead your clients 'up the garden path' in ultimate comfort and security today.

## ✓ Horizontal Tilt



Alternative modules available, please call customer service for RFQ



US: 700 Ensminger Rd. Suite 112  
Tonawanda, NY 14150  
Toll Free Tel: 1.888.433.6818  
Toll Free Fax: 1.888.433.6834

Canada: 84 Citation Dr.,  
Concord, ON L4K 3C1  
Toll Free Tel: 1.866.748.7943  
Tel: 905.695.0134  
Fax: 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)