

Lateral Tilt

### It's time to make a lateral move ...

With all the focus on upward mobility, promotions and job sharing in the corporate world these days, sometimes the benefits of a lateral move can be overlooked and underappreciated.

That's why we created our specialized Lateral Tilt positioning module. Because when it comes to power positioning, regardless of where your clients are, Motion Concepts provides the expert solutions that help them find and maintain their most comfortable, stable position.

The Lateral Tilt enables a Tilt or Tilt/Recline system to tilt up to 15° laterally in either direction. This can assist in pressure reduction, aid in transfers and help alleviate digestive pressure and respiratory issues. From this laterally tilted position, users still retain the ability to tilt posteriorally and/or recline. This system also has a built-in centering function so the system automatically stops when the seat is horizontal to the base.

#### Features include:

- ▶ Lateral tilt up to 15° left and right
- ▶ For use with with UltraLow CG positioning systems on selected power bases:
  - Quantum 6000
  - Invacare 3G Storm
  - Invacare TDX
- ▶ Accommodates most front rigging options (70° swing away footrests are recommended)
- ▶ 20°-25° of Tilt and/or Tilt/Recline available while laterally tilted based on drive lock-out limits
  - ▶ 150lbs weight capacity
  - ▶ Must be ordered through RFQ (Request For Quotation)

Check out the Lateral Tilt from Motion Concepts today and help your clients see how valuable, comfortable and stable a lateral move can really be.



**Motion** Concepts

*US:* 700 Ensminger Rd. Suite 112  
Tonawanda, NY 14150  
Toll Free Tel: 1.888.433.6818  
Toll Free Fax: 1.888.433.6834

*Canada:* 84 Citation Dr.,  
Concord, ON L4K 3C1  
Toll Free Tel: 1.866.748.7943  
Tel: 905.695.0134  
Fax: 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)