

U.S.
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
TEL 1.888.433.6818
FAX 1.888.433.6834

Canada
84 Citation Dr.,
Concord, ON
L4K 3C1
TEL 905.695.0134
FAX 905.695.0138

www.motionconcepts.com

Get the REAL low down ...

Finally ... your clients can get the real low down. Whatever their reasons are for wanting to get lower, the all new, modular CG tilt Low-Pro system from Motion Concepts takes them lower than they ever thought possible — often as much as 2" lower than any other power positioning system available!

The North American made, modular Low-Pro tilts to 55° and provides 174° of recline — without sacrificing the seat-to-floor height. It features a more compact subframe for a sleeker profile and shorter seat depths. And for clients that need extra access, a low STF Power Elevating Seat Module is also available.

Plus, only the Low-Pro, is available with two unprecedented no-charge options: a pair of handy, automotive-style Armrest Storage Compartments and an attractive, stylish Single Post Back! Along with looking great, this Single Post Back offers a multitude of back angle adjustability without tools and folds forward for easier transportation and storage. So you get more for your money!

Like all Motion Concepts' seating systems, the Low-Pro features solid armrests that minimize 'play'. Plus, the tilt actuator is specially situated to provide better linkage position and make lifting easier, so there's less wear on your system — which means you get better mechanical advantage and longer life.

Best of all, the Low-Pro is bundled to fit specifically into workable funding categories. Even the order form has been specially configured to fit into the new HCPCS codes for Power Positioning Systems.

The eye-catching, technologically advanced Low-Pro makes your seating decision so easy, all that's left to consider is which base best suits your clients needs. That's because like every extraordinary seating system created by Motion Concepts, the Low-Pro fits seamlessly on all the major bases from Pride, Sunrise or Invacare.

So, now there's no reason not to give your clients the REAL low down. Check out the new Low-Pro power positioning system from Motion Concepts today.



Rattle-free armrests easily flip back for unobstructed transfers. Mechanical stop prevents armrest from falling all the way back.



Tilt only version offers clean, sleek look with no back canes.



Back easily folds forward, making system compact for easier transportation.

Low-Pro Power Tilt & Recline system shown on Quickie 646 with 17¾" STF Height.

