

August 2002

U.S.  
700 Ensminger Rd.,  
Suite 112  
Tonawanda, New York  
14150  
TEL: 1.888.433.6818  
FAX: 1.888.433.6834

Canada  
84 Citation Dr.,  
Concord, ON  
L4K 3C1  
TEL: 905.695.0134  
FAX: 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)



## Finally ... no need to articulate

**T**here's an old saying that states if you want something done right, do it yourself. Well, here at Motion Concepts, we couldn't agree more.

A few years ago, our chief designer, *Son Ma* — who spends three days each week working in a dealership — noticed a major void in the market for a versatile, yet truly durable legrest that could elevate client's knees and reduce the shearing caused when their knees get pushed upwards and their legs backwards. So, he decided to create one himself.


Now, we're thrilled to present the fruits of *Son's* labor — unique, new Pivot Plus legrests that work by matching the natural pivot point on your client's knees as closely as possible with the pivot point on the legrest.

This system offers precise adjustability — the extension on the hanger receiver is completely adjustable in and out, and up and down — so you get an unlimited combination of adjustments suitable for clients of all shapes and sizes.

As well, angle-adjustable footplates and depth-adjustable calf pads are included. And the Pivot Plus legrests fit all TRx Power Positioning Systems.

Best of all, these legrests are, hands down, the most durable legrests available. Heavy gauge steel construction and sturdy rubber bumpers, combined with rugged foot-plate rollers, help patients deftly handle shocks, door openings and much more.

Pivot Plus is all you could ask for in a legrest and more! Available in power or manual at \$795.00!



\*Color shown for presentation purposes only.