

U.S.
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
TEL 1.888.433.6818
FAX 1.888.433.6834

Canada
84 Citation Dr.,
Concord, ON
L4K 3C1
TEL 905.695.0134
FAX 905.695.0138

www.motionconcepts.com

A powerful back-up plan

Now your clients don't have to get their backs up ... unless they want to. Thanks to the new Power Sliding Back from Motion Concepts, clients can now experience an ideal backup plan that provides extra comfort along with the opportunity for passive movement.

For clients with neurological conditions such as ALS or quadriplegics, it can be difficult to stay comfortable in their chair. Shifting their body weight can often bring much needed relief. With the new Power Sliding Back, all it takes is the flick of a switch, and they're maneuvered backwards, allowing them to passively stretch. This can be an amazing tool for clients that need help to pull themselves up, reposition, or simply shift their body and weight.

Plus, the Power Sliding Back works independently of the recline function, offsets shearing and even allows clients to ease back itch and discomfort through passively repositioning themselves.

But that's not all! It also helps minimize and eliminate certain clients' tendency to go into spasms by offering them more flexible movement, along with the ability to better position alternative drive controls such as Sip and Puff systems or Chin Level Mounting Systems. All this means more comfort and flexibility than ever.

The Power Sliding Back is available for all standard Motion Concepts backs (both ABS and Contour Backs) on all combined tilt and recline systems. It comes in 3", 4" or 5" of slide, and provides much more versatility than the 3 - 1/4" of slide offered with mechanical sliding backs. With the Power Sliding Back, clients' speed is adjustable through the chair's electronics and is available as toggle or non-toggle (single or independent switch).

So go ahead. Help your clients stay more comfortable in their chairs with the help of their own powerful back-up plan. Ask us about the new Power Sliding Back today.

Power Sliding Back works independently of the recline function and can be synchronized with recline



Client can independently adjust the back height to their desired level



Height adjustment is available in 3", 4" or 5" slide options and is available for all standard Motion Concepts backs