

### A Terrific Back-Up Plan

At Motion Concepts, we're setting new standards! Our solutions-oriented products are industry renowned for their innovation and excellence. But we don't rest on our laurels. We constantly evaluate, update and improve every product we create to ensure we maintain these high standards for your clients. And our recently upgraded Standard Rehab Back is no exception.

We've replaced our original standard ABS curved back with an attractive, ergonomically shaped, new model. This new model is symmetrically designed and slightly angled so it's ideal for a wide variety of clients that require solid support.

Our new Rehab Back is not only much better looking than its predecessor, it's made of steel — as opposed to plastic. That means it's thinner, stronger and provides infinite height adjustability. Plus, it's easy to use. You can adjust it up and down even with the client in the chair. No screws, clamps or drills needed. Simply loosen the quick release bolts to slide up or down, find your client's ideal position and tighten. The new Back is now standard on all power positioning systems and can be used in conjunction with most configurations including recline (with and without ESR), and precline.

But that's not all! When your order our new Standard Rehab Back, the back pan comes predrilled to accommodate most standard headrest clamps, multirails and Extended Shear Reduction options. Another new feature is the foam cover. This is designed to extend almost 1" past the top and bottom of the metal pan to provide extra comfort. It also protects clients using tilt from any metal parts digging in to their bodies, and prevent the cushion from sliding down over time. The multirails available on the new Rehab Back can accommodate numerous accessories including lateral supports, chest harnesses and feeding devices.

The Standard Rehab Back is available in sizes 12" W x 12" H - 21" W x 20" H and made from rigid 14 gauge (.075" thick) steel. We've even improved the foam to a high density polyurethane that meets California safety bulletin #117. It also comes complete with a 100% polyester cover that's washable and breathable, surrounding the entire collection with a zipper. (The cover also meets California bulletin #117.)

So go ahead and see for yourself why even our standard items are anything but standard. Check out our new Standard Rehab Back today!

### Seating



New, ergonomically shaped steel back provides angled support and easy adjustability — even while clients are in their chair.

New .075" thick steel replaces thicker ABS plastic for a thinner, more attractive look.

Select either high density polyurethane foam cushion or visco elastic foam insert. Both meet California Safety Bulletin #117.



**US:** 700 Ensminger Rd.,  
Suite 112,  
Tonawanda, NY 14150  
Toll Free Tel: 1.888.433.6818  
Toll Free Fax: 1.888.433.6834

**Canada:** 84 Citation Dr.,  
Concord, ON L4K 3C1  
Tel: 905.695.0134  
Fax: 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)