

U.S.
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
TEL 1.888.433.6818
FAX 1.888.433.6834

Canada
84 Citation Dr.,
Concord, ON
L4K 3C1
TEL 905.695.0134
FAX 905.695.0138

www.motionconcepts.com

Some people just need to vent ...

While it seems there are always people who need to vent their feelings or views, at Motion Concepts we know that some people need to vent just to breathe. And for those who also need mobility assistance, we're proud to present a useful selection of accompanying vent trays.

When you check out our vent trays, the first thing you'll notice is that they stack up with the best. All three of our vent tray styles, The Gimble, The Articulating and the Compact Tray have been sized and created exclusively to meet the varying needs of your vent users.

The options

The Gimble — our most popular vent tray due to its unique ability to gimble, or pivot — is highly effective because it pivots freely under users' movement during both tilt and recline. The Gimble vent tray fits nearly all available vent types and measures in at 14" W x 15" D x 11" H. It's a flat, horizontal style tray that's often used on systems where the vent battery is integral with the base battery tray, as well as in situations where maximum ventilator protection is not necessary. It also features a rubber stopper to help cushion your client's vent movement. Full box protection systems are also available.

The Articulating style vent tray is secured by lower articulating rods, so it ensures that your client's vent box remains vertical during tilt and recline. This style of vent tray is mostly used to reduce your client's overall system length by supporting the vent vertically. Measuring in at 14.6" W x 10" D x 16" H, the Articulating vent tray is used most commonly for Bear and LP style vents. Other valuable features include access cut-outs and a convenient cord holder.

Our exciting, newest offering is the Compact style vent tray. This tray is specifically designed for use with the latest compact style vents. It's smaller size — it measures in at 11" W x 4 1/4" D x 13" H — makes it particularly handy in smaller spaces. The Compact vent tray is held firmly in place with a removable protection cover and features a clear, acrylic view window to monitor vent activity.

There's a suitable Motion Concepts battery tray for most systems. Each one is designed specifically to control and offset the vent's weight distribution and also to be more compact — so they help clients avoid damaging walls during turns and other maneuvers. So, go ahead and do what you can to help your clients vent today with a vent tray from Motion Concepts.



The Gimble style vent tray fits most vent types.



Articulating style tray keeps vent vertical during tilt/recline.



Protective box for newer compact style vents stays within the footprint of the chair.

Example shown of a Quickie P-220 Super Low with compact style vent tray.



Extra batteries always added under the base.