

Motion Concepts gives the letter "f" a whole new meaning!

Introducing the Motion Concepts "F" Series CG Tilt System

A **f**abulous, **f**ully capable and **f**unctional tilt system made especially **f**or today's unique **f**unding situations with all the **f**undamental **f**eatures and **f**lexibility you have **f**ound with Motion Concepts products.

Recent changes in the power mobility device codes have providers looking for economical and reliable solutions to all their power mobility needs. The "F" Series delivers an efficient solution to all their power tilt needs with a system that is long on adjustability and short in price.

"F" Series — The only power tilt system in the industry with lifetime warranties and the ability to interface with every major power base on the current PMD Classification List.

Features:

- ▶ 55° of CG Tilt with 7" CG shift
- ▶ Seat to floor heights starting at 16.5" (depending on power base)
- ▶ Mounting Brackets for base to base interchangeability

System includes:

- ▶ Slotted seat pan and steel back pan
- ▶ Cantilever height and angle adjustable flip back armrests
- ▶ Slotted telescoping legrest receivers offer infinite seat depth adjustment within a 3" range
- ▶ Standard electronics with single push button, dual push button or low amp toggle control
- ▶ Standard fixed center mount foot platform
- ▶ Compatible with all front rigging
- ▶ Elevating seat available



USA: 700 Ensminger Rd. Suite 112
 Tonawanda, NY 14150
 Toll Free Tel: 1.888.433.6818
 Toll Free Fax: 1.888.433.6834

www.motionconcepts.com
info@motionconcepts.com