



Dear MaTRx Customer:

Motion Concepts is excited to introduce the new improved MaTRx/Kid\*ab\*ra Vi cushions/bases. These highly popular products continue to provide the support and skin protection which clients have relied on, as well as many new enhancements:

- 1) **Superior pressure re-distribution** – the new 2 part top layer provides targeted skin protection to the ischial/sacral area, while its unique pocketed design ensures that there are no issues with transition lines.
- 2) **Improved support and stability** – the High Resilient top layer provides better support to the hips and thighs for all types of users, even those who tend to “front load” the cushion.
- 3) **Lighter weight** – the Vi weighs in at approx. 2.5 lbs., almost 1 lb. lighter than the previous V model.
- 4) **Decreased temperature/humidity sensitivity** – the High Resilient top layer provides consistent performance in all types of climatic conditions.
- 5) **Enhanced comfort** – most clients find the new Vi the most comfortable cushion ever!

Due to the ongoing supplier issue with our MaTRx-V foam (pink), we are now filling all current MaTRx-V/Kid\*ab\*ra cushion/base orders with the new and improved Vi model. All dimensions are identical to the MaTRx-V – the only change is the foam composition of the top layer, with its targeted Visco insert pocketed within the High Resilient foam. There are no changes to the pricing or funding numbers for the new Vi models.

Please feel free to contact me at any time with questions or concerns.  
Thank you for your ongoing support of the MaTRx seating line.

Kind Regards,

A handwritten signature in black ink that reads 'Judy Rowley'.

Judy Rowley  
Vice President - MaTRx Seating  
Motion Concepts



For more information call toll free 1.888.433.6818 or [www.motionconcepts.com](http://www.motionconcepts.com) in USA  
1.800.668.5324 or [www.invacare.ca](http://www.invacare.ca) in Canada