



MATRIX SEATING SERIES
ACCESSORIES AND HARDWARE ORDER FORM
EFFECTIVE JANUARY 2010



Rev 05.11.11

Date of Order:	_____	Ship to Address:	_____
Dealer Name:	_____		_____
Dealer Acct#:	_____		_____
Contact:	_____		_____
Telephone:	_____		_____
PO #:	_____		_____
Tag:	_____		_____

Arm Accessories

Quantity	Model #	Accessory	Description	Retail Price	Total
	MPA-D	Multi Position Arm Pad (each)	Desk Length	\$85	
	MPA-F	Multi Position Arm Pad (each)	Full Length	\$85	
	GAP-N	Gel Arm Pad (each)	2" x 12"	\$75	
	GAP-W	Gel Arm Pad (each)	3.5" x 12"	\$75	
	GAP-L	Gel Arm Pad (each)	3.5" x 14"	\$90	
	EAH-R	Ergonomic Arm Trough/ Hand Pad kit* (each) HCPCS E2209	User's Right	\$210	
	EAH-L	Ergonomic Arm Trough/ Hand Pad kit* (each) HCPCS E2209	User's Left	\$210	
	EAT-R	Ergonomic Arm Trough** (each) HCPCS E2209	User's Right	\$100	
	EAT-L	Ergonomic Arm Trough** (each) HCPCS E2209	User's Left	\$100	
	EHP	Ergonomic Hand Pad*** (each)	User's Right / Left	\$85	
	EMH	Arm Trough/Hand Pad Mounting Hardware	For Motion Concepts Hand Pad	\$38	
	EMO		For Otto Bock Hand Pads	\$38	
	MACES-EA	Multi Axis Armrest Mount (each) For use with EAH, EAT or MPA	Allows infinite adjustment of angle and rotation of arm pad	\$225	
	MACES-GA	Multi Axis Armrest Mount (each) For use with GAP-W and GAP-L	Allows infinite adjustment of angle and rotation of arm pad	\$225	
	JOY-MACES-EA	Joystick mounting bar for Multi Axis Armrest Mount with EAH, EAT or MPA armpad		\$110	
	JOY-MACES-GA	Joystick mounting bar for Multi Axis Armrest Mount with Gel armpad		\$110	

* HCPCS K0106 - for Ergonomic Arm Trough Only

** Does not include mounting hardware to connect to hand pad - see EMH or see EAH-L vs. EAH-R.

** For arm trough connecting hardware please order EMH or see EAH-L vs. EAH-R for complete kit.



MATRIX SEATING SERIES
ACCESSORIES AND HARDWARE ORDER FORM
EFFECTIVE JANUARY 2010



Rev 05.11.11

Hardware

Quantity	Model #	Hardware	Description	Retail Price	Total
	BDH0	Base Drop Hooks (per pair)	Flushmount(0") drop for 7/8" or 1" diameter	\$27	
	BDH1-78	Base Drop Hooks (per pair)	1" drop for 7/8" diameter	\$27	
	BDH2-78		2" drop for 7/8" diameter	\$27	
	BDH1-1	Base Drop Hooks (per pair)	1" drop for 1" diameter	\$27	
	BDH2-1		2" drop for 1" diameter	\$27	
	PBMH*	Posture Back Mounting Hardware (per pair)	Complete Mounting assembly for Posture Back or Elite Back	\$225	
	PBMH-Q*	Posture Back Mounting Hardware (per pair)	Complete Mounting assembly for Posture Back or Elite Back 1 1/8" back canes.	\$225	
	SMH*	Split Mounting Hardware (Per pair)	Split Mounting hardware assembly for Posture Back or Elite Back.	\$250	
	HDMH*	Heavy Duty PB Mounting Hardware(per pair)	Mounting Hardware Assembly for Heavy Duty Posture Back	\$325	
* If purchasing replacment hardware and only one side is required please indicate "L" or "R" and split retail price in half					
	PBCR	Posture Back Cane Replacement (per pair)	1" diameter back cane replacement kit	\$115	
	FSSH	Fixed Space Saver Hardware	Allows fixed, non-removable mounting of Elite back or Posture Back**	\$200*	
	FSSH-FT	Fixed Space Saver Hardware - Formula Tilt	Fixed, non-removable mounting of Elite back or Posture back**	\$200*	
	FSSH-FTR	Fixed Space Saver Hardware - Formula Tilt/Recline	Fixed, non-removable mounting of Elite back or Posture back**	\$200*	
	ERMH	Add: Extended Range Mounting Hardware to existing hardware	Allows back to be positioned 2" further forward or behind back canes	\$60	
	ESRMH	ESR Mounting Hardware	Allows PB Back to be mounted to Motion Concepts TRx-Powered Seating with ESR	\$260	
	ESRMH-HD	ESR Mounting Hardware	Allows Heavy Duty PB Back to be mounted to Motion Concepts TRx-Powered Seating with ESR	\$310	
	BDC	DuraClip (per pair)	Top lock retractable hardware for use with drop hooks	\$60	
	BFC	FixClip (per pair)	Top lock hardware for use with drop hooks	\$25	
	EHRA-S	Elite Headrest Adapter - Small	For attachment of headrest to Elite Backs 12"wide	\$35	
	EHRA-M	Elite Headrest Adapter - Medium	For attachment of headrest to Elite Backs 13"-17"wide	\$35	
	EHRA-L	Elite Headrest Adapter - Large	For attachment of headrest to Elite Backs 18"-20"wide	\$35	
	EMP*	Extended Mounting Pins	Allows MaTRx Elite back or Posture Back to be installed on wheelchairs that are 2" wider than the back model width.	\$60	

*Not available on Heavy Duty Backs, weight capacity up to 300 lbs.

**Use this hardware to mount PB and Elite backs to Formula Wheelchairs



MATRIX SEATING SERIES
ACCESSORIES AND HARDWARE ORDER FORM
EFFECTIVE JANUARY 2010



Rev 05.11.11

Calf Straps and Lap Belts

Quantity	Model #	Accessory	Description	Retail Price	Total
	BP-CS-M	Air mesh Calf Strap - Medium HCPSC E0995	3.5" pad height, fits chairs 15"-18" wide	\$75	
	BP-CS-L	Air mesh Calf Strap - Large HCPSC E0995	3.5" pad height, fits chairs 18"-22" wide	\$75	
	BP-CP-MD	Air mesh Calf Panel - Medium HCPSC E0995	35.5"W x 9"H pad height, fits chairs 15"-18" wide	\$125	
	BP-CP-LG	Air mesh Calf Panel - Large HCPSC E0995	38.5"W x 9"H" pad height, fits chairs 18"-22" wide	\$125	
	BP-PLB2	Padded Lap Belt - 2 point HCPSC E0978	Center-pull, push button, padded lap belt fits hip width 13-19"	\$120	
	BP-PLB4	Padded Lap Belt - 4 point HCPSC E0978	Center-pull, push button, padded lap belt fits hip width 13-19"	\$150	