

TRx ULTRA-Low Bariatric Power Positioning System

US PRICE LIST AND ORDER FORM


Effective July 2008



Date of Order: _____
 Dealer Name: _____
 Dealer Account #: _____
 PO #: _____
 Base MFR. PO #: _____
 Tag: _____
 Purchasing Contact: _____
 Phone: _____ Fax: _____
 RTS/Therapist: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____
 Client Info: Gender: _____
 Special Client Conditions: _____
300-500 lbs weight limit. Dependent on weight capacity of the power base.

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Seat Cushion Thickness:	



INVACARE

BASE TYPE SELECTION

TDX SR HD Select STF⁽¹⁾: 18.25"
400 lbs weight capacity

TDX SP* *must order base with 400 lb weight capacity
Select STF⁽¹⁾: 18.25"
400 lbs weight capacity

3G Arrow Select STF⁽¹⁾: 18.5"

3G Torque Select STF⁽¹⁾: 18.5"

3G RangerX Select STF⁽¹⁾: 18.5"
400 lbs weight capacity

M94 Select STF⁽¹⁾: 18.5"
450 lbs weight capacity

ELECTRONICS

Mark VI PSS with ACC SPJ MPJ . . .N/C

Mark V NX™80 with SPJ-80N/C

Select: EX Controller SPJ DPJ MPJ

Mark IVR2 Mark IV Mark IVA . . .N/C

Select Drive: Right Left Specialty Control*

*Please specify: _____

OTHER INFORMATION

Short Base Long Base

4 Pole Motors GB Motors

QUANTUM

BASE TYPE SELECTION

1170 XL Select STF⁽¹⁾: 20.5"
400 lbs weight capacity

6000 XL Select STF⁽¹⁾: 17.25"
400 lbs weight capacity

1420* Select STF⁽¹⁾: 19.75"
450 lbs weight capacity

1650* Select STF⁽¹⁾: 20.5"
(Power Pivot Plus legs will increase STF to 21")
500 lbs weight capacity

614 HD Select STF⁽¹⁾: 18.25" 19.25"
400 lbs weight capacity

*Batteries are required at the time of installation.

ELECTRONICS

VSI 6-Key* VSI 7-Key P&G . . .N/C

VSI Integral 4-key R-NetN/C

Q-LogicN/C

(* Interface is required, see electronic options)

Select Drive: Right Left Specialty Control*

*Please specify: _____

Please note that the HCPCS codes indicated on this order form are our interpretation of the new codes and how we believe our products should be coded. The individual DMERCs will determine if our interpretation is appropriate and whether or not items coded will be accepted or denied. Motion Concepts cannot guarantee that use of the designated codes will ensure coverage or payment. Contact the insurer to confirm compliance with coding and coverage policy.

QUICKIE

BASE TYPE SELECTION

Rhythm Select STF⁽¹⁾: 18.75"
400 lbs weight capacity

S646 Select STF⁽¹⁾: 19.25"
400 lbs weight capacity

ELECTRONICS

Qtronix VSI 6-keyN/C

Delphi Select Joystick: QR* QC . . .N/C

* You must select QR to run the functions through the joystick.

Select Drive: Right Left Specialty

Control*

*Please specify: _____

OTHER INFORMATION

For S646 Bases select:

Standard Width(18") Wide Width (20")

BASES MUST BE ORDERED "TRx" READY FROM MANUFACTURER

*** ON BOARD CHARGERS ARE NOT RECOMMENDED FOR REHAB APPLICATIONS AS THE MAIN CHARGING SOURCE***

SYSTEMS⁽²⁾ MUST SELECT ONE SYSTEM

(Includes seat pan, back pan, armrests, armpads, electronics and standard switch controls)

TRx Ultra-Low Bariatric CG Tilt Only System (45°Tilt)⁽³⁾⁽⁴⁾\$7,395

Note: Tilt range is 45° up to 400 lbs, Tilt range is 30° for 400 lbs and greater

TRx Ultra-Low Bariatric CG Tilt (45°Tilt)⁽³⁾⁽⁴⁾ and Power Recline (135°) System\$11,390

Note: Tilt range is 45° up to 400 lbs, Tilt range is 30° for 400 lbs and greater

TRx Ultra-Low Bariatric CG Tilt (45°Tilt)⁽³⁾⁽⁴⁾ and Power Recline (135°) System with Extended Shear Reduction (ESR)\$11,955

Note: Tilt range is 45° up to 400 lbs, Tilt range is 30° for 400 lbs and greater

TRx Bariatric Recline Only System (135°)\$6,995

TRx Bariatric Recline Only System (135°) with Extended Shear Reduction\$7,560

MODULES - ONLY AVAILABLE ON SYSTEMS WITH RECLINE

TRx 30° Precline Module⁽⁵⁾\$500

TRx 8° Precline Module⁽⁵⁾N/C

STANDARD SEAT SPECIFICATIONS

Seat Width

19" 20" 21" 22" 23" 24" 25"N/C

Seat Depth

16" 17" 18" 19" 20"N/C

21" 22" 23" 24"N/C

CUSTOM SEAT SPECIFICATIONS

Seat Width

26" and overCall for RFQ

SEAT CUSHION OPTIONS

MaTRx-V Seat Cushion 20"w-30"w\$625

MaTRx-Flo-Tech Cushion 20"w-24"w\$595

"ON-CHAIR" Cushion Credit(\$90)

FABRIC OPTIONS FOR MATRX-VI CUSHION

Onyx Fabric Upgrade for MaTRx-Vi Cushion . . . \$100

Vinyl Fabric Upgrade for MaTRx-Vi Cushion . . . \$50

RECLINE BACK OPTIONS

TRx Recline Body Fit System (Minimum back widths of 16") . . . \$995

(This option allows the seat back to be narrower than the seat pan. For example an 18" back pan width with a 20" seat width is accomplished by selecting "2" Narrower" option.)

2" Narrower 4" Narrower

TILT ONLY BACK CANE SPECIFICATIONS

Select Back Cane Height - Measured from Seat pan to top of back cane

22" 24"N/C

Select Back Cane Angle 76° 79° 86° 96° 101° 112° 116°

Custom Cane Height Specify _____\$250

Angled (5°)Tilt Back Canes⁽⁶⁾N/C

Select Back Cane Angle 81° 84° 91° 101° 106° 117° 121°

Endomorph Back Canes Select: 1" 2" 4" Narrower N/C

FINISHED BACK HEIGHT*

Measured from seat pan to top of back shell

18"⁽⁷⁾ 19"⁽⁷⁾ 20" 21" 22" 23" 24" 25"

26" 27" 28" 29" 30" 31" 32" 33"

34" 35" * Comes with standard cane height

BACK OPTIONS

Standard Rehab Back and Back cushion
Finished height from 18" to 25"STD
Finished height from 26" to 35"\$450

MaTRx Fabric Upcharge -for Standard rehab back\$100

Recessed planar interface plate⁽⁸⁾ N/C

OMIT: Rigid Back Pan & Cushion⁽⁹⁾(\$75)

OMIT: Back Cushion Only(\$25)

UPGRADEABLE BACK OPTIONS

MaTRx Contour Back (HCPCS Code E2613)\$495
(For Systems 15"-21" wide under 300lbs)

Specify actual back height dimension: 12" 15" 17" 21"*
* only available in 15, 17, 18 and 20" widths

MaTRx Contour Back (HCPCS Code E2614)\$595
(For Systems 22" wide under 300lbs)

Specify actual back height dimension: 15" 17"

MaTRx High Back (HCPCS Code E2620)\$595
(For Systems 16"-21" wide under 300lbs)

Specify actual back height dimension: 20" 25" 30"

MaTRx PB Heavy Duty Back (5" Contour)\$835
(For Systems 20"-24" wide and over 500lbs)

Specify actual back height dimension: 16" 20"

MaTRx PB Heavy Duty Back (5" Contour)\$1,435
(For Systems 26"-30" wide and over 500lbs)

Specify actual back height dimension: 16" 20"

MaTRx PB Heavy Duty Deep Back (7" Contour)\$935
(For Systems 20"-24" wide and over 500lbs)

Specify actual back pan dimensions: 16" 20"

MaTRx PB Heavy Duty Deep Back (7" Contour) ...\$1,580
(For Systems 26"-30" wide and over 500lbs.)

Specify actual back pan dimensions: 16" 20"

"ON-CHAIR" Matrix Back Credit(\$125)

Onyx Fabric Upgrade\$100

Select Contour Back High Back

NOTE WHEN CHOOSING MaTRx Backs
4" Gap recommended for Tilt systems to allow for the seat cushion.
6" Gap required for systems with Recline/ESR or PSB
* PB Backs are not available with the Precline Module.

HEADRESTS & MOUNTING HARDWARE

- Headrests with Fixed Mounting Hardware** (HCPCS code E0955)
- Concepts Standard Headrest \$295
 - Combination Head/Neck - Small \$295
 - Combination Head/Neck - Large \$295
 - Stealth Comfort Plus 10" (CP250) \$242
 - Stealth Comfort Plus 14" (CP150) \$242
 - Onyx Auto Style Headrest - Straight Mounting available only . . \$350
- Headrest Swing Aways and Mounts - Does not include the Headrest** (HCPCS code E1028)
- Concepts Detachable Mounting Hardware \$200
 - Stealth's World's Best Headrest Hardware (TWB485) . . \$220
 - OMIT: Headrest N/C

ARMRESTS

- Adult Dual Post Adjustable Height (Single Button Release) . . STD
 - Add Quick Height Adjustment Lever \$250
 - Select Size: (9.75"-12.75") (12.75"-15.75") (15.75"-18.75")
 - Select: Flat pouch Glove Boxes*
 - (* Not Available with Mesomorph Back Canes)
- Reclining Adjustable Height STD
 - Select Size: (9.5"-13") (13.5"-18.5")
- 1" Offset Armrest Spacer (pair) \$250
(Allows seat width to increase by 2")
- 1/2" Offset Armrest Spacer (pair) \$250
(Allows seat width to increase by 1")

ARMPADS

- Standard Pads- STD
 - Select: Full (2.25"x14") Left Right
 - Desk (2.25"x10") Left Right
- Waterfall Arm Pads-(Pair) \$35
 - Select: Full (2.25"x14") Left Right
 - Desk (2.25"x10") Left Right
- Flat Pads- (Pair) \$55
 - Select: Full (3"x14") Left Right
 - Desk (3"x10") Left Right
- Gel Pads- (Pair) \$150
 - Select: Narrow (2"x12") Left Right
 - Wide (3.5"x12") Left Right
- Ergonomic Arm Trough - LEFT (4.5"x13") \$85
(HCPCS Code K0106)
- Ergonomic Arm Trough - RIGHT (4.5"x13") \$85
(HCPCS Code K0106)
- Flat Hand Pad - LEFT (4.5"x5") \$55
(HCPCS Code K0108)
- Flat Hand Pad - RIGHT (4.5"x5") \$55
(HCPCS Code K0108)
- OMIT: Armpads (7/8" diameter tube) (\$25)
 - Select: Full size (14") Left Right
 - Desk size(10") Left Right
- Horizontal Armpad Extensions (each) Quantity () . . \$75
 - Select: 3" 4.375" 7"

ARMPADS ACCESSORIES

- TRx Arm Pad Storage/Cup Holder LEFT * \$175
- TRx Arm Pad Storage/Cup Holder RIGHT * \$175
*Not available on the same side as the joystick is mounted
- TRx Elbow Block LEFT (prevents arm from slipping off) \$195
 - Select Pad Size: XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Vertical Extension: S(3") L(5")
 - Horizontal Extension: S(5") L(7")
- TRx Elbow Block RIGHT (prevents arm from slipping off) \$195
 - Select Pad Size: XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Vertical Extension: S(3") L(5")
 - Horizontal Extension: S(5") L(7")
- Multi Purpose Arm Pad Adapter Block (each) Quantity () \$35
(Allows hard mounting & positioning of standard pads or other manufacturer's arm pads. Adapter allows square mounting on round tubes)
- Flip & Fold Away 1/2 Lap Tray (14"W x 10"D) \$795
(Only available with a flat arm pad)
 - Select Mounting*: Left Right
 - *Not available on the same side as the joystick is mounted

LATERALS - FOR ALL OTHER BACKS (EXCLUDING MATRX PB ELITE AND PB BACKS)

- TRx Lateral Trunk Support with Fixed Mounting LEFT⁽¹⁰⁾
 - Pad Size: S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Offset 0"-3" 3"-6"
 - (HCPCS code E0956) \$100
- TRx Lateral Trunk Support with Fixed Mounting RIGHT⁽¹⁰⁾
 - Pad Size: S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Offset 0"-3" 3"-6"
 - (HCPCS code E0956) \$100
- Onyx Fabric Upgrade \$100

**SELECT SWING AWAY HARDWARE TO CONVERT ABOVE:
NOT AVAILABLE ON MATRX PB ELITE, PB ELITE FULL, AND MATRX PB BACKS**

- TRx Swingaway Mounting Hardware⁽¹⁰⁾
 - LEFT (HCPCS code E1028) \$195
 - RIGHT (HCPCS code E1028) \$195
 - Select Lateral Latch plate length:
 - Right: Short (3.5") Long(5") Left: Short (3.5") Long(5")

STEALTH LATERAL OPTIONS

- TWB-L Stealth Laterals⁽¹⁰⁾(pair) \$590
 - Pad Size: S(3"x5") L(4"x6")

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIP SUPPORTS

- Lateral Hip Support Pad with Fixed Mounting LEFT . . \$120
Select Size:
 (5.5"W x 4.5"H) (7"W x 5.75"H) (9"W x 4.5"H) (13"W x 4.5"H)
(HCPCS code E0956)
- Lateral Hip Support Pad with Fixed Mounting RIGHT . \$120
Select Size:
 (5.5"W x 4.5"H) (7"W x 5.75"H) (9"W x 4.5"H) (13"W x 4.5"H)
(HCPCS code E0956)
- Onyx Fabric Upgrade \$100

HARDWARE TO CONVERT ABOVE HIP SUPPORT TO REMOVABLE/SWINGAWAY HIP SUPPORTS

- TRx Lift Off, Removable Hardware for Hip Supports
 - LEFT (HCPCS code E1028) \$195
 - RIGHT (HCPCS code E1028) \$195
- TRx Swingaway Removable Hardware for Hip Supports
 - LEFT (HCPCS code E1028) \$215
 - RIGHT (HCPCS code E1028) \$215
 - Select Mounting Position:
 Hip Knee or Both... Add \$215 per side

STANDARD FRONT RIGGING

CALL TO CONFIRM ON SEAT SIZES BELOW 16"W X 16"D

- Legrest Receivers Only STD

****Please Specify Front rigging to be used on TRx System ⁽¹⁾**

POWER CENTER MOUNT FOOT PLATFORM OPTIONS

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Immersion: (Add)	
Seat Pan to Footplate: (Total)	

*****CENTER MOUNTS ARE RESTRICTED TO SEAT DEPTHS 16" AND GREATER******

- CMB Belt Drive Center Mount Articulating Foot Platform ⁽¹²⁾⁽¹³⁾⁽¹⁴⁾
(HCPCS code K0108) \$2,595
Seat pan to footplate: - select one:
 Short 9"-11"- 3.5"-5" adjustable articulation
 Medium 11"-15"- 3.5"-7" adjustable articulation
 Long 14.5"-21"- 5"-7" adjustable articulation
- Standard with one piece calf pad
- Comes with a 12" x 10" footplate made of durable aluminum
1/4" thick with flexible heel loops

POWER FRONT RIGGING OPTIONS

- TRx Power Elevating Pivot Plus Legrests ⁽¹⁵⁾ \$2,490
(Includes Pivot Plus Legrests & Actuator)(HCPCS code E1010)
Select Function: Individual ⁽¹⁶⁾ Combined Synchronized W/Recline
Seat pan to footplate: S (12"-16") M (15"-19") L (18"-21")
Select Footplate: Regular (7.5"W x 9"D) Large (9"W x 11"D)
ADD: Power Actuator Quick Release \$350

ELEVATING LEGRESTS

- TRx Pivot Plus Manual Elevating Legrests (Pair) \$720
Seat pan to footplate: S (12"-16") M (15"-19") L (18"-21")
Select Footplate: Regular (7.5"W x 9"D) Large (9"W x 11"D)
- AT5544 Standard Elevating Legrests (Pair) \$495
- Jazzy Elevating Legrests (Pair) \$495
- Quickie Standard Elevating Legrests (Pair) \$495

ADDITIONAL FRONT RIGGING OPTIONS

90° FRONT RIGGING/CENTER MOUNTED OPTIONS NOT AVAILABLE ON STORM SERIES BASES

- TRx-Heavy Duty 70° Swingaway Footrests (Pair) ⁽¹⁷⁾ . . . \$495
Select Receivers: Straight Flared inward 1" Flared outward 1"
Seat pan to footplate: S (12"-16") M (15"-19") L (18"-22")
Select Footplate: S (4.25"W x 8"D) M (5.25"W x 8"D) L (6.25"W x 8"D)
- Dynamic Coil Springs (Pair) for TRx-Heavy Duty 70° Swingaway Footrests \$165
- Action 70° Swingaway Footrests (Pair) \$395
- Jazzy 70° Swingaway Footrests (Pair) \$395
- Quickie 70° Swingaway Footrests (Pair) \$395
- Adjustable 70°-90° Fixed Center Mount Foot platform ⁽¹⁸⁾ \$595
(Depth and Height Adjustable -- Comes with a 12" x 10" footplate made of durable aluminum and calf pads are included)
Seat pan to footplate: S (12"-16") M (15"-19") L (18"-22")
- Invacare Split Foot Platform Interface Bracket \$595
(Invacare CTMT-T Foot Platform is included)

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FOOT PLATFORM & FOOTPLATE OPTIONS

Upgraded Flip-up Foot Platform* (12"W x 10"D) . . . \$180

Bariatric Flip-up Foot Platform*(15"W x 10"D) \$250

Split Footplates* (6" x 10" each) pair \$275

(*For use with Adjustable 70°-90° Fixed Center Mount Foot platform -Made of durable aluminum 1/4" thick with flexible heel loops.)

Adjustable Angle Footplates(Each)(HCPCS code K0040)

Left:

S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) . . . \$125

Right:

S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) . . . \$125

Adjustable Angle Footplates with Cutouts for Caster Clearance (Each)(HCPCS code K0040)

Left Regular (7.5"W x 9"D) Large (9"W x 11"D) . . . \$125

Right Regular (7.5"W x 9"D) Large (9"W x 11"D) . . . \$125

Heel Loops (Each) QUANTITY () \$20
(HCPCS code E0951)

CALF PAD AND STRAP OPTIONS

Flat Calf pads \$275

Select: TRx HD Swing Away Footrests

Fixed Center Mount Foot platform

Curved Calf pads \$275

Select: TRx HD Swing Away Footrests

TRx Pivot Plus Legrests

Fixed Center Mount Foot platform

One piece calf pad for Fixed Center Mount Footplate . \$275

Custom Bariatric Calf Pads Call for RFQ

Airmesh Calf Strap for chairs 18"-22" wide-3.5" pad height
(HCPCS code E0978) BP-CS-L \$50

SWING AWAY FOOTREST MECHANISM

- Assists user to swing away 70° footrests
- Power version requires enhanced electronics
- Compatible with Motion Concepts TRx HD 70° Swing away footrests or Invacare Action Swing away footrests
- Footrests NOT included

TRx Power Swing Away Mechanism . .Suggested Retail \$1,995

.RFQ Required

Specify : Switch or Joystick (See page 7 for electronics)

TRx Manual Swing Away Mechanism . . .Suggested Retail \$995

.RFQ Required

Specify Lever Mounting Right or Left

ADDITIONAL OPTIONS AND ACCESSORIES

TRx Concepts Back Pack - Single \$125

TRx Concepts Back Pack - Double \$195

Transfer Handles (Height Adjustable) \$350

Select: 4"-6.5"

9"-11.5"

D-Loops (set of 2)* \$250

-D-Loops increase the seat to floor height by 3/8"

BODY POINT LAP BELTS

Padded Lap Belt 2 point - BP-PLB2 (HCPCS code E0978) . . \$120

Padded Lap Belt 4 point - BP-PLB4 (HCPCS code E0978) . . \$150

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SPECIFY DRIVE

Select: Left hand Right hand

ELECTRONICS SWITCH MOUNTING

Select: Left Right Do Not Mount

Specialty Controls Please call for specification sheet 1-888-433-6818

STANDARD ELECTRONICS - N/A WITH VENT TRAYS OR HEAD ARRAYS

(These electronics are equipped with drive lock-out. They do not include attendant control, tilt limits, and variable speed control(s).)

TRx Single-Function Control BoxSTD

Please Select One: Single Push Button Dual Push Button

Single Toggle Omit Switch*

(*Single-Function Control Box uses a phono jack connection. Compatible with other phono jack switches.)

HIGH END ELECTRONICS

(These electronics are equipped with drive lock-out, attendant control, tilt limits and variable speed control(s).)

SINGLE-FUNCTION SYSTEM - (IE: TILT ONLY)

TRx Enhanced Single-Function Control Box \$500

Select: Single Push Button Dual Push Button⁽¹⁹⁾

Single Toggle

Enhanced Single Function through Joystick* \$1,150

(HCPCS code E2310)

Additional Switch: Single Push Button Dual Push Button⁽¹⁹⁾ \$97

Single Toggle

DUAL-FUNCTION SYSTEM - (IE: TILT/RECLINE)

TRx Enhanced Dual-Function Control Box \$700

Select: Dual Push Button Quad Push Button⁽¹⁹⁾

Dual Toggle Four-Way Toggle

Enhanced Dual Function through Joystick* \$1,550

(HCPCS code E2311)

Additional Switch: Dual Push Button Quad Push Button⁽¹⁹⁾ . . \$97

Dual Toggle Four-Way Toggle

MULTI-FUNCTION SYSTEM - (IE: TILT/RECLINE/POWER LEGS)

TRx Enhanced 2-5 Function Control Box \$1,150

Select: Quad Push Button Four-Way Toggle

Add: Function Display* \$250

Add: Function Display/Step Scanner* \$400

* The display is operated through left/right motion on 4-Way Toggle or buttons 3 and 4 on the Quad Push Button.

Enhanced 2-5 Function through Joystick \$2,350

(HCPCS code E2311)

Includes: Function Display Yes No

Add: Function Display/Step Scanner* \$150

Additional Switch: Quad Push Button Four-Way Toggle . . . \$97

(*Motion Concepts provides the communication interface for Invacare Electronics and includes a function display)

ELECTRONICS OPTIONS - REQUIRES HIGH END ELECTRONICS

TRx M16 Splitter(HCPCS code E2399) \$125

(Allows for multiple switch options with phono plug connection)

TRx Feather Touch Switch - Single function \$97

Dime Size Quarter Size

TRx Feather Touch Switch Multi function \$135

Dual Dime Size Quad Dime Size

Stealth Green Egg (ESGRN) \$77

(For use with Single function electronics only)

TRx Power Shut Off Switch⁽²⁰⁾ - Toggle Switch Only . . . \$175

(HCPCS code E2399)

TRx Lazarus Jump Starter \$150

(HCPCS code A9270)

Accessory Port Power Supply \$195

(3 amp capacity)

TRx Lighting Effects - Through Switch \$350

TRx Front Safety Lights (pair)

Operated through switches \$495

Operated through Joystick \$845

TRx Rear Safety Lights (pair)

Operated through switches \$495

Operated through Joystick \$845

TRX Front & Rear Light Combo Package Credit* . . . (\$150)

*Purchase both sets and receive a credit

MEMORY SEATING OPTIONS

Single Function Memory Seating* Suggested Retail \$495

.RFQ Required

(For use with Enhanced electronics only - Not available to control Elevating seat modules)

*Please call for specification sheet 1-888-433-6818

Multiple Function Memory Seating* . . Suggested Retail \$995

.RFQ Required

(For use with Enhanced electronics only - Not available to control Elevating seat modules)

*Please call for specification sheet 1-888-433-6818

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

O2 HOLDERS

- Non-Articulating O² Gas Holder (Cylinder Style) \$895
 - Articulating O² Holder⁽²¹⁾ \$1,250
- Please Select One:
- O² Gas (Cylinder Style) Liquid O² (6.5" x 4.5")

BATTERIES

- MK U1 Battery (MU-1 SLD G) Each \$140
Quantity () (HCPCS code E2365)
- MK GROUP 22NF Battery (M22NF SLD G) Each
Quantity () (HCPCS code E2361) \$173
- MK GROUP 24 Battery (M24 SLD G) Each
Quantity () (HCPCS code E2363) \$215
- MK GROUP 34 Battery (M34 SLD G) Each
Quantity () \$215

INSTALLATION

- Installation by Motion Concepts⁽²²⁾ N/C
- Dealer Installation (Subject to approval) (\$100)

Total Retail Value \$ _____

NOTES

- 1 Seat-to-floor height is measured at the front of the seat pan based on a standard seat size with a 0° pre-tilt. Seat to floor heights may vary +/- .25" dependent on conditions of the power base. Suspension castors may affect seat-to-floor height up to an additional .5"-1".
- 2 Systems are designed for use with a 2"-3" cushion.
- 3 Tilt range is 45° up to 400 lbs and 30° for 400 lbs or greater.
- 4 Back angle is adjustable from -76° to 116° degrees in increments.
- 5 Total degrees of recline is 174 ie: if 30 degrees of precline is used, then 144 degrees of recline will be available.
- 6 This option allows the back angle to adjust from 81° to 121° in increments.
- 7 This back height is not recommended with Recline Systems, please call for further details. Power Sliding Back (PSB) is not available on these back heights.
- 8 Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
- 9 Not available with Recline/ESR or Recline/PSB.
- 10 Laterals are not available with the MaTRx PB Backs.
- 11 Standard front rigging will only be provided when specified. Please provide or order separately.
- 12 Includes a flip up foot platform. Increases seat-to-floor height.
- 13 This item is only warrantied for client weights up to 375 lbs.
- 14 90° front rigging will collide with front casters on the Storm 3G, S626/S646.
- 15 Under normal usage Power Legrest actuators tend to be louder than Tilt/Recline Actuators.
- 16 Individual legs count as two functions.
- 17 Heavy duty durable construction (4.5lbs ea.). Lift and swing away setup.
- 18 90° front rigging will collide with front casters on the Storm 3G.
- 19 Non-reversing relays.
- 20 Must order an Enhanced control box.
- 21 Articulating O2 holders are not standard configuration on Precline Systems. Please call for quotation.
- 22 Batteries may be required at time of installation, please call to confirm.