




Date of Order: \_\_\_\_\_  
 Dealer Name: \_\_\_\_\_  
 Dealer Account #: \_\_\_\_\_  
 PO #: \_\_\_\_\_  
 Base MFR. PO #: \_\_\_\_\_  
 Tag: \_\_\_\_\_  
 Purchasing Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 RTS/Therapist: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 Client Info: Gender: \_\_\_\_\_  
 Special Client Conditions: \_\_\_\_\_  
**250-300 lbs weight limit. Dependent on weight capacity of the power base.**

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
<b>A.</b> Seat to Shoulder:	
<b>B.</b> Trunk Depth:	
<b>C.</b> Chest Width:	
<b>D.</b> Knee to Back:	
<b>E.</b> Seat to Top of Head:	
<b>F.</b> Elbow to Hand:	
<b>G.</b> Seat to Elbow:	
<b>H.</b> Hip Width:	
<b>I.</b> Knee to Heel:	
Seat Cushion Thickness:	



**BASE TYPE SELECTION**

TDX SI Select STF<sup>(1)</sup>: 17.25" 18.25" 19.25"  
 TDX SR Select STF<sup>(1)</sup>: 17.25" 18.25" 19.25"  
 TDX SC Select STF<sup>(1)</sup>: 16.25" 17.25" 18.25"  
 TDX SP Select STF<sup>(1)</sup>: 17.25" 18.25" 19.25"  
 Select TDX-SP Battery Tray Group 22 Group 24  
 TDX 3, 4, 5  
 Select STF<sup>(1)</sup>: 17.25"\* 18.25" 19.25"  
 \*This STF is not available when using Invacare's Transport Ready kit.  
 3G Storm Select STF<sup>(1)</sup>: 18" 19" 20"  
 Arrow  Torque  RangerX  
 \* Requires custom battery tray modification and is not available with elevating seat. ^ This STF is not available when using a vent tray due to battery terminal clearance  
 M91 Select STF<sup>(1)</sup>: 17.75" 18.75" 19.75"  
 M71<sup>(2)</sup> Select STF<sup>(1)</sup>: 17.75" 18.75" 19.75"  
 M51 Select STF<sup>(1)</sup>: 18.255" 19.25" 20.255"  
**NEW** Select:  4 Post  Single Post

**NOTE:**  
The F-Series is available on all base models offered by the base manufacturers for Group 2, 3 and 4 base funding codes.

**ELECTRONICS**

Mark 6 Select:  SPJAP ^  SPJ-ACC  MPJ  
 Mark 5  NX™80/SPJ-80 ^ or  EX Controller  
 Select:  DPJ  MPJ  
 Mark 4  ^ For separate switch control only  
**Select Drive:**  Right  Left  Specialty Control\*  
 \*Please specify: \_\_\_\_\_

**OTHER INFORMATION**

Short Base  Long Base  Heavy Duty  
 G-Trac Module  4 Pole Motors  GB Motors  
 Transport Ready Option:  Yes  No

**BASE TYPE SELECTION**

6000/Z Select STF<sup>(1)</sup>: 17.5" 18.5" 19.5"  
 600 Select STF<sup>(1)</sup>: 17.25" 18.25" 19.25"  
 1400 Select STF<sup>(1)</sup>: 17.5" 18.5" 19.25"  
 1402 Select STF<sup>(1)</sup>: 17.5" 18.5" 19.25"  
 610 Select STF<sup>(1)</sup>: 16.75" 17.75"  
 614 Select STF<sup>(1)</sup>: 17.25" 18.25" 19.25"  
 1107<sup>(3)</sup> Select STF<sup>(1)</sup>: 15.75" 16.75" 17.75"  
 \*VSI electronics are not compatible you must upgrade to P & G or Q-Logic.  
 1113 ATS<sup>(2)</sup> Select STF<sup>(1)</sup>: 16.5" 16.75" 17.5"  
 1103 Ultra Select STF<sup>(1)</sup>: 16.5" 17.5"  
 R-4000 Select STF<sup>(1)</sup>: 18" 18.5" 19"

**NOTE:**  
The F-Series is available on all base models offered by the base manufacturers for Group 2, 3 and 4 base funding codes.

**ELECTRONICS**

VSI 6-Key  VSI 7-Key  P&G  
 VSI Integral 4-key  R-Net  Q-Logic  
 VR2 4-Key  VR2 6-Key  
 Q-Logic NE+  Q-Logic NEX  
**Select Drive:**  Right  Left  Specialty Control\*  
 \*Please specify: \_\_\_\_\_

Please note that the HCPCS codes indicated on this order form are our interpretation of the new codes and how we believe our products should be coded. The individual DMERCs will determine if our interpretation is appropriate and whether or not items coded will be accepted or denied. Motion Concepts cannot guarantee that use of the designated codes will ensure coverage or payment. Contact the insurer to confirm compliance with coding and coverage policy.

**BASE TYPE SELECTION**

Pulse Select STF<sup>(1)</sup>: 17" 18" 19" **NEW**  
 Rhythm\* Select STF<sup>(1)</sup>: 18.25" 19.25" 20.25"  
 F11 Select STF<sup>(1)</sup>: 17.75" 18.75" 19.75"  
 M11<sup>(3)</sup> Select STF<sup>(1)</sup>: 16.25" 17.25" 18.25"  
 S646 Select STF<sup>(1)</sup>: 17.75" 18" 19"  
 S525 Select STF<sup>(1)</sup>: 18" 19" 20"  
 For S525 Bases select:  Standard Width (18")  
 Wide Width (20")  
 P222  P220 Select STF<sup>(1)</sup>: 18.75"  
 For P222/P220 and S646 Bases select:  
 Frame Size (14"-18")\* Order this size for Motion 16"-18" seat widths  
 Frame Size (19"-20")\* Order this size for Motion 16"-20" seat widths  
 Please specify the frame size to Sunrise when ordering\*\* On P220/P222 Base - The harness (#139460) is required with the base for installation of the seating system when using VSI electronics.  
 On S525 Base - The harness (#910149) is required with the base for installation of the seating system when using VSI electronics.

**NOTE:**  
The F-Series is available on all base models offered by the base manufacturers for Group 2, 3 and 4 base funding codes.

**ELECTRONICS**

Qtronix  VSI 6-key  Delphi  
 Select Joystick:  QR\*  QC  
 \* You must select QR to run the functions through the joystick.

**BASE TYPE SELECTION**

CompassTRO Select STF<sup>(1)</sup>: 17.5" 18.5" 19.5"  
**Select Drive:**  Right  Left  Specialty Control\*  
 \*Please specify: \_\_\_\_\_

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**"RFQ" = Request For Quotation** - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

**SINGLE FUNCTION SYSTEM<sup>(4)</sup>**

*(Systems include a seat pan, back pan, armrests, armpads and front rigging)*

- F-Series CG Tilt Only System (55°Tilt)<sup>(6)</sup>**  
(HCPCS code E1002) .....\$4,700
- 5° Pre-Tilt\*\*** ..... N/C  
\*\*This option increases seat-to-floor height by approximately 1.5"-2" minimum  
Please call for confirmation.

**ELEVATING SEAT OPTION**

- TRx Elevating Seat Module <sup>(6)(7)</sup>**  
**NOT AVAILABLE ON M71, M51, 1107, 614, M11 AND COMPASS POWER BASES**  
Weight capacity is 175lbs on 1113, 610, 1103 Ultra  
Select Elevation:  
  - 7" Elevation - 225lbs max. ....\$2,550
  - 6"HD Elevation - 280lbs max. ....\$2,995**AVAILABLE ON F11, RYTHYM, TDX SP, STORM, Q610, Q600 AND Q6000** (lower STF elevating module)  
  - 6" Elevation - 225lbs max. ....\$2,995

\* The seat-to-floor height may increase .5" - 1.5" on the lowest available STF, see page 1 or call for specific details.

- Battery Tray Modification for TDX-SP, TDX-SR, TDX 3,4,5 & Storm 3, Q6000/Z to allow for Two(2) Group 24 Batteries** .....\$295
- Battery Tray for Two(2) Group 24 Batteries P220/P222 Only N/C**  
(Group 22 Batteries are required on elevating modules unless an alternate tray is selected)

**STANDARD SEAT SPECIFICATIONS - FIXED FRAME**

- Seat Width**  
 14"  16"  18"  20"<sup>(8)</sup> .....N/C  
 22"<sup>(9)</sup> .....\$595
- Seat Depth** (depth adjustable seat pan)  
 16"-19"  17-19"  18-20"<sup>(10)</sup>  19-21"<sup>(8)(11)</sup> ....N/C  
 22"<sup>(9)(11)</sup> .....\$595  
 Specify seat depth: \_\_\_\_\_

**SEAT CUSHION OPTIONS**

- MaTRx-Vi Seat Cushion 14"w-20"w** (HCPCS Code E2607) ..\$425
- MaTRx-Vi Seat Cushion 22"w** (HCPCS Code E2608) .....\$475
- MaTRx-Flo-Tech Cushion 14"w-22"w** (HCPCS CodeK0734) ..\$475
- "ON-CHAIR" Cushion Credit** .....(\$20)

**FABRIC OPTIONS FOR MATRX-VI CUSHION**

- Onyx Fabric Upgrade for MaTRx-Vi Cushion** ..\$100
- Vinyl Fabric Upgrade for MaTRx-Vi Cushion** ...\$50

**BACK OPTIONS**

- Standard Rehab Back and Back cushion**  
finished height for 18" to 25" .....STD  
finished height from 26" to 35" .....\$450
- MaTRx Fabric Upcharge -for Standard rehab back** ...\$100
- OMIT: Rigid Back Pan & Cushion** .....(\$75)
- OMIT: Back Cushion Only** .....(\$25)

**UPGRADEABLE BACK OPTIONS**

- MaTRx Contour Back** (HCPCS Code E2613) .....\$495  
(For Systems 15"-21" wide under 300lbs)  
Specify actual back height dimension:  12"  15"  17"  21"
- MaTRx Contour Back** (HCPCS Code E2614) .....\$595  
(For Systems 22" wide under 300lbs)  
Specify actual back height dimension:  15"  17"
- MaTRx High Back** (HCPCS Code E2620) .....\$595  
(For Systems 16"-20" wide under 300lbs)  
Specify actual back height dimension:  20"  25"  30"
- MaTRx PB Elite Back**(HCPCS Code E2620) .....\$625  
(For Systems 14"-20" wide under 300lbs)  
Specify actual back height dimension:  12"  16"  20"
- MaTRx PB Elite "Full" Back**(HCPCS Code E2620) .....\$625  
(For Systems 16"-20" wide under 300lbs)  
Specify actual back height dimension:  20"
- MaTRx PB Back (5" Contour)** (HCPCS Code E2615) ...\$625  
(For Systems 14"-20" wide under 300lbs)  
Specify actual back height dimension:  12"  16"  20"
- MaTRx PB Deep Back (7" Contour)** (HCPCS Code E2615) \$645  
(For Systems 14"-20" wide under 300lbs)  
Specify actual back height dimension:  12"  16"  20"
- MaTRx PB Heavy Duty Back (5" Contour)** .....\$835  
(For Systems 20"-22" wide and over 300lbs)  
(HCPCS Code E2615 for 20" width & HCPCS Code E2616 for 22" width)  
Specify actual back height dimension:  16"  20"
- MaTRx PB Heavy Duty Deep Back (7" Contour)** ...\$935  
(For Systems 20"-22" wide and over 300lbs)  
(HCPCS Code E2615 for 20" width & HCPCS Code E2616 for 22" width)  
Specify actual back pan dimensions:  16"  20"
- "ON-CHAIR" Matrx Back Credit** .....(\$50)
- Onyx Fabric Upgrade** .....\$100  
Select  Contour Back  High Back

**NOTE WHEN CHOOSING MaTRx Backs**  
4" Gap recommended for Tilt systems to allow for the seat cushion.

**For Options not on the F-Series system, please see the Ultra low order form or call Customer Service at 1-888-433-6818.**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**"RFQ" = Request For Quotation** - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

**TILT ONLY BACK CANE SPECIFICATIONS**

- Select Back Cane Height - Measured from Seat pan to top of back cane  
 22"  24" .....N/C  
 Select Back Cane Angle  76°  79°  86°  96°  101°  112°  116°  
 Custom Cane Height Specify \_\_\_\_\_ .\$.250  
 Angled (5°) Tilt Back Canes<sup>(13)</sup> .....N/C  
 Select Back Cane Angle  81°  84°  91°  101°  106°  117°  121°  
 Endomorph Back Canes Select:  1"  2"  4" Narrower .\$.250  
 Mesomorph Back Canes\* Select:  1"  2"  4" Wider .\$.265

\*Not available with Standard Dual Post height adjustable armrests. Must use Dual Post height adjustable armrests with Outback option or Trx Height & Angle Adjustable Cantilever Flip back armrests.

**FINISHED BACK HEIGHT\***

- Measured from seat pan to top of back shell  
 18"  19"  20"  21"  22"  23"  24"  25"  26"  
 27"  28"  29"  30"  31"  32"  33"  34"  35"  
 \* Comes with standard cane height

**HEADRESTS & MOUNTING HARDWARE**

- Headrests with Fixed Mounting Hardware (HCPCS code E0955)  
 Concepts Standard Headrest .....\$345  
 Stealth All Positioning Headrest Large (AP485) .....\$315  
 Stealth All Positioning Headrest X-Large (AP585) .....\$315  
 Nino Comfort Plus 6" (CPN350) .....\$315  
 Stealth Comfort Plus 8" (CP450) .....\$315  
 Stealth Comfort Plus 10" (CP250) .....\$315  
 Stealth Comfort Plus 14" (CP150) .....\$315  
 Onyx Auto Style Headrest - Straight Mounting available only ..\$435

Headrest Swing Aways and Mounts - Does not include the Headrest (HCPCS code E1028)

- Concepts Detachable Mounting Hardware .....\$250  
 Stealth's World's Best Headrest Hardware (TWB485) ..\$285  
 OMIT: Headrest .....N/C

**ARMRESTS**

- TRx Height & Angle Adjustable Cantilever Flip Back  
 Armrests (Minimum armrest height is 9.5")<sup>(14)</sup> .....STD  
 Add: Quick Lock Mechanism (each) .....\$100  
 Adult Dual Post Adjustable Height (Single Button Release) \$265  
 Add  Quick Height Adjustment Lever .....\$330  
 Select Size:  (9.75"-12.75")  (12.75"-15.75")  (15.75"-18.75")  
 Select:  Flat pouch  Glove Boxes\*  
 (\* Not Available with Mesomorph Back Canes)

**NEW** Add  Outback Armrest (each) Specify: Right  Left  .....\$285  
 Allows outward rotation while flipping back armrest

- 1" Offset Armrest Spacer (pair) .....\$250  
 (Allows seat width to increase by 2")  
 1/2" Offset Armrest Spacer (pair) .....\$250  
 (Allows seat width to increase by 1")

**ARMPADS**

- Standard Pads- .....STD  
 Select: Full (2.25"x14")  Left  Right  
 Desk (2.25"x10")  Left  Right

**ARMPADS - CONTINUED**

- Waterfall Arm Pads-(Pair) .....\$45  
 Select: Full (2.25"x14")  Left  Right  
 Desk (2.25"x10")  Left  Right  
 Flat Pads- (Pair) .....\$75  
 Select: Full (3"x14")  Left  Right  
 Desk (3"x10")  Left  Right  
 Gel Pads- (Pair) .....\$200  
 Select: Narrow (2"x12")  Left  Right  
 Wide (3.5"x12")  Left  Right  
 Long (3.5"x14")  Left  Right  
 Ergonomic Arm Trough - LEFT (4.5"x13") .....\$110  
 (HCPCS Code E2209)  
 Ergonomic Arm Trough - RIGHT (4.5"x13") .....\$110  
 (HCPCS Code E2209)  
 Flat Hand Pad - LEFT (4.5"x5") (HCPCS Code K0108) ...\$75  
 Flat Hand Pad - RIGHT (4.5"x5") (HCPCS Code K0108) ...\$75  
 OMIT: Armpads (7/8" diameter tube) .....(\$25)  
 Select: Full size (14")  Left  Right  
 Desk size(10")  Left  Right  
 Horizontal Armpad Extensions (each) Quantity ( ) ..\$95  
 Select:  3"  4.375"  7"

**ARMPADS ACCESSORIES**

- TRx Arm Pad Storage/Cup Holder LEFT \* .....\$225  
 TRx Arm Pad Storage/Cup Holder RIGHT \* .....\$225  
 \*Not available on the same side as the joystick is mounted  
 TRx Elbow Block RIGHT (prevents arm from slipping off) .....\$260  
 Select Pad Size:  
 XS(3.25"x4.25")  S(3.5"x5.25")  M(3.75"x5.75")  L(4.25"x6.5")  
 Select Extension Vertical:  S(3")  L(5") Horizontal:  S(5")  L(7")  
 TRx Elbow Block LEFT (prevents arm from slipping off) .....\$260  
 Select Pad Size:  
 XS(3.25"x4.25")  S(3.5"x5.25")  M(3.75"x5.75")  L(4.25"x6.5")  
 Select Extension Vertical:  S(3")  L(5") Horizontal:  S(5")  L(7")  
 Stealth Swing-Away Elbow Support - RIGHT .....\$365  
 Select Pad Size:  S(3" x 3")  M(4"x4")  L(3"x 10")  
 Stealth Swing-Away Elbow Support - LEFT .....\$365  
 Select Pad Size:  S(3" x 3")  M(4"x4")  L(3"x 10")  
 Multi Purpose Arm Pad Adapter  
 Block (each) Quantity ( ) .....\$35  
 (Allows hard mounting & positioning of standard pads or other manufacturer's arm pads. Adapter allows square mounting on round tubes)  
**NEW**  Multi Axis Upper Extremity Support (each) .....\$225  
 Mechanical locking bolt .....STD  
 Quick Adjustment Lever .....N/C  
 Allows rotation of the armpad

- Flip & Fold Away 1/2 Lap Tray (14"W x 10"D) .....\$995  
 (Only available with a flat arm pad)  
 Select Mounting\*:  Left  Right

\*Not available on the same side as the joystick is mounted

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

**FIXED LATERALS - FOR MATRX PB ELITE AND MATRX PB ELITE FULL BACKS**

- Standard Fixed LEFT - PBL-S-L (HCPCS code E0956) . . . . . \$100  
Select pad :  Standard (7.5" x 5")  Small (5.5"x3")
- Standard Fixed RIGHT - PBL-S-R (HCPCS code E0956) . . . . . \$100  
Select pad :  Standard (7.5" x 5")  Small (5.5"x3")
- Offset Fixed LEFT - PBOL-L (HCPCS code E0956) . . . . . \$110
- Offset Fixed RIGHT - PBOL-R (HCPCS code E0956) . . . . . \$110

**SWING AWAY LATERALS - FOR MATRX PB ELITE AND MATRX PB ELITE FULL BACKS**

- Lateral Trunk Support Left pad (each) (HCPCS code E0956) . . \$115  
Select Size:  
 XS(3.25"x4.25")  S(3.5"x5.25")  M(3.75"x5.75")  L(4.25"x6.5")
- Lateral Trunk Support Right pad (each) (HCPCS code E0956) \$115  
Select Size:  
 XS(3.25"x4.25")  S(3.5"x5.25")  M(3.75"x5.75")  L(4.25"x6.5")

**NEW**

**SELECT SWING AWAY HARDWARE TO CONVERT ABOVE:**

- Elite Back Swing Away Hardware
  - LEFT (HCPCS code E1028) . . . . . \$215
  - RIGHT (HCPCS code E1028) . . . . . \$215
- Offset Elite Swing Away Hardware
  - LEFT (HCPCS code E1028) . . . . . \$215
  - RIGHT (HCPCS code E1028) . . . . . \$215

**LATERALS - FOR ALL OTHER BACKS (EXCLUDING MATRX PB ELITE AND PB BACKS)**

- TRx Lateral Trunk Support with Fixed Mounting LEFT<sup>(15)</sup>  
Pad Size:  S(3.5"x5.25")  M(3.75"x5.75")  L(4.25"x6.5")  
Select Offset  0"-3"  3"-6"  
(HCPCS code E0956) . . . . . \$115
- TRx Lateral Trunk Support with Fixed Mounting RIGHT<sup>(15)</sup>  
Pad Size:  S(3.5"x5.25")  M(3.75"x5.75")  L(4.25"x6.5")  
Select Offset  0"-3"  3"-6"  
(HCPCS code E0956) . . . . . \$115
- Onyx Fabric Upgrade . . . . . \$100

**SELECT SWING AWAY HARDWARE TO CONVERT ABOVE:**

- TRx Swingaway Mounting Hardware<sup>(15)</sup>
  - LEFT (HCPCS code E1028) . . . . . \$215
  - RIGHT (HCPCS code E1028) . . . . . \$215
- Select Lateral Latch plate length:  
Right:  Short (3.5")  Long(5") Left:  Short (3.5")  Long(5")

**STEALTH LATERALS**

- FMTWB-L Stealth Laterals<sup>(15)</sup> -Fixed (pair) . . . . . \$230  
Pad Size:  S(3"X5")  L(4"X6") (HCPCS code E0956)

**SELECT SWING AWAY HARDWARE TO CONVERT ABOVE:**

- TWB-L Stealth Adjustable Swingaway Hardware(pair)  
(HCPCS code E1028) . . . . . \$440

**HIP SUPPORTS**

- Lateral Hip Support Pad with Fixed Mounting LEFT . . \$145  
Select Size:  
 (5.5"W x 4.5"H)  (7"W x 5.75"H)  (9"W x 4.5"H)  (13"W x 4.5"H)  
(HCPCS code E0956)
- Lateral Hip Support Pad with Fixed Mounting RIGHT . \$145  
Select Size:  
 (5.5"W x 4.5"H)  (7"W x 5.75"H)  (9"W x 4.5"H)  (13"W x 4.5"H)  
(HCPCS code E0956)
- Onyx Fabric Upgrade . . . . . \$100

**HARDWARE TO CONVERT ABOVE HIP SUPPORT TO REMOVABLE/SWINGAWAY HIP OR KNEE SUPPORTS**

- TRx Lift Off, Removable Hardware for Hip Supports
  - LEFT (HCPCS code E1028) . . . . . \$205
  - RIGHT (HCPCS code E1028) . . . . . \$205
- TRx Swingaway Removable Hardware for Hip Supports
  - LEFT (HCPCS code E1028) . . . . . \$265
  - RIGHT (HCPCS code E1028) . . . . . \$265

Select Mounting Position:  Hip  Knee or  Both... Add \$215 per side

Please note that the HCPCS codes indicated on this order form are our interpretation of the codes and how we believe our products should be coded. The individual DMERCs will determine if our interpretation is appropriate and whether or not items coded will be accepted or denied. Motion Concepts cannot guarantee that use of the designated codes will ensure coverage or payment. Contact the insurer to confirm compliance with coding and coverage policy.

For Options not on the F-Series system, please see the Ultra low order form or call Customer Service at 1-888-433-6818.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**"RFQ" = Request For Quotation** - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

**STANDARD FRONT RIGGING**

CALL TO CONFIRM ON SEAT SIZES BELOW 16" W X 16" D

Legrest Receivers Only .....STD

**\*\*Please Specify Front rigging to be used on TRx System (16)**

Fixed Center Mount Foot Platform  70°  90° .....N/C  
 (Comes with tapered footplate 12" x 10.75" made of steel 1/8" thick with built in heel wells)

Seat pan to footplate:  S (9"-15")  L (16"-21") Set at \_\_\_\_\_

OR

Action 70° Swingaway Footrests (Pair)\* .....N/C

**POWER CENTER MOUNT FOOT PLATFORM OPTIONS**

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Immersion: (Add)	
Seat Pan to Footplate: (Total)	

\*\*\*CENTER MOUNTS ARE RESTRICTED TO SEAT DEPTHS . . .16" AND GREATER\*\*\*\*\*

CMB Belt Drive Center Mount Articulating Foot Platform<sup>(17)(18)(19)</sup>  
 (HCPCS code K0108) .....\$2,595

- Short 9"-11"- 3.5"-5" adjustable articulation
- Medium 11"-15"- 3.5"-7" adjustable articulation
- Long 14.5"- 21"- 5"-7" adjustable articulation
- Standard with one piece calf pad
- Comes with a 12" x 10.75" footplate made of durable aluminum 1/4" thick with flexible heel loops

Add:  Individual Calf pads for CMB Belt Drive .....\$275

CMX Center Mount Articulating Foot Platform<sup>(17)(18)(19)</sup>  
 (HCPCS code K0108) .....\$2,995

- Seat pan to footplate: - select one:
- Short 14"-17"- 3"-5"
  - Medium 15"-19"- 5"-7.5"
  - Long 17"-21"- 5"-7.5"
  - Standard with individual calf pads
  - 6 3/4" extension
  - Plantar flexion and dynamic calf support
  - N/A with one piece calf pad
  - Comes with a 12" x 10" footplate made of durable aluminum 1/4" thick with flexible heel loops

**POWER FOOT PLATFORM OPTIONS**

Power Telescoping Foot Platform . . .Suggested Retail \$2,250  
 .....RFQ Required  
 - Comes with Adjustable 70°-90° Fixed Center Mount Foot platform<sup>(18)</sup>  
 - Does not elevate  
 - Powers lower leg length vertically

**POWER FRONT RIGGING OPTIONS**

TRx Swing Away Power Elevating Pivot Plus Legrests<sup>(20)</sup> \$2,750  
 (Includes Swing Away Receivers, Pivot Plus Legrests & Actuators - legs swing away manually)(HCPCS code E1010 )

Select Receivers:  Straight  Flared outward 1"

Select Function:  Individual<sup>(21)</sup>  Combined

Seat pan to footplate:  S (12"-16")  M (15"-19")  L (18"-21")

Footplates:  S(7.25"W x 9"D)  M(7.5"W x 11"D)  L(10"W x 14"D)

TRx Power Elevating Pivot Plus Legrests<sup>(20)</sup> ..... \$2,490  
 (Includes Pivot Plus Legrests & Actuator)(HCPCS code E1010 )

Select Function:  Individual<sup>(21)</sup>  Combined

Seat pan to footplate:  S (12"-16")  M (15"-19")  L (18"-21")

Footplates:  S(7.25"W x 9"D)  M(7.5"W x 11"D)  L(10"W x 14"D)

ADD: Power Actuator Quick Release .....\$350

TRx Power Elevating AT5544 Elevating Legrests<sup>(20)</sup> . \$1,990  
 (Includes AT5544 Legrests & Actuators)(HCPCS code E1010 )

Select Function:  Individual<sup>(21)</sup>  Combined

TRx Power Legrest Actuators (PAIR)<sup>(16)(20)</sup> .....\$1,695  
 (Legrests not included - For use with Elevating Legrests)

Select Function:  Individual<sup>(21)</sup>  Combined

**ELEVATING LEGRESTS**

TRx Pivot Plus Manual Elevating Legrests (Pair) .....\$720

Select Receivers:  Straight  Flared outward 1"

Seat pan to footplate:  S (12"-16")  M (15"-19")  L (18"-21")

Footplates:  S(7.25"W x 9"D)  M(7.5"W x 11"D)  L(10"W x 14"D)

ADD: Optional Swing away receiver .....\$150

AT5544 Standard Elevating Legrests (Pair) .....\$495

Jazzy Elevating Legrests (Pair) .....\$625

Quickie Standard Elevating Legrests (Pair) .....\$625

**90° FRONT RIGGING/CENTER MOUNTED OPTIONS NOT RECOMMENDED ON STORM SERIES, R-4000, 646 AND P220/222 BASES**

**For Options not on the F-Series system, please see the Ultra low order form or call Customer Service at 1-888-433-6818.**

**"RFQ" = Request For Quotation** - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

Dealer Name: \_\_\_\_\_  
 Dealer Account Number: \_\_\_\_\_  
 PO #: \_\_\_\_\_  
 Tag: \_\_\_\_\_

**ADDITIONAL FRONT RIGGING OPTIONS**  
**90° FRONT RIGGING/CENTER MOUNTED OPTIONS NOT RECOMMENDED ON STORM SERIES, R-4000, 646 AND P220/222 BASES**

- NEW**  TRx-Heavy Duty 70° Swingaway Footrests (Pair)<sup>(22)</sup> . . . \$495  
 Select Receivers:  Straight  Flared outward 1"  
 Seat pan to footplate:  S (12"-16")  M (15"-19")  L (18"-22")  
 Select Footplate:  S(4.25"W x 8"D)  M(5.25"Wx8"D)  
 L(6.25"Wx8"D)  XL(9"Wx11"D)
- Dynamic Coil Springs (Pair) for TRx-Heavy Duty 70° Swingaway Footrests . . . . . \$205
- Action 70° Swingaway Footrests (Pair) . . . . . \$515
- Jazzy 70° Swingaway Footrests (Pair) . . . . . \$515
- Quickie 70° Swingaway Footrests (Pair) . . . . . \$515
- Adjustable 70°-90° Fixed Center Mount Foot platform<sup>(18)</sup> \$595  
 (Depth and Height Adjustable -- Comes with a 12" x 10" footplate made of durable aluminum and calf pads are included)  
 Seat pan to footplate:  S (9"-15")  L (16"-21") Set at \_\_\_\_\_
- Invacare Split Foot Platform Interface Bracket<sup>(18)</sup> . . . . . \$795  
 (Invacare CTMT-T Foot Platform is included)  
 Seat pan to footplate:  S (9.25"-12")  M (12"-16")  L (16"-18.75")
- Tapered Manual Elevating Foot cage (Range 12"-18") \$1,275  
 Select side for adjustment lever:  Left  Right  
 Seat pan to footplate :

(For seat widths 16"-20")

**FOOT PLATFORM & FOOTPLATE OPTIONS**

- Small Flip Foot Platform\*(12"W x 8.5"D)..... \$235
- Medium Flip-up Foot Platform\* (12"W x 10.75"D) . . . \$240
- Large Flip-up Foot Platform\*(17"W x 13.75"D) . . . . . \$325
- Extra Large Flip-up Foot Platform\* ^ (21"W x 14.75"D) \$325
- Split Footplates\* . . . . . \$345

In cast aluminum with built in heel cups and rubber mat.

Select:  S(5"W x 7.5"D)  M(5.5"Wx9.5"D)  L(6"Wx11.5"D)

(\*For use with Adjustable 70°-90° Fixed Center Mount Foot platform or CMX Center Mount Power Articulating & Elevating Foot Platforms)  
 ^ May cause interference with certain bases and seat widths call to confirm configuration.

- Adjustable Angle Footplates (HCPCS code K0040)  
 Left:  S(4.25"W x 8"D)  M(5.25"Wx8"D)  
 L(6.25"Wx8"D)  XL(9"Wx11"D) . . . . . \$125  
 Right:  S(4.25"W x 8"D)  M(5.25"Wx8"D)  
 L(6.25"Wx8"D)  XL(9"Wx11"D) . . . . . \$125

**FOOT PLATFORM & FOOTPLATE OPTIONS**

- Adjustable Angle Footplates with Cutouts for Caster Clearance (HCPCS code K0040)  
 Left:  S(7.25"W x 9"D)  M(7.5"W x 11"D)  L(10"W x 14"D) . \$125  
 Right:  S(7.25"W x 9"D)  M(7.5"W x 11"D)  L(10"W x 14"D) . \$125
- Heel Loops (Each) QUANTITY ( ) . . . . . \$20  
 (HCPCS code E0951)
- NEW**  Multi-Axis Adjustable Angle Footplates (HCPCS code K0040)  
 Left:  S(4.25"W x 8"D)  M(5.25"Wx8"D)  
 L(6.25"Wx8"D)  XL(9"Wx11"D) . . . . . \$225  
 Right:  S(4.25"W x 8"D)  M(5.25"Wx8"D)  
 L(6.25"Wx8"D)  XL(9"Wx11"D) . . . . . \$225
- Multi-Axis Adjustable Angle Footplates Cutout for Caster Clearance (HCPCS code K0040 )  
 Left:  S(7.25"W x 9"D)  M(7.5"W x 11"D)  L(10"W x 14"D) . \$225  
 Right:  S(7.25"W x 9"D)  M(7.5"W x 11"D)  L(10"W x 14"D) . \$225

**CALF PAD AND STRAP OPTIONS**

- Flat Calf pads . . . . . \$290  
 Select:  TRx HD Swing Away Footrests  
 Fixed Center Mount Foot platform
- Curved Calf pads (inside width is 6") . . . . . \$290  
 Select:  TRx HD Swing Away Footrests  
 TRx Pivot Plus Legrests  
 Fixed Center Mount Foot platform
- One piece calf pad for Fixed Center Mount Foot platform . . \$290
- Airmesh Calf Strap for chairs 15"-18" wide - 3.5" pad height (HCPCS code E0978) BP-CS-M . . . . . \$70
- Airmesh Calf Strap for chairs 18"-22" wide-3.5" pad height (HCPCS code E0978) BP-CS-L . . . . . \$70
- Airmesh Calf Panel for chairs 15"-18" wide 35.5"w x 9"h pad BP-CP-MD . . . . . \$155
- Airmesh Calf Panel for chairs 18"-22" wide 38.5"w x 9"h pad BP-CP-LG . . . . . \$155

**BODY POINT LAP BELTS**

- Padded Lap Belt 2 point - BP-PLB2 (HCPCS code E0978) . . \$155
- Padded Lap Belt 4 point - BP-PLB4 (HCPCS code E0978) . . \$195

**ADDITIONAL OPTIONS AND ACCESSORIES**

- TRx Concepts Back Pack - Single . . . . . \$160
- TRx Concepts Back Pack - Double . . . . . \$255
- Transfer Handles (Height Adjustable) . . . . . \$455  
 Select:  4"-6.5"  
 9"-11.5"
- D-Loops (set of 2)\* . . . . . \$325  
 -D-Loops increase the seat to floor height by 3/8"

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**"RFQ" = Request For Quotation** - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

**ELECTRONICS SWITCH MOUNTING**

Select:  Left hand  Right hand  Provide Mounting Hardware  
 Do Not Mount

**SINGLE FUNCTION ELECTRONICS**

Single Function through Switch . . . . .STD

Select:  Single Push Button  Single Toggle

Enhanced Single Functions through Switches . . . . . \$500

Select:  Single Push Button  Dual Push Button  
 Low Amp Toggle

Single Function Control Module . . . . . \$1,200

(HCPCS code E2310)

Through Joystick

Through Specialty Control-Fill out Specialty control configuration sheet on page 8

Additional Switch:  Single Push Button  Dual Push Button . . . \$97

**MULTIPLE ELECTRONIC FUNCTIONS**

Multiple Functions through Switches/Toggle . . . . . \$1,550

Select:  Dual Push Button  Dual Toggle  
 Quad Push Button  4-way toggle

Add:  Reduced Drive with Elevating Seat . . . . . \$400

Add:  Function Display\* . . . . . \$250

Add:  Function Display/Step Scanner\* . . . . . \$400

\* The display is operated through left/right motion on 4-Way Toggle or buttons 3 and 4 on the Quad Push Button.

Multiple Function Control Module ^ . . . . . \$2,600

(HCPCS code E2311)

Through Joystick

Through Specialty Control-Fill out Specialty control configuration sheet on page 8

Includes: Reduced Drive with Elevating Seat . . . . .  Yes  No

Includes: Function Display . . . . .  Yes  No

Add:  Function Display/Step Scanner\* . . . . . \$150

Additional Switch:  Quad Push Button  Four-Way Toggle . . . \$97

\* The display is operated through left/right motion on 4-Way Toggle or buttons 3 and 4 on the Quad Push Button.

(^ Motion Concepts provides the communication interface for Invacare Electronics and includes a function display)

**ELECTRONICS OPTIONS**

**REQUIRE ENHANCED OR MULTIPLE ELECTRONIC FUNCTIONS**

TRx M16 Splitter . . . . . \$165

(Allows for multiple switch options with phono plug connection)

TRx Feather Touch Switch - Single function . . . . . \$130

Dime Size  Quarter Size

TRx Feather Touch Switch Multi function . . . . . \$180

Dual Dime Size  Quad Dime Size

Stealth Egg Switch  Green  Black . . . . . \$105

(For use with Single function & Enhanced Single Function electronics only)

TRx Power Shut Off Switch - Toggle Switch Only . . . . . \$230

TRx Lazarus Jump Starter . . . . . \$195

Accessory Port Power Supply . . . . . \$265

(3 amp capacity)

TRx Lighting Effects - Through Switch . . . . . \$500

TRx Front Safety Lights (pair)

Operated through switches . . . . . \$650

Operated through Joystick . . . . . \$1,150

TRx Rear Safety Lights (pair)

Operated through switches . . . . . \$650

Operated through Joystick . . . . . \$1,150

TRX Front & Rear Light Combo Package Credit\* . . . (\$150)

\*Purchase both sets and receive a credit

**MEMORY SEATING OPTIONS**

Single Function Memory Seating\* . . . .Suggested Retail \$645

. . . . .Call for RFQ

(For use with Enhanced electronics only - Not available to control Elevating seat modules)

\*Please call for specification sheet 1-888-433-6818

Multiple Function Memory Seating\* .Suggested Retail \$1,325

. . . . .Call for RFQ

(For use with Enhanced electronics only - Not available to control Elevating seat modules)

\*Please call for specification sheet 1-888-433-6818

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_


PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834


## SPECIALTY CONTROL CONFIG SHEET

### Mult-Channel Joystick Control Configuration Sheet



Notes: \_\_\_\_\_

### 3-Axis Head Array Configuration Sheet



Notes: \_\_\_\_\_

### SIP and PUFF Configuration Sheet

ECU 1 Action:	Equivalent Action:	DB9 Pin out:	Function:
Hard Puff _____	Forward _____	Pin1 _____	
Hard Sip _____	Reverse _____	Pin2 _____	
Soft Sip _____	Left _____	Pin3 _____	
Soft Puff _____	Right _____	Pin4 _____	

ECU 2 Action:	Equivalent Action:	DB9 Pin out:	Function:
Hard Puff _____	Forward _____	Pin1 _____	
Hard Sip _____	Reverse _____	Pin2 _____	
Soft Sip _____	Left _____	Pin3 _____	
Soft Puff _____	Right _____	Pin4 _____	

Notes: \_\_\_\_\_

**"RFQ" = Request For Quotation** - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**O2 HOLDERS - ADULT SEAT SIZES ONLY**

- Non-Articulating O<sup>2</sup> Gas Holder (Cylinder Style) . . \$1,175
  - Articulating O<sup>2</sup> Holder . . . . . \$1,650
- Please Select One:
- O<sup>2</sup> Gas (Cylinder Style)       Liquid O<sup>2</sup> (6.5" x 4.5")

**BATTERIES**

- MK U1 Battery (MU-1 SLD G) Each . . . . . \$185  
Quantity ( ) (HCPCS code E2365)
- MK GROUP 22NF Battery (M22NF SLD G) Each  
Quantity ( ) (HCPCS code E2361) . . . . . \$230
- MK GROUP 24 Battery (M24 SLD G) Each  
Quantity ( ) (HCPCS code E2363) . . . . . \$285
- MK GROUP 34 Battery (M34 SLD G) Each  
Quantity ( ) . . . . . \$285

**INSTALLATION**

- Installation by Motion Concepts<sup>(23)</sup> . . . . . N/C
- Not Required<sup>(24)</sup> . . . . . (\$100)

**Total Retail Value \$** \_\_\_\_\_

**NOTES**

- 1 Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18". Suspension castors may affect seat-to-floor height up to an additional .5"-1".
- 2 Weight limit is 225lbs.
- 3 Uses 2 U1 Batteries
- 4 Systems are designed for use with a 2"-3" cushion.
- 5 Back angle is adjustable from 76° to 116° degrees in increments.
- 6 These systems are not available over 20" wide by 20" deep and increase seat to floor height. Elevating systems use 2-22NF batteries. Group 24 batteries are available for the TDX and TDX-SP upon request.
- 7 TDX-SP base must have the TRO (Transport ready option) if utilizing a vent and elevating seat module.
- 8 RFQ required for this seat size on M51 and M71.
- 9 RFQ required for this seat size on M51, M71, 1103, 1107, 1113 ATS.
- 10 We highly recommend a 3G - Long base for optimal maneuverability
- 11 3G - Long base must be used for these seat depths.
- 12 May be used as an interface for after market backs.
- 13 This option allows the back angle to adjust from 81° to 121° in increments.
- 14 When omitting the back there may be compatibility issues when using certain aftermarket back and/or lateral options. When using these arms with endomorph or mesomorph canes the lowest arm setting is 13.5".
- 15 Laterals are not available with the MaTRx PB Backs.
- 16 Standard front rigging will only be provided when specified. Please provide or order separately.
- 17 Includes a standard flip up foot platform. May increase seat-to-floor height.
- 18 90° front rigging will collide with front casters on the P220, S525, S646 and 3G Storm.
- 19 The increased weight of this option reduces the weight capacity by 20 lbs on the M71, 1107, 1103, 1113 ATS.
- 20 Under normal usage Power Legrest actuators tend to be louder than Tilt/Recline Actuators.
- 21 Individual legs count as two functions.
- 22 Heavy duty durable construction (4.5lbs ea.).
- 23 Batteries may be required at time of installation, please call to confirm.
- 24 Elevating seat systems must be installed by Motion Concepts.