



# Low-PRO

## TRx-Pediatric CG System for Quickie Z500

### US PRICE LIST AND ORDER FORM


Effective February 2005



To ensure system is accurately configured please fill in all required measurements.

Date of Order: \_\_\_\_\_  
 Dealer Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_  
 Tag: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 Client Info: Gender: \_\_\_\_\_ Special Client Conditions: \_\_\_\_\_

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Seat Cushion Thickness:	



\* 125 lbs. maximum weight limit. Dependent on weight capacity of the power base.

#### BASE TYPE - MUST BE ORDERED "TRx READY"

Z500 Seat-to-floor height is 18.75"

Please Specify:

- Left Hand Drive     Right Hand Drive

Please Specify Electronics:

- Qtronix     VSI

#### SYSTEMS<sup>(2)</sup> MUST SELECT ONE SYSTEM

TRx Low-PRO CG Tilt Only System (55° Tilt)<sup>(3)</sup>  
 (Includes seat pan, back pan, armrests, armpads, fixed front rigging, electronics and standard switch controls )

(HCPCS code E1002) ..... \$5,495

TRx Low-PRO CG Tilt (55° Tilt) and Power Recline (174°) System  
 (Includes seat pan, back pan, armrests, armpads, fixed front rigging, electronics and standard switch controls )

(HCPCS code E1006) ..... \$8,490

TRx Low-PRO CG Tilt (55° Tilt) and Power Recline (174°) System with Extended Shear Reduction (ESR)  
 (Includes seat pan, back pan, armrests, armpads, fixed front rigging, electronics and standard switch controls )

(HCPCS code E1007) ..... \$9,055

TRx Elevate Only System  
 (Includes seat pan, back pan, armrests, armpads, fixed front rigging, electronics and standard switch controls ) (HCPCS code E2300)

- 8" Elevation     7" Elevation

Seat to Floor on Low-PRO CG Tilt with Elevating Module<sup>(1)</sup> 19.5"

(HCPCS code E1002) ..... \$5,495

#### MODULES MUST SELECT TRx SYSTEM-TO ADD MODULES

TRx Precline Module (up to 30° Precline)<sup>(4)</sup> ..... \$500  
 (HCPCS code K0108)

TRx Precline Module (up to 8° Precline)<sup>(4)</sup> ..... N/C  
 (HCPCS code K0108)

TRx Manual Recline (135°) ..... \$2,000  
 (HCPCS code K0108)

TRx Elevating Seat Module ..... \$2,995  
 (HCPCS code E2300)

- 7" Elevation     6" Elevation

Seat to Floor on Low-PRO CG Tilt with Elevating Module<sup>(1)</sup> 19.5"

#### SEAT SPECIFICATIONS

Pediatric Seat Width  
 11"     12"     13"     14"     15" ..... N/C

Pediatric Seat Depth  
 11"     12"     13"     14"     15" ..... N/C

#### SEAT CUSHION OPTION

TRx Kid-ab-ra Matrix-V Seat Cushion 11"w-14"w ..... \$425  
 (HCPCS Code K0656)

TRx Matrix-V Seat Cushion 15"w ..... \$425

Please note that the HCPCS codes indicated on this order form are our interpretation of the new codes and how we believe our products should be coded. The individual DMERCs will determine if our interpretation is appropriate and whether or not items coded will be accepted or denied. Motion Concepts cannot guarantee that use of the designated codes will ensure coverage or payment. Contact the insurer to confirm compliance with coding and coverage policy.

#### NOTES

- |   |  |  |
|---|--|--|
| 1 Seat to floor height is measured at the front of the seat pan based on a standard seat size with a 0 degree pre-tilt. The Z500 uses 2-22NF or 2-U1 Batteries. | 3 Back angle is adjustable from -5° to 15° degrees in 5 degree increments. | 4 Total degrees of recline is 174° ie: if 30° degrees of precline is used, then 144° degrees of recline will be available. |
| 2 Systems are designed for use with a 2"-3" cushion.  |  |  |

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**BACK OPTIONS**

- Curved Rigid Back and Back cushion with 1" foam . . . . .STD
- Recessed planar interface plate<sup>(5)</sup> . . . . . N/C

**BACK SPECIFICATIONS**

**Back Height - Seat pan to top of back shell**

- 14"    15"    16"    17"    18"    19"
- 20"    21"    22"    23"    24"    25"
- Custom Back Height Upcharge . . . . . \$450  
Specify height \_\_\_\_\_  
(Comes with standard cane height unless specified by RFQ)

**HEADRESTS & MOUNTING HARDWARE**

- Headrests with Fixed Mounting Hardware (HCPCS code E0955)**
- Concepts Standard Headrest - with Concepts fixed hardware . . \$250
- Combination Head/Neck - Small - with Concepts fixed hardware \$295
- Combination Head/Neck - Large - with Concepts fixed hardware \$295
- Stealth Comfort Plus 10" (CP250) - with Concepts fixed hardware \$242
- Stealth Comfort Plus 14" (CP150) - with Concepts fixed hardware \$242
- Headrest Swing Aways and Mounts - Does not include the Headrest (HCPCS code E1028)**
- Concepts Detachable Mounting Hardware . . . . . \$200
- Stealth's World's Best Headrest Hardware (TWB485) . . . . . \$220

**ARMRESTS**

- Pediatric Adj. Height Select:  Dual Post  Reclining . . . . . STD
- S (6"-7.5")    M (6.75"-8.25")    L (8.25"-9")

**ARMPADS**

- Standard Pads-   Full (2.25"x14")    Left    Right . . . . . N/C
- Desk (2.25"x10")    Left    Right . . . . . N/C
- Flat Pads- (Pair) . . . . . \$35
- Select:    Full (3"x14")    Left    Right
- Desk (3"x10")    Left    Right

**STANDARD LEGRESTS**

CALL TO CONFIRM ON SEAT SIZES BELOW 16"W X 16"D

- Legrest Hangers Only . . . . . STD
- \*\*Specify Leg/Footrest style to be used on TRx System <sup>(1)</sup>**

OR \_\_\_\_\_

- 70° Fixed Footplate Mount . . . . . N/C
- Select Extension:  S (4"-6")    M (7"-10")    L (11"-14")

**ELEVATING LEGRESTS**

- Quickie Standard Elevating Legrests (Pair) . . . . . \$495
- Quickie P210 Elevating Legrests (Pair) . . . . . \$495
- Action AT5544 Elevating Legrests (Pair) . . . . . \$495

**LEGREST OPTIONS**

- Quickie 70° Swingaway Footrests (Pair) . . . . . \$395
- Action 70° Swingaway Footrests (Pair) . . . . . \$395
- Flip-up Foot Platform (for use with above fixed footplate mount) . . . . . \$180
- Adjustable Angle Footplates (Each)(HCPCS code K0040)
- Left  S(4.25" x 8")  M(5.25"x8")  L(6.25"x8")    \$125
- Right  S(4.25" x 8")  M(5.25"x8")  L(6.25"x8")    \$125
- Heel Loops (Each) **QUANTITY (    )** . . . . . \$20  
(HCPCS code E0951)

**NOTES**

5 Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.

\*\* For Custom variations - please call for quotation

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**LATERAL OPTIONS**

LEFT Lateral Support with Fixed Mounting  
 Pad Size: S(3.25"x5.25") M(4"x6.25") L(4.5"x7.5")  
 (HCPCS code E0956) .....\$100

RIGHT Lateral Support with Fixed Mounting  
 Pad Size: S(3.25"x5.25") M(4"x6.25") L(4.5"x7.5")  
 (HCPCS code E0956) .....\$100

TRx Swingaway Mounting Hardware  
 (to convert above to Swingaway Laterals)  
 LEFT (HCPCS code E1028) .....\$195  
 RIGHT (HCPCS code E1028) .....\$195

**Select Offset for Fixed or Swingaway Lateral Support:**

Right: 0"-3" 3-6"

Left: 0"-3" 3-6"

**Select Lateral Latch plate length for Fixed or Swingaway Lateral Supports:**

Right: Short (3.5") Long(5")

Left: Short (3.5") Long(5")

LEFT Integrated Arm Block\* .....\$125  
 Pad Size: S(3.25"x5.25") M(4"x6.25") L(4.5"x7.5")

RIGHT Integrated Arm Block\* .....\$125  
 Pad Size: S(3.25"x5.25") M(4"x6.25") L(4.5"x7.5")

(\*Prevents arm from sliding back off armrest while tilted. Not available with Recline.)

**ELECTRONICS SWITCH MOUNTING**

Select:  Left Side Mounting  Right Side Mounting

**STANDARD ELECTRONICS**

(These electronics are equipped with drive lock-out. They do not include attendant control, tilt limits, and variable speed control(s).)

TRx Single-Function Control Box .....STD  
 Please Select One:

Single Push Button  Dual Push Button .....  
 Single Toggle  Omit Switch\*

(\*Single-Function Control Box uses a phono jack connection. Compatible with other phono jack switches.)

Through Qtronix Joystick (HCPCS code E2310) .....\$750

TRx Dual-Function Control Box - Tilt/Recline Only\* .....STD  
 (\*ESR or PSB included - PSB is not operated as separate function - order 3-5 function Enhanced electronics for separate control)

Please Select One:  
 Dual Push Button  Quad Push Button<sup>(6)</sup>  
 Dual Toggle  Four-Way Toggle

Through Qtronix Joystick (HCPCS code E2311) .....\$950

**HIGH END ELECTRONICS**

(These electronics are equipped with drive lock-out, attendant control, tilt limits and variable speed control(s).)

**Single-Function System - (ie: Tilt only)**

TRx Enhanced Single-Function Control Box .....\$500

Please Select One:

Single Push Button  Dual Push Button<sup>(6)</sup>  
 Single Toggle

Enhanced Single Function through Qtronix Joystick .....\$1,150  
 (HCPCS code E2310)

**Dual-Function System - (ie: Tilt/Recline)**

TRx Enhanced Dual-Function Control Box .....\$700

Please Select One:

Dual Push Button  Quad Push Button<sup>(6)</sup>  
 Dual Toggle  Four-Way Toggle

Enhanced Dual Function through Qtronix Joystick .....\$1,550  
 (HCPCS code E2311)

**Multi-Function System - (ie: Tilt/Recline/Future options)**

TRx Enhanced 3-5 Function Control Box .....\$1,150

Please Select One:

Quad Push Button  Four-Way Toggle  
 Eight-Way Toggle

Enhanced 3-5 Function through Qtronix Joystick .....\$2,350  
 (HCPCS code E2311)

**ELECTRONICS OPTIONS**

Stealth Green Egg (ESGRN) .....\$77

(For use with Single function & Enhanced Single Function electronics only)

TRx Lazarus Jump Starter .....\$150  
 (HCPCS code A9270)

TRx M16 Splitter<sup>(7)</sup> .....\$125  
 (Allows for multiple switch options) (HCPCS code E2399)

TRx M06 Attendant Power Switch<sup>(7)</sup> .....\$175  
 (HCPCS code E2399)

**OMISSIONS**

Rigid Back Pan<sup>(8)</sup> .....(\$75)

Arms .....(\$75)

Arm Pads .....(\$25)

"ON-CHAIR" Matrix-V Cushion Credit .....(\$90)

**INSTALLATION**

Installation by Motion Concepts<sup>(9)</sup> .....N/C

Not Required .....(\$100)

MK GROUP 22NF Battery (M22NF SLD G) Each  
 Quantity ( ) (HCPCS code E2361) .....\$163

MK U1 Battery (MU-1 SLD G) Each  
 Quantity ( ) (HCPCS code E2365) .....\$175

**Total Retail Value \$** \_\_\_\_\_

**NOTES**

6 Non-reversing relays.  
 7 Must order an Enhanced control box.

8 N/A with Recline/ESR .

9 Batteries may be required at time of installation, please call to confirm.

\*\* For Custom variations - please call for quotation