

Easy Does it

TRx-CG Power Positioning System for Quantum Q600

US PRICE LIST AND ORDER FORM

Effective October 2005




Date of Order: _____
 Dealer Name: _____
 Dealer Account #: _____
 PO #: _____
 Tag: _____
 Purchasing Contact: _____
 Phone: _____ Fax: _____
 RTS/Therapist: _____
 Ship to Address: _____

 City: _____ State: _____ Zip: _____

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Seat Cushion Thickness:	



300lbs. maximum weight limit. Dependent on weight capacity of the power base.

BASE SELECTION

- Quantum 600 2SP (Group 2 Single Power)
(HCPCS code K0835) \$5,595
- Quantum 600 3SP (Group 3 Single Power)
(HCPCS code K0835) \$5,595
- Quantum 600 2SP (Group 3 Multiple Power)
(HCPCS code K0835) \$5,595
Programmable speed up to 5 mph

COLOR

- Champagne Candy Apple Red Viper Blue
- Onyx Black Sunburst Orange

DRIVE WHEELS

- 14" Pneumatic (pair) STD
- 14" Flat-Free inserts(pair) \$213
(HCPCS code K0093)

CONTROLLER

P&G Electronics

- VSI Integral 70 amp⁽¹⁾ N/C
- VSI Integral 6-Key 70 amp⁽²⁾ (HCPCS code E2399) \$200
- 3-Key Remote Plus 70 amp⁽³⁾⁽⁴⁾ (HCPCS code E2320) \$1,025
Select one 1 Profile Multiple Profile
- 7-Key Access Remote Plus 70 amp (HCPCS code E2320) \$1,300
- Remote Plus 70 amp with Speed Dial⁽¹⁾⁽⁴⁾ (HCPCS code E2320) . \$1,300
- Remote Plus 70 amp with Toggles⁽³⁾⁽⁵⁾ HCPCS code E2320) \$1,300
- Remote Plus 70 amp with Switch Jacks⁽³⁾⁽⁵⁾ (HCPCS code E2320) \$1,300
- Microdrive, 70 amp with Switch Jacks⁽³⁾⁽⁶⁾ (HCPCS code E2320) . \$1,300

JOYSTICK HANDLE OPTIONS

- P&G Controls (HCPCS code E2323)
- Carrot Handle \$20 Mushroom Handle . . \$71
 - T-Bar \$45 Chin Cup \$71
 - Extended I-Handle . . . \$71 Soft Ball Handle . . . \$50
 - Medium Ball Handle . \$50

PROGRAMMER

- Remote Plus/Pilot Programmer \$423

CONTROLLER MOUNT

- Fixed, In-line, Joystick Mount Select: Left Right STD
- Quick Disconnect, Joystick Mount Select: Left Right . . N/C
- Swing-Away Joystick Mount Select: Left Right \$196
(HCPCS code E1028)
- Multi-axis Joystick Mount Select: Left Right \$100
(HCPCS code E1028)

BATTERY CHARGER

- 8 amp Off-Board Charger \$360

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SELECT SEAT-TO-FLOOR HEIGHT

Select Seat-to-Floor (where applicable) ⁽⁷⁾:

Q600

Select STF⁽¹⁾: 17.5"* 18.5" 19.5" 20.5"

*This STF not available with D-Loops

(If a specific seat to floor height is required please request STF confirmation)

BATTERIES - REQUIRED⁽²⁹⁾

MK GROUP 22NF Battery (M22NF SLD G) Each
Quantity () (HCPCS code E2361)\$163

SYSTEMS⁽⁸⁾ MUST SELECT ONE SYSTEM

(Includes seat pan, back pan, armrests, arm pads, fixed front rigging, electronics and standard switch controls)

TRx Ultra Low CG Tilt Only System (55° Tilt)⁽⁹⁾
(HCPCS code E1002)\$5,120

TRx Ultra Low CG Tilt (55° Tilt) and Power Recline (174°) System
(HCPCS code E1006)\$8,115

TRx Ultra Low CG Tilt (55° Tilt) and Power Recline (174°) System with Extended Shear Reduction (ESR)
(HCPCS code E1007)\$8,680

TRx Ultra Low CG Tilt(55° Tilt) and Power Recline (174°) System with Synchronized Power Sliding Back (PSB)
(HCPCS code E1008)\$9,115

TRx Recline Only System (174°)
(HCPCS code E1003)\$5,695

TRx Recline Only System (174°) with Extended Shear Reduction (ESR)
(HCPCS code E1004)\$6,260

TRx Recline Only System (174°) with Synchronized Power Sliding Back (PSB)
(HCPCS code E1005)\$7,460

TRx Elevate Only System⁽¹⁰⁾
(HCPCS code E2300)\$4,495

You must select TRx Enhanced 3-5 Function electronics to drive while elevated.

PSB is also available as an independent function - Please select 3-5 function Enhanced electronics.

MODULES MUST SELECT TRx SYSTEM-TO ADD MODULES

TRx Recline Body Fit System (Fits 16" - 20" widths) . \$995
(This option allows the seat back to be narrower than the seat pan.
For example an 18" back pan width with a 20" seat width is accomplished by selecting "2" Narrower" option.)

2" Narrower 4" Narrower

TRx Precline Module (up to 30° Precline)⁽¹¹⁾\$500
(HCPCS code K0108)

TRx Precline Module (up to 8° Precline)⁽¹¹⁾N/C
(HCPCS code K0108)

TRx Manual Recline (135°)\$2,000
(HCPCS code K0108)

TRx Elevating Seat Module⁽¹⁰⁾\$2,250
(HCPCS code E2300)

7" Elevation - 225 lbs max. 6" Elevation - 225 lbs max.

Seat-to-floor height is 19.313"

You must select TRx Enhanced 3-5 Function electronics to drive while elevated.

STANDARD SEAT SPECIFICATIONS

Seat Width

16" 17" 18" 19" 20"N/C

Seat Depth

15" 16" 17" 18" 19" 20"N/C

CUSTOM SEAT SPECIFICATIONS

Pediatric Seat Width (Only available with Pediatric Seat Depths)

11" 12" 13" 14" 15"\$250

Pediatric Seat Depth (Only available with Pediatric Seat Widths)

11" 12" 13" 14"\$250

BACK SPECIFICATIONS

Back Height - Seat pan to top of back shell

18"⁽¹²⁾ 19"⁽¹²⁾ 20" 21"STD

22" 23" 24" 25"STD

Custom Back Height Upcharge Specify height _____ . . . \$450
(Comes with standard cane height unless specified)

TRx Recline Body Fit System (Minimum back width of 16") . \$995
(This option allows the seat back to be narrower than the seat pan.
For example an 18" back pan width with a 20" seat width is accomplished by selecting "2" Narrower" option.)

2" Narrower 4" Narrower

BACK CANE SPECIFICATIONS - TILT/ELEVATE ONLY

Select Back Cane Height - Measured from Seat pan to top of back cane

22" 24"N/C

Custom Cane Height(Tilt Only)Specify height _____ . . . \$250

Angled (5°)Tilt Back Canes⁽¹³⁾N/C

Endomorph Back Canes Select: 1" 2" 4" NarrowerN/C

Mesomorph Back Canes* Select: 1" 2" 4" Wider \$250

* Comes with Single Post Flip Back Arms

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

BACK OPTIONS

- Standard Rehab Back and Back cushionSTD
- Recessed planar interface plate⁽¹⁴⁾N/C
- TRx Single Post Back interface plate* (Tilt Only)N/C
(*Single Post Back is angle adjustable in 6° increments from 24° forward to 18° back. (The back folds down and is only available for 16"-20" seat widths & 21"-23" back height)
- OMIT: Rigid Back Pan & Cushion⁽¹⁵⁾(\$75)
- OMIT: Back Cushion Only(\$25)

UPGRADEABLE BACK OPTIONS

- MaTRx Contour Back (HCPCS Code E2613)\$495
(For Systems 15"-21" wide under 300lbs)
 Specify actual back height dimension: 12" 15" 17"
- MaTRx Contour Back (HCPCS Code E2614)\$595
(For Systems 22" wide under 300lbs)
 Specify actual back height dimension: 15" 17"
- MaTRx High Back (HCPCS Code E2620) \$595
(For Systems 16"-21" wide under 300lbs)
 Specify actual back height dimension: 20" 25" 30"
- MaTRx PB Elite Back(HCPCS Code E2620) \$625
(For Systems 12"-20" wide under 300lbs)
 Specify actual back height dimension: 12" 16" 20"
- MaTRx PB Elite Back with Laterals (HCPCS Code E2620) . . \$775
(For Systems 12"-20" wide under 300lbs)
 Specify actual back height dimension: 12" 16" 20"
- MaTRx PB Back (5" Contour) (HCPCS Code E2615) \$625
(For Systems 14"-20" wide under 300lbs)
 Specify actual back height dimension: 12" 16" 20"
- MaTRx PB Deep Back (7" Contour) (HCPCS Code E2615) . \$695
(For Systems 14"-20" wide under 300lbs)
 Specify actual back height dimension: 12" 16" 20"
- MaTRx PB Heavy Duty Back (5" Contour) \$835
(For Systems 20"-22" wide and over 300lbs)
 (HCPCS Code E2615 for 20" width & HCPCS Code E2616 for 22" width)
 Specify actual back height dimension: 16" 20"
- MaTRx PB Heavy Duty Deep Back (7" Contour) \$935
(For Systems 20"-22" wide and over 300lbs)
 (HCPCS Code E2615 for 20" width & HCPCS Code E2616 for 22" width)
 Specify actual back pan dimensions: 16" 20"
- "ON-CHAIR" Matrix Back Credit(\$125)

NOTE WHEN CHOOSING MaTRx Backs

4" Gap required for Tilt
 6" Gap required for systems with Recline/ESR
 * PB Backs are not available with the Precline Module.

SEAT CUSHION OPTIONS

- Kid-ab-ra MaTRx-V Seat Cushion 11"w-14"w \$425
- MaTRx-V Seat Cushion 14"w-20"w (HCPCS Code E2607) . \$425
- MaTRx-V Seat Cushion 21"w-22"w (HCPCS Code E2608) . \$475
- MaTRx-Flo-Tech Cushion 15"w-20"w (HCPCS Code K0734) \$475
- "ON-CHAIR" Matrix-V Cushion Credit (\$90)

CUSHION AND BACK FABRIC SELECTION

- Onyx Fabric Upcharge - Select below
- Contour Back...\$100 High Back...\$100 MaTRx Cushion...\$100

ARMRESTS

- Dual Post Adj. HeightSTD
 Select Size: (9.75"-12.75") (12.75"-15.75") (15.75"-18.75")
- Dual Post with Quick Height Adjustment Lever \$250
 Select Size: (9.75"-12.75") (12.75"-15.75") (15.75"-18.75")
- Reclining Adj. HeightSTD
 Select Size: (9.5"-13") (13.5"-18.5")
- Pediatric Adj.Height Select: Dual Post Reclining . . \$250
 S (6"-7.5") M (6.75"-8.25") L (8.25"-9")
- TRx Height & Angle Adjustable Cantilever
 Flip Back Armrests⁽¹⁰⁾ Tilt/Elevate Recline
(*Not Available with a Single Post Back) \$250
- 1" Offset Armrest Spacer (each) \$125
- 1/2" Offset Armrest Spacer (each) \$125
- TRx Armrest Storage Compartment RIGHT*N/C
- TRx Armrest Storage Compartment LEFT*N/C
(*N/A with Recline/ Pediatric or Single Pivot Flip Back Armrests or Mesomorph Back Cones)
- OMIT: Armrests (\$75)

ARMPADS

- Standard Pads-STD
 Select: Full (2.25"x14") Left Right
 Desk (2.25"x10") Left Right
- Waterfall Arm Pads-(Pair) \$35
 Select: Full (2.25"x14") Left Right
 Desk (2.25"x10") Left Right
- Flat Pads- (Pair) \$35
 Select: Full (3"x14") Left Right
 Desk (3"x10") Left Right
- Gel Pads- (Pair) \$150
 Select: Narrow (2"x12") Left Right
 Wide (3.5"x12") Left Right
- TRx Contour Arm Trough - LEFT (4.5"x13") \$85
(HCPCS Code K0106)
- TRx Contour Arm Trough - RIGHT (4.5"x13") \$85
(HCPCS Code K0106)
- TRx Flat Hand Pad - LEFT (4.5"x5") \$55
(HCPCS Code K0108)
- TRx Flat Hand Pad - RIGHT (4.5"x5") \$55
(HCPCS Code K0108)
- TRx Arm Pad Storage/Cup Holder LEFT * \$175
- TRx Arm Pad Storage/Cup Holder RIGHT * \$175
*When using a joystick this feature is only available on the armrest opposite the joystick.
- TRx Elbow Block LEFT* \$195
 Pad Size: XS(3"x4") S(3.5"x5") M(3.75"x5.75") L(4.5"x6.25")
- TRx Elbow Block RIGHT* \$195
 Pad Size: XS(3"x4") S(3.5"x5") M(3.75"x5.75") L(4.5"x6.25")
(*Prevents arm from sliding back off armrest while tilted. Not recommended with Recline)
- Multi Purpose Arm Pad Adapter Block (each) \$35
(Allows hard mounting & positioning of standard pads or other manufacturer's arm pads)
- OMIT: Armpads (\$25)

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HEADRESTS & MOUNTING HARDWARE

- Headrests with Fixed Mounting Hardware** (HCPCS code E0955)
- Concepts Standard Headrest - with Concepts fixed hardware . . .\$295
 - Combination Head/Neck - Small - with Concepts fixed hardware . . .\$295
 - Combination Head/Neck - Large - with Concepts fixed hardware . . .\$295
 - Stealth Comfort Plus 10" (CP250) - with Concepts fixed hardware . . . \$242
 - Stealth Comfort Plus 14" (CP150) - with Concepts fixed hardware . . . \$242
 - Onyx Auto Style Headrest - Fixed Mounting only available \$350
- Headrest Swing Aways and Mounts - Does not include the Headrest** (HCPCS code E1028)
- Concepts Detachable Mounting Hardware \$200
 - Stealth's World's Best Headrest Hardware (TWB485) \$220
 - OMIT: Headrest** N/C

LATERAL OPTIONS

- TRx Lateral Trunk Support with Fixed Mounting LEFT⁽¹⁴⁾**
 Pad Size: S(3.5"x5") M(3.75"x5.75") L(4.5"x6.25")
 (HCPCS code E0956) \$100
 - TRx Lateral Trunk Support with Fixed Mounting RIGHT⁽¹⁴⁾**
 Pad Size: S(3.5"x5") M(3.75"x5.75") L(4.5"x6.25")
 (HCPCS code E0956) \$100
 - TRx Swingaway Mounting Hardware⁽¹⁴⁾**
 (to convert above to Swingaway Laterals)
 - LEFT (HCPCS code E1028) \$195
 - RIGHT (HCPCS code E1028) \$195
- Select Offset for Fixed or Swingaway Lateral Trunk Support:
- Right: 0"-3" 3-6"
 Left: 0"-3" 3-6"
- Select Lateral Latch plate length:
- Right: Short (3.5") Long(5")
 Left: Short (3.5") Long(5")

STEALTH LATERAL OPTIONS⁽¹⁷⁾

- TWB-L** Pad Size: S(3" x 5") L(4"x6") \$590
- Note:** Due to the thickness of the Stealth Lateral pads the armrests will not flip back. Should you require this function, please order the 1" offset armrest spacers.

STANDARD FRONT RIGGING

CALL TO CONFIRM ON SEAT SIZES BELOW 16"W X 16"D

- Legrest Hangers Only STD
 - **Please Specify Front rigging to be used on TRx System⁽¹⁸⁾**
 - Fixed Center Mount Footplate Mount 70° 90° N/C
 (Comes with Footplate)
 Select Extension:
 S (12"-16") M (15"-19") L (18"-22")
 - Action 70° Swingaway Footrests (Pair)* N/C
- *Not covered by Motion Concepts Lifetime warranty

POWER FRONT RIGGING OPTIONS

- Center Mount Power Articulating & Elevating Foot Platform⁽¹⁹⁾ . . . \$2,995
 Extensions: (12"-16") (16"-21") (HCPCS code K0108)
 (Available on 16"-22" widths 18"-20" depths)
-
- Independent Elevating & Articulating Legrests
 (Center Mounted)⁽¹⁹⁾⁽²⁰⁾ \$2,995
 -can be mounted in multiple positions across front of seat frame
 Extensions: (12"-16") (16"-21") (HCPCS code K0108)
 Select Footplate: S(4.25" x 8") M(5.25"x8") L(6.25"x8")
 (Available on 16"-22" widths 16"-22" depths)
- TRx Power Elevating Pivot Plus Legrests⁽²⁰⁾ \$2,490
 (Includes Pivot Plus legrest & Actuator)
 (HCPCS code E1010)
 Select Function: Individual⁽²¹⁾ Combined Synchronized
 Extensions: S (12"-16") M (15"-19") L (18"-21")
 Select Footplate: Regular (7.5" x 9") Large (9 x 11")
- TRx Power Elevating Jazzy Legrests⁽²⁰⁾ \$2,490
 (Includes legrest & Actuator) (HCPCS code E1010)
 Select Function: Individual⁽²¹⁾ Combined Synchronized
- TRx Power Legrest Actuators (PAIR)⁽¹¹⁾⁽¹²⁾ \$1,695
 (Legrests not included - For use with Elevating Legrests)
 Select Function: Individual⁽²¹⁾ Combined Synchronized

ELEVATING LEGRESTS

- TRx Pivot Plus Manual Elevating Legrests (Pair) \$720
 Extensions: S (12"-16") M (15"-19") L (18"-21")
 Select Footplate: Regular (7.5" x 9") Large (9 x 11")
- Jazzy Standard Elevating Legrests (Pair) \$495

ADDITIONAL FRONT RIGGING OPTIONS

- TRx-Heavy Duty 70° Swingaway Footrests (Pair)⁽²¹⁾ \$420
 Select Receivers: Straight Flared(1")
 - Reversible inward/outward
 Select Extensions: S (12"-16") M (15"-19") L (18"-22")
 Select Footplate: S(4.25" x 8") M(5.25"x8") L(6.25"x8")
- Action 70° Swingaway Footrests (Pair) \$395
- Adjustable 70°-90° Center Mount Footplatform⁽¹⁴⁾ \$595
 (Depth and Height Adjustable - Standard Foot Platform and calf pads are included)
 S (12"-16") M (15"-19") L (18"-22")
- Foot cage fits 16"-20" widths \$895

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FOOTPLATES

Upgraded Flip-up Foot Platform (11" x 10") \$180
(for use with Adjustable 70°-90° Center Mount Footplatform)

Bariatric Flip-up Foot Platform(15.5" x 10.75") \$250
(for use with Adjustable 70°-90° Center Mount Footplatform)

Pediatric Flip-up Foot Platform(8" x 10") \$250
(for use with Adjustable 70°-90° Center Mount Footplatform)

Split Footplates (6" x 10" each) pair \$275
(For use with Adjustable 70°-90° Center Mount Footplatform or Center Mount Power Articulating & Elevating Foot Platform)

Adjustable Angle Footplates(Each)(HCPCS code K0040)
Left S(4.25" x 8") M(5.25"x8") L(6.25"x8") \$125
Right S(4.25" x 8") M(5.25"x8") L(6.25"x8") \$125

Adjustable Angle Footplates with Cutouts for
Caster Clearance (Each)(HCPCS code K0040)
Left Regular (7.5" x 9") Large (9 x 11") \$125
Right Regular (7.5" x 9") Large (9 x 11") \$125

Heel Loops (Each) QUANTITY () \$20
(HCPCS code E0951)

OTHER OPTIONS

Calf pads for the Fixed Center Mount Footplate \$220

Calf pads for TRx HD 70 Swing Away Legrests \$220

ELECTRONICS SWITCH MOUNTING

Select: Left Right Do Not Mount

Specialty Controls _____

Please call for specification sheet 1-888-433-6818

STANDARD ELECTRONICS - N/A WITH VENT TRAYS OR HEAD ARRAYS

(These electronics are equipped with drive lock-out. They do not include attendant control, tilt limits, and variable speed control(s).)

TRx Single-Function Control Box \$STD

Please Select One:

Single Push Button Dual Push Button

Single Toggle Omit Switch*

(*Single-Function Control Box uses a phono jack connection. Compatible with other phono jack switches.)

Through P&G Joystick \$750
(HCPCS code E2310)

STANDARD ELECTRONICS - N/A WITH VENT TRAYS OR HEAD ARRAYS

TRx Dual-Function Control Box - Tilt/Recline Only* \$STD
(*ESR or PSB included - PSB is not operated as separate function - order 3-5 function Enhanced electronics for separate control)

Please Select One:

Dual Push Button Dual Toggle

Dual Function Through P&G Joystick \$950
(HCPCS code E2311)

HIGH END ELECTRONICS

(These electronics are equipped with drive lock-out, attendant control, tilt limits and variable speed control(s).)

SINGLE-FUNCTION SYSTEM - (IE: TILT ONLY)

TRx Enhanced Single-Function Control Box \$500

Select: Single Push Button Dual Push Button⁽²²⁾

Single Toggle

Add: Reduced Drive with Elevating Seat \$400

Enhanced Single Function through Joystick \$1,150

(HCPCS code E2310)

Add: Reduced Drive with Elevating Seat \$400

DUAL-FUNCTION SYSTEM - (IE: TILT/RECLINE)

TRx Enhanced Dual-Function Control Box \$700

Select: Dual Push Button Quad Push Button⁽²²⁾

Dual Toggle Four-Way Toggle

Enhanced Dual Function through Joystick \$1,550

(HCPCS code E2311)

MULTI-FUNCTION SYSTEM - (IE: TILT/RECLINE/POWER LEGS)

TRx Enhanced 3-5 Function Control Box \$1,150

Select: Quad Push Button Four-Way Toggle

Add: Reduced Drive with Elevating Seat \$400

Add: Function Display* \$250

* The display is operated through left/right motion on 4-Way Toggle or buttons 3 and 4 on the Quad Push Button.

Enhanced 2-5 Function through Joystick \$2,350

(HCPCS code E2311)

Add: Reduced Drive with Elevating SeatN/C

ELECTRONICS OPTIONS

Stealth Green Egg (ESGRN) \$77
(For use with Single function & Enhanced Single Function electronics only)

TRx Lazarus Jump Starter \$150
(HCPCS code A9270)

TRx M16 Splitter⁽²⁵⁾(HCPCS code E2399) \$125
(Allows for multiple switch options with phono plug connection)

TRx M06 Attendant Power Switch⁽²⁵⁾ \$175
(HCPCS code E2399)

Accessory Port Power Supply \$195
(3 amp capacity)

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

VENTS AND O² HOLDERS (CALL TO CONFIRM OPTIONS ON SEAT SIZES BELOW 16"W X 16"D)

TILT WILL BE LIMITED TO 45° WHEN USING VENT APPLICATIONS

- TRx Vent Tray - Articulating⁽²⁶⁾ Horizontal Vertical
 - With Battery Tray⁽²⁷⁾\$1,750
 - Without Battery Tray \$1,350

- TRx Vent Tray -Gimbal (Pendulum Style)⁽²⁶⁾
 - With Battery Tray⁽²⁷⁾\$1,750
 - Without Battery Tray \$1,350

- TRx Vent Tray - Laptop Style⁽²⁶⁾
 - With Battery Tray⁽²⁷⁾\$1,750
 - Without Battery Tray \$1,350

- TRx Vent Bag Holder ⁽²⁶⁾
 - With Battery Tray⁽²⁷⁾\$1,540
 - Without Battery Tray \$990

- TRx Vent Bag Holder - Non Articulating⁽²⁶⁾
 - With Battery Tray⁽²⁷⁾\$1,075
 - Without Battery Tray \$525

- Non-Articulating O² Gas Holder (Cylinder Style)\$895

- Articulating O² Holder ⁽²⁸⁾\$1,250

Please Select One:

 - O² Gas (Cylinder Style) Liquid O² (6.5" x 4.5")

OPTIONS

- TRx Lighting Effects\$350
- TRx Headlights\$495
- TRx Concepts Back Pack - Single\$125
- TRx Concepts Back Pack - Double\$195
- Transfer Handles (Height Adjustable from 3.5"-6") ...\$350

CALL TO CONFIRM OPTIONS ON SEAT SIZES BELOW 16"W X 16"D AND FOR OTHER AVAILABLE OPTIONS

INSTALLATION

- Installation by Motion Concepts⁽²⁹⁾N/C
- D-Loops (set of 2)\$250

*D-Loops increase the seat to floor height by 3/8"

NOTES

1. A separate switch is required for power positioning functions. Functions cannot be run through this VSI controller.
2. For use with single function power positioning system. Function may be controlled through a switch or the VSI 6-key controller. Select single function control through joystick on page 6.
3. Needed for controlling power positioning seating systems with more than one function through controller.
4. Standard with one programmable drive profile.
5. Available with five individually programmable drive profiles.
6. Requires multi-axis joystick mount; No charge for mount.
7. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18". Suspension castors may affect seat-to-floor height up to an additional .5"-1".
8. Systems are designed for use with a 2"-3" cushion.
9. Back angle is adjustable from -85° to 105° in 5° increments.
10. These systems are not available over 20" wide by 20" deep and increase seat to floor height. Please call for actual seat to floor height. Elevating systems use 2-22NF batteries.
11. Total degrees of recline is 174° ie: if 30° degrees of precline is used, then 144° degrees of recline will be available.
12. This back height is not recommended with Recline Systems, please call for further details. Power Sliding Back (PSB) is not available on these back heights.
13. This option allows the back angle to adjust from 90° to 110° in 5° increments.
14. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
15. Not available with Recline/ESR or Recline/PSB.
16. When omitting the back there may be compatibility issues when using certain aftermarket back and/or lateral options.
17. Laterals are not available with the MaTRx PB Backs.
18. Standard front rigging will only be provided when specified. Please provide or order separately.
19. Includes a flip up foot platform. Increases seat-to-floor height.
20. Under normal usage Power Legrest actuators tend to be louder than Tilt/Recline Actuators.
21. Power Legs - Individual count as two functions.
22. Heavy duty durable construction (4.5lbs ea.). Lift and swing away setup. Please call to confirm.
23. Non-reversing relays.
25. Must order an Enhanced control box.
26. Vent trays are not standard configuration on Elevating or Precline Systems or with Single Post Backs-Please request a quotation.
27. Comes with battery tray. Uses 3-22NF Batteries
28. Articulating O² holders are not standard configuration on Precline Systems, please call for quotation.
29. Batteries may be required at time of installation, please call to confirm.

Total Retail Value \$ _____