



ALL-TERRAIN MANUAL WHEELCHAIRS

Date of Order: _____
 Dealer Name: _____
 P.O. #: _____
 Tag: _____
 Contact: _____
 Ph: _____
 Client Info: M F Weight _____ Height _____

SHIP-TO ADDRESS

<input type="checkbox"/> TerraTrek Outlook 02-TKLT <i>(Handles rugged terrain - formerly the "LT")</i> Complete with 26" all-terrain wheels with aluminum push rims, 12.5" casters, seat back, wheel locks, axle assembly, adjustable foot rests and bolt-on clothing guards. <p style="text-align: right;">\$ 3,150</p>	<input type="checkbox"/> TerraTrek Urban 03-TKCP <i>(Versatile for indoors & outdoors - formerly the "Compact")</i> Complete with 26" all-terrain wheels with aluminum push rims, 8" x 2" pneumatic casters, seat back, wheel locks, axle assembly, adjustable foot rest and bolt-on clothing guards. <p style="text-align: right;">\$ 3,150</p>
--	--

Specifications

Frame Colour

<input type="checkbox"/> Red	01-001	N/C
<input type="checkbox"/> Hunter Green	01-002	N/C
<input type="checkbox"/> Black	01-003	\$ 270
<input type="checkbox"/> Yellow*	01-004	\$ 270
<input type="checkbox"/> Midnight Blue*	01-005	\$ 270
<input type="checkbox"/> White*	01-006	\$ 270
<input type="checkbox"/> Lime Green*	01-007	\$ 270
<input type="checkbox"/> Pumpkin Orange*	01-008	\$ 270

Seat Width

Measured from outside to outside of seat tubes

<input type="checkbox"/> 12 inches*	01-009	N/C
<input type="checkbox"/> 13 inches*	01-010	N/C
<input type="checkbox"/> 14 inches*	01-011	N/C
<input type="checkbox"/> 15 inches	01-012	N/C
<input type="checkbox"/> 16 inches	01-013	N/C
<input type="checkbox"/> 17 inches	01-014	N/C
<input type="checkbox"/> 18 inches	01-015	N/C
<input type="checkbox"/> 19 <input type="checkbox"/> 20 <input type="checkbox"/> 21 inches*	01-078	\$ 190
<input type="checkbox"/> 22+ inches*	01-079	Quote

Seat Depth

Note: 18"-20" seat depths require a 1" axle extender at no charge.

<input type="checkbox"/> 13 inches*	01-016	N/C
<input type="checkbox"/> 14 inches*	01-017	N/C
<input type="checkbox"/> 15 inches	01-018	N/C
<input type="checkbox"/> 16 inches	01-019	N/C
<input type="checkbox"/> 17 inches	01-020	N/C
<input type="checkbox"/> 18 inches	01-021	\$ 130
<input type="checkbox"/> 19 <input type="checkbox"/> 20 inches	01-080	\$ 190
<input type="checkbox"/> 21+ inches*	01-081	Quote

Knee to Heel Measurement (Required)

<input type="checkbox"/> 16 <input type="checkbox"/> 17 <input type="checkbox"/> 18 <input type="checkbox"/> 19 <input type="checkbox"/> 20	01-084	STD
<input type="checkbox"/> Pediatric 8" to 15" knee to heel	01-088	\$ 150

Note: Measure from seat pan to sole of shoe.

Axle Extenders

<input type="checkbox"/> Standard	01-022	STD
<input type="checkbox"/> 1 inch	01-023	\$ 80

Rear Seat Height

Note: This is the height of the wheelchair measured from the ground to the seat rail at the back

<input type="checkbox"/> 15 inches	01-026	N/C
<input type="checkbox"/> 16 inches	01-027	N/C
<input type="checkbox"/> 17 inches	01-028	N/C
<input type="checkbox"/> 18 inches	01-029	N/C
<input type="checkbox"/> 19 inches	01-030	N/C

Back Upholstery

<input type="checkbox"/> Sling, Black Ballistic Nylon	01-031	STD
<input type="checkbox"/> Lumbar Support, Black Nylon*	01-032	\$ 60
<input type="checkbox"/> Omit Back Upholstery	01-033	(\$25)

Back Height

<input type="checkbox"/> 8" to 12" Adjustable*	01-034	N/C
<input type="checkbox"/> 13" to 17" Adjustable	01-035	N/C
<input type="checkbox"/> 18" to 20" Adjustable*	01-036	\$ 150

Push Handles

<input type="checkbox"/> No Push Handles	01-037	STD
<input type="checkbox"/> Push Handles	01-038	\$ 80
<input type="checkbox"/> Stroller Handles	01-082	\$ 200

* These items may increase delivery times from 1-3 weeks, please confirm.

On-Chair Options

This section lists the standard components and all of the options which can be ordered with the chair. Charges are as indicated and must be added to the base price.

Armrests

<input type="checkbox"/> No Arms	01-089	STD
<input type="checkbox"/> L-Style Armrests with receivers (Pair right/left)	01-039	\$ 199

Rear Wheels

<input type="checkbox"/> 26" All-Terrain Wheels	01-041	N/C
<input type="checkbox"/> Spoke Guards	01-042	\$ 100
<input type="checkbox"/> Omit Wheels	01-043	-\$ 100

Push Rims

<input type="checkbox"/> Anodized Aluminum	01-046	STD
<input type="checkbox"/> Vinyl Coated	01-047	\$ 99

Quick Release Axles

<input type="checkbox"/> Quick Release Axles	01-048	STD
<input type="checkbox"/> Quad Quick Release Axles*	01-049	\$ 75

Rear Wheel Camber

<input type="checkbox"/> 0°	01-050	N/C
<input type="checkbox"/> 3°	01-051	N/C
<input type="checkbox"/> 6°	01-052	N/C
<input type="checkbox"/> 9°	01-053	\$ 55

Wheel Locks

<input type="checkbox"/> Retractable (Pair)	01-054	STD
<input type="checkbox"/> Push to Lock (Pair)	01-055	N/C

Casters & Forks

<input type="checkbox"/> Bolt On Caster Forks	01-090	STD
<input type="checkbox"/> Quick Release Casters (Pair)	01-059	\$ 135

Anti Tippers

<input type="checkbox"/> Anti Tipper Bars with Adaptor Bushings (Pair)	01-061	\$ 90
---	--------	--------------

Cushions

<input type="checkbox"/> 2" Foam Cushion with Nylon Cover (specify size) Size: _____	01-062	\$ 95
--	--------	--------------

Seat Belts

<input type="checkbox"/> Velcro™	01-063	\$ 25
<input type="checkbox"/> Nylon Buckle	01-064	\$ 25
<input type="checkbox"/> Positioning Belt	01-065	\$ 40

Clothing Guards

<input type="checkbox"/> Bolt-on (Black ABS Plastic)	01-066	N/C
<input type="checkbox"/> Fender (Black ABS Plastic)	01-067	\$ 125

Extra Rear Wheel Camber Inserts

<input type="checkbox"/> 0° <input type="checkbox"/> 3° <input type="checkbox"/> 6°	01-068	\$ 150
<input type="checkbox"/> 9°	01-083	\$ 220

Extra Axle Extenders (Pair)

<input type="checkbox"/> 1"+ Axle Extenders	01-070	\$ 175
---	--------	---------------

Other Accessories*

<input type="checkbox"/> Fishing Rod Holder (specify left or right) L <input type="checkbox"/> R <input type="checkbox"/>	01-072	\$ 145
<input type="checkbox"/> Rifle/Shot Gun Holder (specify left or right) L <input type="checkbox"/> R <input type="checkbox"/>	01-073	\$ 145
<input type="checkbox"/> Nylon Back Pack	01-074	\$ 70
<input type="checkbox"/> Nylon Under Seat Pouch	01-075	\$ 40
<input type="checkbox"/> Crutch/Cane Holder with Adaptor Bushing (Universal right/left)	01-076	\$ 50
<input type="checkbox"/> Black Neoprene Frame Tube Protectors (Pair)	01-077	\$ 70

Two Position Legrest

(Comes with adjustable angle flip back footplate)

<input type="checkbox"/> Small - for seat widths 17"-18" (Footplate measures 10 1/2" inner width)	01-091	\$ 395
<input type="checkbox"/> Large - for seat widths 19"-20" (Footplate measures 12 1/2" inner width)	01-092	\$ 395

Footrests/ Footplates

<input type="checkbox"/> Standard (8.5" W x 12.75" D)	01-085	STD
<input type="checkbox"/> Wide (10.5"W x 11 7/8" D)	01-086	N/C
<input type="checkbox"/> Pediatric (8.5"W x 12.75"D) (knee to heel 8"-15")	01-087	N/C

Overall Dimensions

<i>Front Seat Height:</i>	21"
<i>Rear Seat Height:</i>	16" - 19"
<i>Overall Length:</i>	Standard 36" Long 36"
<i>Overall Width with 6° axle camber:</i>	
<i>Seat Width:</i>	16" 17" 18"
<i>Overall Width:</i>	28" 29" 30"
<i>Weight average based on 16" x 16" seat 30lbs.</i>	



ALL-TERRAIN MANUAL WHEELCHAIRS



700 Ensminger Road, Suite 112,
Tonawanda, New York 14150
Tel: (888) 433-6818 Fax: (888) 433-6834

* These items may increase delivery times from 1-3 weeks, please confirm.