



ALL-TERRAIN MANUAL WHEELCHAIRS

Date of Order: _____
 Dealer Name: _____
 P.O. #: _____
 Tag: _____
 Contact: _____
 Ph: _____
 Client Info: M F Weight _____ Height _____

SHIP-TO ADDRESS

<input type="checkbox"/> TerraTrek Outlook (Handles rugged terrain - formerly the "LT") Complete with 26" all-terrain wheels with aluminium push rims, 12.5" casters, seat back, wheel locks, axle assembly, adjustable foot rests and bolt-on clothing guards. \$ 4,295	02-TKLT	<input type="checkbox"/> TerraTrek Urban (Versatile for indoors & outdoors - formerly the "Compact") Complete with 26" all-terrain wheels with aluminium push rims, 8" x 2" pneumatic casters, seat back, wheel locks, axle assembly, adjustable foot rest and bolt-on clothing guards. \$ 4,295	03-TKCP
--	---------	--	---------

Specifications

Frame Colour

<input type="checkbox"/> Red	01-001	N/C
<input type="checkbox"/> Hunter Green	01-002	N/C
<input type="checkbox"/> Black	01-003	\$ 350
<input type="checkbox"/> Yellow*	01-004	\$ 350
<input type="checkbox"/> Midnight Blue*	01-005	\$ 350
<input type="checkbox"/> White*	01-006	\$ 350
<input type="checkbox"/> Lime Green*	01-007	\$ 350
<input type="checkbox"/> Pumpkin Orange*	01-008	\$ 350

Seat Width

Measured from outside to outside of seat tubes

<input type="checkbox"/> 12 inches*	01-009	N/C
<input type="checkbox"/> 13 inches*	01-010	N/C
<input type="checkbox"/> 14 inches*	01-011	N/C
<input type="checkbox"/> 15 inches	01-012	N/C
<input type="checkbox"/> 16 inches	01-013	N/C
<input type="checkbox"/> 17 inches	01-014	N/C
<input type="checkbox"/> 18 inches	01-015	N/C
<input type="checkbox"/> 19 <input type="checkbox"/> 20 <input type="checkbox"/> 21 inches*	01-078	\$ 250
<input type="checkbox"/> 22+ inches*	01-079	Quote

Seat Depth

Note: 18"-20" seat depths require a 1" axle extender at no charge.

<input type="checkbox"/> 13 inches*	01-016	N/C
<input type="checkbox"/> 14 inches*	01-017	N/C
<input type="checkbox"/> 15 inches	01-018	N/C
<input type="checkbox"/> 16 inches	01-019	N/C
<input type="checkbox"/> 17 inches	01-020	N/C
<input type="checkbox"/> 18 inches	01-021	\$ 160
<input type="checkbox"/> 19 <input type="checkbox"/> 20 inches	01-080	\$ 250
<input type="checkbox"/> 21+ inches*	01-081	Quote

Knee to Heel Measurement

(Required)

<input type="checkbox"/> 16 <input type="checkbox"/> 17 <input type="checkbox"/> 18 <input type="checkbox"/> 19 <input type="checkbox"/> 20	01-084	STD
<input type="checkbox"/> Pediatric 8" to 15" knee to heel	01-088	\$ 150

Note: Measure from seat pan to sole of shoe.

Axle Extenders

<input type="checkbox"/> Standard	01-022	STD
<input type="checkbox"/> 1 inch	01-023	\$ 80

Rear Seat Height

Note: This is the height of the wheelchair measured from the ground to the seat rail at the back

<input type="checkbox"/> 15 inches	01-026	N/C
<input type="checkbox"/> 16 inches	01-027	N/C
<input type="checkbox"/> 17 inches	01-028	N/C
<input type="checkbox"/> 18 inches	01-029	N/C
<input type="checkbox"/> 19 inches	01-030	N/C

Back Upholstery

<input type="checkbox"/> Sling, Black Ballistic Nylon	01-031	STD
<input type="checkbox"/> Lumbar Support, Black Nylon*	01-032	\$ 60
<input type="checkbox"/> Omit Back Upholstery	01-033	(\$25)

Back Height

<input type="checkbox"/> 8" to 12" Adjustable*	01-034	N/C
<input type="checkbox"/> 13" to 17" Adjustable	01-035	N/C
<input type="checkbox"/> 18" to 20" Adjustable*	01-036	\$ 150

Push Handles

<input type="checkbox"/> No Push Handles	01-037	STD
<input type="checkbox"/> Push Handles	01-038	\$ 80
<input type="checkbox"/> Stroller Handles	01-082	\$ 290

* These items may increase delivery times from 1-3 weeks, please confirm.

On-Chair Options

This section lists the standard components and all of the options which can be ordered with the chair. Charges are as indicated and must be added to the base price.

Armrests

- No Arms 01-089 STD
- L-Style Armrests with receivers (Pair right/left) 01-039 \$ 240

Rear Wheels

- 26" All-Terrain Wheels 01-041 N/C
- Spoke Guards 01-042 \$ 150
- Omit Wheels 01-043 (\$150)

Push Rims

- Anodized Aluminium 01-046 STD
- Vinyl Coated 01-047 \$ 115

Quick Release Axles

- Quick Release Axles 01-048 STD
- Quad Quick Release Axles* 01-049 \$ 100

Rear Wheel Camber

- 0° 01-050 N/C
- 3° 01-051 N/C
- 6° 01-052 N/C
- 9° 01-053 \$ 70

Wheel Locks

- Retractable (Pair) 01-054 STD
- Push to Lock (Pair) 01-055 N/C

Casters & Forks

- Bolt On Caster Forks 01-090 STD
- Quick Release Casters (Pair) 01-059 \$ 160

Anti Tipppers

- Anti Tipper Bars with Adaptor Bushings (Pair) 01-061 \$ 120

Cushions

- 2" Foam Cushion with Nylon Cover (specify size) 01-062 \$ 130
- Size: _____

Seat Belts

- Velcro™ 01-063 \$ 30
- Nylon Buckle 01-064 \$ 30
- Auto Style 01-065 \$ 45

Clothing Guards

- Bolt-on (Black ABS Plastic) 01-066 N/C
- Fender (Black ABS Plastic) 01-067 \$ 145

Extra Rear Wheel Camber Inserts

- 0° 3° 6° 01-068 \$ 180
- 9° 01-083 \$ 250

Extra Axle Extenders (Pair)

- 1"+ Axle Extenders 01-070 \$ 250

Other Accessories*

- Fishing Rod Holder (specify left or right) L R 01-072 \$ 180
- Rifle/Shot Gun Holder (specify left or right) L R 01-073 \$ 180
- Nylon Back Pack 01-074 \$ 96
- Nylon Under Seat Pouch 01-075 \$ 50
- Crutch/Cane Holder with Adaptor Bushing (Universal right/left) 01-076 \$ 90
- Black Neoprene Frame Tube Protectors (Pair) 01-077 \$ 98

Two Position Legrest

(Comes with adjustable angle flip back footplate)

- Small - for seat widths 17"-18" (Footplate measures 10 1/2" inner width) 01-091 \$ 495
- Large - for seat widths 19"-20" (Footplate measures 12 1/2" inner width) 01-092 \$ 495

Footrests/ Footplates

- Standard (8.5" W x 12.75" D) 01-085 STD
- Wide (10.5"W x 11 7/8" D) 01-086 N/C
- Pediatric (8.5" W x 12.75" D) (knee to heel 8"-15") 01-087 N/C

Overall Dimensions

- Front Seat Height:** 21"
 - Rear Seat Height:** 16" - 19"
 - Overall Length:** Standard 36"
Long 36"
 - Overall Width with 6° axle camber:**
 - Seat Width:** 16" 17" 18"
 - Overall Width:** 28" 29" 30"
- Weight average based on 16" x 16" seat 30lbs.



ALL-TERRAIN MANUAL WHEELCHAIRS



700 Ensminger Road, Suite 112,
Tonawanda, New York 14150
Tel: (888) 433-6818 Fax: (888) 433-6834
website: www.motionconcepts.com

* These items may increase delivery times from 1-3 weeks, please confirm.