


# TRI138 - INSTALLATION OVERVIEW, S626, S646, ELEVATE-ONLY, VAN SEAT



ITEM NO.	QTY.	PART NO.	DESCRIPTION	ITEM NO.	QTY.	PART NO.	DESCRIPTION
1	1	TRC0420	MAIN ASS'Y	17	8 (10)	TRH0101	1/4" FLAT WASHER
2	1	TRC0475	MANUAL ANTI-TIP ASS'Y	18	2	TRH0104	3/8" FLAT WASHER
3	1	TRC0581	PES ACTUATOR HARDWARE KIT, STD	19	10	TRH0107	5/16" FLAT WASHER
4	2 (0)	TRX0291	SPACER, ANTI-TIPPER	20	6 (10)	TRH0115	1/4" LOCK WASHER
5	1	TRX0631	ACTUATOR SUPPORT ASS'Y	21	6 (8)	TRH0201	1/4-20 NYLOC NUT
6	1	TRX0634	FRONT SUPPORT	22	2	TRH0204	3/8-16 NYLOC NUT
7	1	TRX0635	CLAMP	23	4	TRH0205	3/8-16 SHALLOW NYLOC NUT
8A	1	TRX0641	ANTI-TIP ACTUATOR SUPPORT	24	4	TRH0206	5/16-18 NYLOC NUT
8B	1	TRX2661	MOUNTING BRACKET WLD	25	6	TRH0207	5/16-18 SHALLOW NYLOC NUT
9 N/S	2	TRX0645	SPACER	26	4	TRH0303	1/4-20 x 1" HEX HEAD CAP SCREW
10	2	TRX0692	EXTENSION BLOCK	27	2	TRH0316	5/16-18 x 1" HEX HEAD CAP SCREW
11 N/S	1	TRX1852	BATTERY COVER	28A	6	TRH0335	1/4-20 x 1" SOCKET HEAD CAP SCREW GR8
12	2	TRX2001	SPACER	28B	2 (12)	TRH0354	1/4-20 x 2" SOCKET HEAD CAP SCREW GR8
13	1	TRX2489	VAN SEAT FOOTPLATE MNT WLD	29	2	TRH0378	5/16 x 2 1/4" SOCKET HEAD SHOULDER BOLT
14	2	TRM0074	SPRING	30	2	TRH0398	3/8-16 x 3" SOCKET HEAD CAP SCREW GR8
15	2	TRM0075	ANTI-TIP LOCKOUT LATCH	31	4	TRH0416	3/8-16 x 2 1/2" BUTTON SOCKET CAP SCREW
16	4	TRM0076	ID 0.32 X OD 0.38 BUSHING	32	2	TRH0493	5/16-18 x 3/4" HEX HEAD CAP SCREW