



**MATRX SEATING SERIES**  
**Effective February 2006**



**MaTRx PB Elite Back with Laterals**

Back includes outer cover and hardware to fit both 7/8" or 1" diameter back canes (3/4" inserts available upon request) and is drilled for head support attachment. Contour Depth 5". Includes formable, width adjustable lateral supports - 4" high x 8" long. Weight Capacity up to 300 lbs.

Model #	Model #	Fits Wheelchair	
Back	Extra Outer Cover*	Width	Height
HCPCS code pending	HCPCS code pending		
PBE1212-LS	PBEC1212	12"-13"	12"
PBE1312-LS	PBEC1312	12"-14"	12"
PBE1412-LS	PBEC1412	13"-15"	12"
PBE1416-LS	PBEC1416	13"-15"	16"
PBE1512-LS	PBEC1512	14"-16"	12"
PBE1516-LS	PBEC1516	14"-16"	16"
PBE1612-LS	PBEC1612	15"-17"	12"
PBE1616-LS	PBEC1616	15"-17"	16"
PBE1620-LS	PBEC1620	15"-17"	20"
PBE1712-LS	PBEC1712	16"-18"	12"
PBE1716-LS	PBEC1716	16"-18"	16"
PBE1720-LS	PBEC1720	16"-18"	20"
PBE1812-LS	PBEC1812	17"-19"	12"
PBE1816-LS	PBEC1816	17"-19"	16"
PBE1820-LS	PBEC1820	17"-19"	20"
PBE1912-LS	PBEC1912	18"-20"	12"
PBE1916-LS	PBEC1916	18"-20"	16"
PBE1920-LS	PBEC1920	18"-20"	20"
PBE2012-LS	PBEC2012	19",20"	12"
PBE2016-LS	PBEC2016	19",20"	16"
PBE2020-LS	PBEC2020	19",20"	20"

\* Includes foam