



SPECIFICATIONS: AS AN ASSEMBLY, CUSHION MUST MEET MINIMUM REQUIREMENTS OF THE CALIFORNIA TECHNICAL BULLETIN #117



ITEM	QTY	PART NO.	DESCRIPTION
1	1	TRS0060	ELITE POSTURE BACK COVER, MESHTEX, VC-1297
2	1	TRS0061	ELITE POSTURE BACK FOAM PAD
3	1	TRS0070	LUMBAR FOAM INSERT
4A	1	TRX3604	ELITE BACK PAN WLD, 16W x 16H
4B N/S	1	TRX3603	ELITE BACK PAN WLD, 16W x 12H
4C N/S	1	TRX3605	ELITE BACK PAN WLD, 16W x 20H
4D N/S	1	TRX3622	ELITE BACK PAN WLD, 17W x 12H
4E N/S	1	TRX3623	ELITE BACK PAN WLD, 17W x 16H
4F N/S	1	TRX3624	ELITE BACK PAN WLD, 17W x 20H
4G N/S	1	TRX3628	ELITE BACK PAN WLD, 18W x 12H
4H N/S	1	TRX3629	ELITE BACK PAN WLD, 18W x 16H
4J N/S	1	TRX3630	ELITE BACK PAN WLD, 18W x 20H
4K N/S	1	TRX3634	ELITE BACK PAN WLD, 19W x 12H
4L N/S	1	TRX3635	ELITE BACK PAN WLD, 19W x 16H
4M N/S	1	TRX3636	ELITE BACK PAN WLD, 19W x 20H
4N N/S	1	TRX3640	ELITE BACK PAN WLD, 20W x 12H
4P N/S	1	TRX3641	ELITE BACK PAN WLD, 20W x 16H
4Q N/S	1	TRX3642	ELITE BACK PAN WLD, 20W x 20H
4R N/S	1	TRX3646	ELITE BACK PAN WLD, 15W x 12H
4S N/S	1	TRX3647	ELITE BACK PAN WLD, 15W x 16H
4T N/S	1	TRX3648	ELITE BACK PAN WLD, 15W x 20H
4V N/S	1	TRX3652	ELITE BACK PAN WLD, 14W x 12H
4W N/S	1	TRX3653	ELITE BACK PAN WLD, 14W x 16H
4X N/S	1	TRX3654	ELITE BACK PAN WLD, 14W x 20H
4Y N/S	1	TRX3735	ELITE BACK PAN WLD, 13W x 12H
4Z N/S	1	TRX3736	ELITE BACK PAN WLD, 13W x 16H
4AA N/S	1	TRX3731	ELITE BACK PAN WLD, 12W x 12H
4BB N/S	1	TRX3732	ELITE BACK PAN WLD, 12W x 16H
5 N/S	5ft	TRM0568	RUNNER EDGE TRIM, SMALL PROFILE
6 N/S	1	TRM0577	ELITE RAISED DECAL

PART NUMBER	SIZE CONFIGURATION	LUMBAR PAD
TRS0062-1212	12 X 12	TRS0070-040
TRS0062-1216	12 X 16	TRS0070-040
TRS0062-1312	13 X 12	TRS0070-050
TRS0062-1316	13 X 16	TRS0070-050
TRS0062-1412	14 X 12	TRS0070-060
TRS0062-1416	14 X 16	TRS0070-060
TRS0062-1420	14 X 20	TRS0070-060
TRS0062-1512	15 X 12	TRS0070-070
TRS0062-1516	15 X 16	TRS0070-070
TRS0062-1520	15 X 20	TRS0070-070
TRS0062-1612	16 X 12	TRS0070-080
TRS0062-1616	16 X 16	TRS0070-080
TRS0062-1620	16 X 20	TRS0070-080
TRS0062-1712	17 X 12	TRS0070-090
TRS0062-1716	17 X 16	TRS0070-090
TRS0062-1720	17 X 20	TRS0070-090
TRS0062-1812	18 X 12	TRS0070-100
TRS0062-1816	18 X 16	TRS0070-100
TRS0062-1820	18 X 20	TRS0070-100
TRS0062-1912	19 X 12	TRS0070-110
TRS0062-1916	19 X 16	TRS0070-110
TRS0062-1920	19 X 20	TRS0070-110
TRS0062-2012	20 X 12	TRS0070-120
TRS0062-2016	20 X 16	TRS0070-120
TRS0062-2020	20 X 20	TRS0070-120

	TORONTO, ONTARIO, CANADA	SCALE: NTS
	TOLERANCE (UNLESS OTHERWISE SPECIFIED) .X ±.03 .XX ±.015 .XXX ±.005 ANGLES ±0.5°	DIMENSIONS IN: INCHES
DESCRIPTION: ELITE POSTURE BACK CUSHION ASS'Y	DRAWN BY: DC	DATE: JANUARY 25, 2006
APPROVED BY:	DWG NUMBER: TRS0062	REV 1

1	ADD 12", 13" SIZES. ADD LUMBAR FOAM	DC	02/17/06
REV	DESCRIPTION	REV. BY	DATE