





INVACARE® MATRX® ELITE FULL BACK SPECIFICATIONS



	BACK SHELL DIMENSIONS				OVERALL DIMENSIONS (CONTACT AREA)									
	A WIDTH		B HEIGHT		C INNER MAXIMUM WIDTH		D OUTER MAXIMUM WIDTH		E HEIGHT		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
MODEL	in	cm	in	cm	in	cm	in	cm	in	cm	lb	kg	lb	kg
PBEFB1620	14½	37	19½	49.5	15	38	15	38	21	53	5.4	2.4	300	136
PBEFB1720	15½	39.5	19½	49.5	16	41	16	41	21	53	5.6	2.5	300	136
PBEFB1820	16½	42	19½	49.5	17	43	17	43	21	53	5.8	2.6	300	136
PBEFB1920	17½	44.5	19½	49.5	18	46	18	46	21	53	6.0	2.7	300	136
PBEFB2020	18½	47	19½	49.5	19	48	19	48	21	53	6.2	2.8	300	136

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols™ and®.
All trademarks are owned by or licensed to Invacare
Corporation unless otherwise noted. Specifications
are subject to change without notification.
Form 12-573 Rev. 06/14 140443