





INVACARE® MATRX® GENERA BACK SPECIFICATIONS



MODEL	BACK SHELL DIMENSIONS						OVERALL DIMENSIONS (CONTACT AREA)									
	A WIDTH		B HEIGHT		C LATERAL CONTOUR HEIGHT		D INNER MINIMUM WIDTH		E OUTER MAXIMUM WIDTH		F HEIGHT		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lb	kg	lb	kg
GB1616	15	38	15	38	11	28	13	33	15½	39.5	16	41	3.8	1.7	250	136
GB1620	15	38	19	48	12½	32	13	33	15½	39.5	20	51	4.3	1.9	250	136
GB1816	17	43	15	38	11	28	15	38	17½	44.5	16	41	4.2	1.9	250	136
GB1820	17	43	19	48	12½	32	15	38	17½	44.5	20	51	4.7	2.1	250	136
GB2016	19	48	15	38	11	28	17	43	19½	49.5	16	41	4.5	2.0	250	136
GB2020	19	48	19	48	12½	32	17	43	19½	49.5	20	51	5.0	2.3	250	136

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols ™ and ©.
All trademarks are owned by or licensed to Invacare
Corporation unless otherwise noted. Specifications
are subject to change without notification.
Form 12-576 Rev. 06/14 140552