

INVACARE® MATRX®

MX1

BACK



The latest in back support technology combines aesthetic appeal with superior positioning and function. The full range of patent-pending hardware options ensures optimal adjustment for the widest range of applications, while the lightweight carbon fiber shell provides exceptional support and durability. The **Invacare® Matrx® MX1 Back** is designed for active individuals who demand peak performance and great looks.

Features and Functional Benefits

Carbon Fiber Back Shell - Durable solid carbon fiber shell provides controlled flex to reduce vibration. Contoured to provide hip/trunk stability and enhance upper extremity mobility and function. Very light weight and aesthetically appealing

Cover - Attractive combination of breathable space mesh fabric and durable reinforced vinyl panels. Cover and foam overlap the shell edge to protect the user

Adjustable Sacral/Lumbar Support - Includes lumbo-sacral support that is easily adjustable to meet individual needs

Mounting Hardware

Mono Mount (Patent Pending) - Attaches to wheelchair rigidizer bar providing complete clearance from back canes for optimal propulsion and reach. Fully angle adjustable, remains securely fastened in desired position and included track

provides height adjustment. Fits rigidizer bars 3/4" optional 7/8" and 1" band clamps available

Fixed Dual Mount (Patent Pending) - Attaches to wheelchair back canes (requires tools for removal). Fits back canes 1" in diameter - optional 3/4" and 7/8" diameter band clamps available. Attractive, low profile hardware with no sharp edges. Adjustable width to fit wheelchairs 1" wide or 1" narrower than back width, providing optimal fit and support to user*

Quick Release Dual Mount (Patent Pending) - Back cane mount hardware with quick release tabs for easy removal from wheelchair. Includes band clamps to fit 3/4", 7/8" or 1" diameter back canes. Attractive, low profile hardware and 2" adjustable width*

SuperLite - Minimalist hardware option when depth or angle adjustment are not required. Provides lightest weight option and fits back canes 1" in diameter

HCPCS Code: E2615



- 1. Mono Mount (Patent Pending)
- 2. Fixed Dual Mount (Patent Pending)
- 3. Quick Release Dual Mount (Patent Pending)
- 4. SuperLite
- 5. Adjustable Sacral/Lumbar Support
- 6. Optional Privacy Flap
- 7. 3" Contour Depth



Model #	Height	Width	Depth Adjustment	Angle Adjustment	Product Weight	Weight Capacity	Limited Warranty
IMX1-SL SuperLite	9", 12", 16"	14", 15", 16", 17", 18"	N/A	N/A	1.5 lb. and up	250 lb.	2 Years
IMX1-RM Mono Mount	9", 12", 16"	14", 15", 16", 17", 18"	N/A	30° range* +/- 15°	1.9 lb. and up	250 lb.	2 Years
IMX1-CF Fixed Dual Mount	9", 12", 16"	14", 15", 16"	2** 17", 18"	30° range*	1.7 lb. and up +/- 15°	250 lb.	2 Years
IMX1-CR Quick Release Dual Mount	9", 12", 16"	14", 15", 16", 17", 18"	2**	30° range* +/- 15°	1.7 lb. and up	250 lb.	2 Years

Products listed below are NOW Exclusively Available from Motion Concepts USA
 IMX1-SL SuperLite

*Limitation of depth/angle adjustment if used on chairs narrower than model width

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor to discuss specific situations in further detail.