

# matrx™ GENERA

MaTRx Genera back provides support and comfort at a price so affordable that no wheelchair user needs to make do with sling upholstery.



### GENERA BACK HIGHLIGHTS:

- > Attractive, water resistant and breathable cover, vinyl backed
- > Hardware fits 3/4", 7/8" and 1" back canes
- > Provides 1" depth and 20° recline adjustability
- > Quick release hardware for easy portability
- > Sturdy contoured shell, durable plated steel hardware
- > Built-in T-nuts for headrest attachment
- > Hardware adjusts to fit chairs 1" wider than back model width

### STANDARD SIZES

width	16"	18"	20"
height			
16"	●	●	●
20"	●	●	●

Weight capacity 300 lbs.

# matrx™ HIGH BACK

MaTRx High back is a stylish option with a built-in head support - great for tilt/recline applications

### HIGH BACK HIGHLIGHTS:

- > Removable, water resistant, breathable cover
- > Adjustable lumbar and neck support
- > Adjustable molded foam lateral wedges
- > Sacral/lumbar pad
- > Fits 7/8" or 1" diameter back canes

### STANDARD SIZES

width	16"	17"	18"	19"	20"	21"
height						
20"	●	●	●	●	●	●
22"	●	●	●	●	●	●
25"	●	●	●	●	●	●

Weight capacity 300 lbs.



# matrx™ CONTOUR BACK

MaTRx Contour back provides comfortable back support with a limited amount of adjustability for users with basic positioning needs.

### CONTOUR BACK HIGHLIGHTS:

- > Moisture resistant, breathable cover
- > Mildly contoured for gentle midline support
- > Ergonomically contoured aluminum back pan
- > Fixed hardware on back canes, small range of depth recline adjustment
- > Fits 7/8" or 1" diameter back canes

### STANDARD SIZES

width	15"	16"	17"	18"	19"	20"	21"	22"
height								
12"	●	●	●	●				
15"	●	●	●	●	●	●	●	●
17"	●	●	●	●	●	●	●	●
21"		●	●	●	●	●	●	●

Weight capacity 300 lbs.



### > PELVIC/THORACIC FIXED LATERAL SUPPORTS



- > Designed especially for MaTRx Contour or High Backs providing infinite/easy adjustability up and down
- > 3 pad sizes: S, M & L

### > SWINGAWAY LATERAL SUPPORTS



- > Designed especially for MaTRx Contour or High backs providing infinite/easy adjustability and obstruction-free transfers
- > 3 pad sizes: S, M & L

### > LATERAL PADS

- Small 
- Medium 
- Large 

- > Lateral Pads are available in 3 pad sizes: S, M & L
- S: 3.5" high x 5.25" long
- M: 3.75" high x 5.75" long
- L: 4.25" high x 6.5" long

### > EXTENDED SHEAR REDUCTION MOUNTING HARDWARE




- > Allows MaTRx backs to be mounted directly to Motion Concepts TRx power seating systems with recline

### > ONYX UPGRADE



- > Faux Leather automotive look
- > Available on MaTRx Contour back and High back

## PRODUCT REFERENCE



MATRIX-PB CONTOUR SHORT 12" H    MATRIX-PB CONTOUR SHORT 15" H    MATRIX-PB CONTOUR SHORT 17" H    MATRIX-PB CONTOUR SHORT 21" H

Weight capacity	300 lb	300 lb	300 lb	300 lb
Back width	15" to 18"	15" to 22"	15" to 22"	15" to 22"
Fits w/c widths	15" to 18"	15" to 22"	15" to 22"	15" to 22"
Lateral Option	✓	✓	✓	✓
Contour depth	3"	3"	3"	3"

Please see Price Lists/Order Forms for ADP/Medicare Codes

## CONTOUR BACK WEIGHTS -

Weights include mounting hardware

MATRIX-PB CONTOUR SHORT 12" H			MATRIX-PB CONTOUR SHORT 15" H			MATRIX-PB CONTOUR SHORT 17" H			MATRIX-PB CONTOUR SHORT 17" H		
SIZE	KG	LB	SIZE	KG	LB	SIZE	KG	LB	SIZE	KG	LB
1512	1.0	2.25	1515	1.2	2.75	1517	1.5	3.25	1621	3.3	7.2
1612	1.1	2.5	1615	1.4	3.0	1617	1.6	3.5	1721	3.4	7.4
1712	1.2	2.75	1715	1.5	3.25	1717	1.7	3.75	1821	3.5	7.8
1812	1.4	3.0	1815	1.6	3.5	1817	1.8	4.0	2021	3.8	8.4
1912	1.5	3.25	1915	1.7	3.75	1917	1.9	4.25			
2012	1.6	3.5	2015	1.8	4.0	2017	2.0	4.5			
2112	1.7	3.75	2115	1.9	4.25	2117	2.2	4.75			
2212	1.8	4.0	2215	2.0	4.5	2217	2.3	5.0			

## CONTOUR BACK CUSHION WEIGHTS

MATRIX-PB CONTOUR SHORT 12" H			MATRIX-PB CONTOUR SHORT 15" H			MATRIX-PB CONTOUR SHORT 17" H		
SIZE	KG	LB	SIZE	KG	LB	SIZE	KG	LB
1512	0.58	1.2	1515	0.59	1.3	1517	0.7	1.5
1612	0.6	1.3	1615	0.61	1.3	1617	0.72	1.6
1712	0.62	1.4	1715	0.62	1.4	1717	0.74	1.6
1812	0.64	1.4	1815	0.63	1.4	1817	0.76	1.7
1912	0.66	1.5	1915	0.65	1.4	1917	0.78	1.7
2012	0.68	1.5	2015	0.67	1.5	2017	0.8	1.8
2112	0.7	1.5	2115	0.71	1.6	2117	0.82	1.8
2212	0.72	1.6	2215	0.73	1.6	2217	0.84	1.9

## COVER WEIGHTS

MATRIX-PB CONTOUR SHORT 12" H			MATRIX-PB CONTOUR SHORT 15" H			MATRIX-PB CONTOUR SHORT 17" H		
SIZE	KG	LB	SIZE	KG	LB	SIZE	KG	LB
1512	0.15	0.3	1515	0.17	0.4	1517	0.19	0.42
1612	0.15	0.3	1615	0.17	0.4	1617	0.19	0.42
1712	0.15	0.3	1715	0.17	0.4	1717	0.19	0.42
1812	0.15	0.3	1815	0.17	0.4	1817	0.19	0.42
1912	0.15	0.3	1915	0.17	0.4	1917	0.19	0.42
2012	0.15	0.3	2015	0.17	0.4	2017	0.19	0.42
2112	0.15	0.3	2115	0.17	0.4	2117	0.19	0.42
2212	0.15	0.3	2215	0.17	0.4	2217	0.19	0.42