

INVACARE®

# Matrix® Elite Back

The Invacare® Matrix® Elite Back is designed to provide precise orientation within the wheelchair for optimal postural support. The attractive, lightweight design appeals to users with active lifestyles.

**HCPCS Code**  
**E2620**

## Invacare® Matrix® Elite Back

Model no. PBE

### Features and Functional Benefits

#### Outer cover:

- Oversized cover/foam for extra comfort and protection
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available

#### Hardware:

- "Easy Set" hardware installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair
- Refer to hardware page

#### Back shell:

- Stylish and breathable back shell
- Lightweight - most models under 3 lb.

#### Sacral/lumbar support:

- Includes removable positioning pad



### STANDARD

Width	12"– 20"
Height	10"
	12"– 20"
	12"– 20"
	14"– 14"
	12"– 20"
	12"– 20"
	14"– 20"
20" Full Back*	16"– 20"
Weight capacity	300 lb.

\* Optional Full back style - no scapular cut-outs

**FIXED**

2.5" of vertical adjustment on back pan  
Mount inside back shell without interfering with depth/angle adjustment  
Available in standard 5" h x 7.5" l and pediatric 3" h x 5.5" l sizes

**OFFSET FIXED**

Provides an additional 3" of vertical position (higher or lower) relative to mounting slots  
Available 5" h x 7.5" l

**SWING-AWAY**

2.5" of vertical adjustment on back pan  
Mount inside back shell without interfering w/depth/angle adjustment  
Easy push button release 4 pad sizes

**OFFSET SWING-AWAY**

Provides an additional 3" of vertical position (higher or lower)  
Mount inside back shell without interfering w/depth/angle adjustment  
Easy push button release, 4 pad sizes

Swing-Away and Offset Swing-Away sizes:  
**XS:** 3.25" h x 4.25" l, **S:** 3.5" h x 5.25" l, **M:** 3.75" h x 5.75" l, **L:** 4.25" h x 6.5" l

INVACARE®

# Matrix® Elite Deep Back

The Invacare® Matrix® Elite Deep Back is designed with lateral supports to provide optional trunk support while allowing clearance for the shoulders and hips. The attractive, lightweight design appeals to users with active lifestyle.

**HCPCS Code**  
**E2620**

**NEW** **Invacare® Matrix® Elite Deep Back**

Model no. EDB

**Features and Functional Benefits**

**Outer cover:**

- Oversized cover/foam for extra comfort and protection
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available

**Hardware:**

- "Easy Set" hardware installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair
- Refer to hardware page

**Back shell:**

- Stylish and breathable back shell
- Lightweight - most models under 4 lb.

**Sacral/lumbar support:**

- Includes removable positioning pad



**6" CONTOUR DEPTH-DEEP**



**SACRAL/  
LUMBAR  
SUPPORT**

**LATERAL PADS**



Lateral Pad sizes: (for Elite Back)

**XS:** 3.25" high x 4.25" long

**S:** 3.5" high x 5.25" long

**M:** 3.75" high x 5.75" long

**L:** 4.25" high x 6.5" long

**HEADREST ADAPTER**



Allows the headrest to be mounted directly onto Invacare Matrix Elite backs

Available in 3 sizes

**PRIVACY FLAP**



**DEEP**

Width	14" – 20"
Height	12", 16", 18", 20"
Weight capacity	300 lb.