

INVACARE®

Matrix® Flo-tech™ Cushion

The Invacare® Matrix® Flo-tech™ Cushion uses a combination of molded foam and dual layer fluid sacs to provide the highest level of skin protection for users with high risk of skin breakdown. The Fit Kit components are available to address individual postural requirements.



Invacare® Matrix® Flo-tech Cushion

Model no. FT

Features and Functional Benefits

Outer cover:

- Moisture resistant and breathable
- Includes non-slip base, hook and loop fastener strips and a lifting strap
- Expandable cover facilitates use with Fit Kit components

Platilon liner:

- An additional layer of moisture protection

Fluid sac:

- Unique dual-layer fluid in ischial/sacral area to help prevent bottoming out
- Non-temperature sensitive fluid resists freezing
- No kneading or maintenance required

Foam base:

- High quality, comfortable, resilient foam base
- Molded in hook and loop fastener strips for secure attachment of Fit Kit components
- Pre-ischial contour, trochanteric support and sacral support



STANDARD

Width	14" – 20"
Depth	14" – 21"
A Height at adductor	3.375
B Height at abductor	3.375
C Height at leg trough	2.75
Weight Capacity	300 lb.

HEAVY-DUTY

Width	20" – 24"
Depth	18", 20", 22"
Weight capacity	500 lb.

PELVIC OBLIQUITY BUILD-UP

Accommodates obliquities up to 1". For correctable obliquities use under low side of pelvis and for fixed obliquities use under high side of pelvis

ADDUCTORS

FRONT - Provides increased adductor contour to thighs
REAR - Provides increased adductor contour to hips

SEAT WEDGES

Increases degree of pre-ischial contour angle to help prevent sliding

HEMI WEDGE

May be used to provide increased support to one lower extremity – useful for clients who use foot propulsion or with hip flexion limitations

ABDUCTOR

SMALL - Provides .75" increased height of abductor contour
LARGE - Provides 2" increased height of abductor contour

TAPERED WEDGES

Accommodate mild pelvic rotation, mild hip flexion discrepancies or increase adductor contour

FLUID SUPPLEMENT

May be placed to provide pelvic obliquity accommodation or increased sacral/ischial immersion

FIT KIT COMPONENTS:

