

ELITE LATERAL OPTIONS



SWING-AWAY LATERAL PAD


- Swing-Away lateral Pads are sold as an assembly with Offset Swing-Away and Swing-Away Hardware
- 4 pad sizes: XS, S, M & L

Extra Small Pad: 3.25"/8.2cm high x 4.25"/11cm long

Small Pad: 3.5"/8.9cm high x 5.25"/13.5cm long

Medium Pad: 3.75"/9.5cm high x 5.75"/14.5cm long

Large Pad: 4.25"/11cm high x 6.5"/16.5cm long



OFFSET SWING-AWAY LATERALS

- OSALEB-XS-R (XS-RIGHT)
- OSALEB-XS-L (XS-LEFT)
- OSALEB-S-R (S-RIGHT)
- OSALEB-S-L (S-LEFT)
- OSALEB-M-R (M-RIGHT)
- OSALEB-M-L (M-LEFT)
- OSALEB-L-R (L-RIGHT)
- OSALEB-L-L (L-LEFT)

- Offset Swing-Away assembly Left or Right side (each) - for user's left when mounted up or the user's right when mounted down
- Provides an additional 3"/7.2cm of vertical adjustment (higher or lower) relative to mounting slots
- Mount inside back shell without interfering with depth/angle adjustment
- Easy push button release
- 4 sizes: XS, S, M & L

SWING-AWAY LATERALS

- SALEB-XS-R (XS-RIGHT)
- SALEB-XS-L (XS-LEFT)
- SALEB-S-R (S-RIGHT)
- SALEB-S-L (S-LEFT)
- SALEB-M-R (M-RIGHT)
- SALEB-M-L (M-LEFT)
- SALEB-L-R (L-RIGHT)
- SALEB-L-L (L-LEFT)

- Swing- Away assembly Left or Right side (each)
- 2½"/6.3cm of vertical adjustment on back pan
- Mount inside back shell without interfering with depth/angle adjustment
- Easy push button release
- 4 pad sizes: XS, S, M & L



OFFSET LATERALS

- PBO-L (LEFT);
- PBOL-R (RIGHT)

- Offset Lateral Pad with Fixed Hardware(each) - for user's left when mounted up or the user's right when mounted down
- Provides an additional 3"/7.2cm of vertical adjustment (higher or lower) relative to mounting slots



FIXED LATERALS

- PBLS-L (LEFT STANDARD);
- PBLS-R (RIGHT STANDARD);
- PBSS-L (LEFT SMALL);
- PBSS-R (RIGHT SMALL)

- Small or Standard Left or Right Lateral with Fixed Hardware (each)
- 2½"/6.3cm of vertical adjustment on back pan
- Mount inside back shell without interfering with depth/angle adjustment
- Available in Standard and Pediatric pad sizes:
Small Pad Size: 5.5"/14cm long x 3"/7.6cm high
Standard Pad Size: 7.5"/19cm long x 5"/12.7cm high