

INVACARE®

Matrix® PS Cushion

The Invacare® Matrix® PS Cushion is designed to provide superior positioning, stability and comfort. The design features an anatomically contoured shape with a unique waffled ischial relief recess and an innovative reversible outer cover.

HCPCS Code
E2605/(E2606 pending)



Invacare® Matrix® PS Cushion

Model no. PS

Features and Functional Benefits

Reversible outer cover:

- Innovative reversible Dartex™ outer cover may be used smooth or fabric side out
- Dartex™ fabric is moisture resistant and breathable
- Hook and loop fastener strips

Inner cover:

- Moisture resistant
- Easily wiped clean

Cushion:

- Durable high resilient foam prevents "bottoming out"
- Anatomically molded shape for optimal positioning & comfort
- Rear radius for sacral support
- Waterfall front edge
- Made with Ultra-Fresh™ for antibacterial and odor protection

Ischial relief area:

- Waffled recess provides immersion of bony prominences
- Helps maintain pelvic position and prevent sliding

PEDIATRIC

Width	10", 12", 14"
Depth	12", 14", 16"
A Height at adductor	3.75"
B Height at abductor	3.5"
C Height at leg trough	2.5"
Weight capacity	150 lb.

STANDARD

Width	16", 18", 20"
Depth	16", 18", 20"
A Height at adductor	4"
B Height at abductor	3.5"
C Height at leg trough	2.5"
Weight capacity	300 lb.

HEAVY-DUTY

Width	21" – 30"
Depth	21" – 30"
A Height at adductor	4"
B Height at abductor	3.5"
C Height at leg trough	2.5"
Weight capacity	600 lb.

CUSTOM MODIFICATIONS



Leg length discrepancies, width/depth reductions
Foam, inner and outer covers modified to your specifications

CUSHION RIGIDIZER



Lightweight, slim profile solid seat insert
Contoured to eliminate "hammock effect" of sling upholstery

PELVIC OBLIQUITY KIT



Easily installed on Invacare Matrix PS cushions
Accommodates .5" or 1" obliquities