

INVACARE®

# Matrix® PB Back

The Invacare® Matrix® PB Posture Back is contoured to provide midline support and positioning. The extensive range of sizes provides exceptional support and comfort for a variety of customers.

**HCPCS Code**  
**E2615/E2616**



**BACK SHELL**

**5" CONTOUR DEPTH-STANDARD**



**SACRAL/  
LUMBAR  
SUPPORT**

## Invacare® Matrix® PB Back

Model no. MPB

### Features and Functional Benefits

#### Foam:

- Molded foam wraps around shell for maximum comfort and protection

#### Outer cover:

- Moisture resistant and breathable
- Available in black or Kid\*ab\*ra print
- Optional privacy flap is available

#### Sacral/lumbar support:

- Includes removable positioning pad

#### Back shell:

- Contoured shell provides exceptional midline support and positioning
- Durable, lightweight ABS back shell includes built-in plate for headrest attachment

#### Hardware:

- "Easy Set" hardware is easily installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair
- Refer to Easy Set Hardware page

### STANDARD

Width 14", 16", 17", 18", 20"

Height

12" 14"–20"

16" 14"–20"

20" 16"–20"

Weight capacity 300 lb.



Kid\*ab\*ra fabric optional

### PRIVACY FLAP



### FOAM-IN-PLACE KIT



Individually molded insert produced in minutes with user seated in back  
Optional high stretch cover available

### VISCO FOAM UPGRADE



Laminate of .5" high resilient foam and .625" Visco Foam  
Available on Invacare Matrix PB and Elite Backs  
Molded foam is not available on Invacare Matrix PB Back with this option

### HEAVY-DUTY

Width 20", 22", 24"

Height

16" 20", 22", 24"

20" 20", 22", 24"

Weight capacity 500 lb.

### EXTRA WIDE

Width 24"–30"

Height

16" 24"–30"

20" 24"–30"

Weight capacity 500 lb.

INVACARE®


# Matrix® PB Deep Back

The Invacare® Matrix® PB Deep back is contoured to provide increased midline support and positioning. The extensive range of sizes provides exceptional comfort and support for a wide range of user needs.

**HCPCS Code**  
E2615/E2616



**BACK SHELL**



**7" CONTOUR DEPTH-DEEP**



**SACRAL/  
LUMBAR  
SUPPORT**

## Invacare® Matrix® PB Deep Back

Model no. PBD (Deep)

### Features and Functional Benefits

#### Foam:

- Molded foam wraps around shell for maximum comfort and protection

#### Outer cover:

- Moisture resistant and breathable
- Available in black or Kid\*ab\*ra print
- Optional privacy flap is available

#### Sacral/lumbar support:

- Includes removable positioning pad

#### Back shell:

- Contoured shell (standard/deep) provides exceptional midline support and positioning
- Durable, lightweight ABS back shell includes built-in plate for headrest attachment

#### Hardware:

- "Easy Set" hardware is easily installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair
- Refer to Easy Set Hardware page

### DEEP

Width	14", 16", 17", 18", 20"
Height	
12"	14" – 20"
16"	14" – 20"
20"	16" – 20"

Weight capacity 300 lb.



Kid\*ab\*ra fabric optional

### PRIVACY FLAP



### FOAM-IN-PLACE KIT



Individually molded insert produced in minutes with user seated in back  
Optional high stretch cover available

### VISCO FOAM UPGRADE



Laminate of .5" high resilient foam and .625" Visco Foam  
Available on Invacare Matrix PB and Elite Backs  
Molded foam is not available on Invacare Matrix PB Back with this option

### HEAVY-DUTY

Width	20", 22", 24"
Height	
16"	20", 22", 24"
20"	20", 22", 24"

Weight capacity 500 lb.

### EXTRA WIDE

Width	24" – 30"
Height	
16"	24" – 30"
20"	24" – 30"

Weight capacity 500 lb.