

INVACARE®

# Matrix® PS Cushion

The Invacare® Matrix® PS Cushion is designed to provide superior positioning, stability and comfort. The design features an anatomically contoured shape with a unique waffled ischial relief recess and an innovative reversible outer cover.

HCPCS Code  
E2605/E2606



## Invacare® Matrix® PS Cushion

Model no. PS

### Features and Functional Benefits

#### Reversible outer cover:

- Innovative reversible Dartex™ outer cover may be used smooth or fabric side out
- Dartex™ fabric is moisture resistant and breathable
- Hook and loop fastener strips

#### Inner cover:

- Moisture resistant
- Easily wiped clean
- Zippered cover allows access for cushion modifications

#### Cushion:

- Durable high resilient foam prevents "bottoming out"
- Anatomically molded shape for optimal positioning & comfort
- Rear radius for sacral support
- Waterfall front edge
- Made with Ultra-Fresh™ for antibacterial and odor protection

#### Ischial relief area:

- Waffled recess provides immersion of bony prominences
- Helps maintain pelvic position and prevent sliding

### PEDIATRIC

Width	10", 12", 14"
Depth	10", 12", 14", 16"
A Height at adductor	3.75"
B Height at abductor	3.5"
C Height at leg trough	2.5"
Weight capacity	150 lb.

### STANDARD

Width	16", 18", 20"
Depth	16", 18", 20"
A Height at adductor	4"
B Height at abductor	3.5"
C Height at leg trough	2.5"
Weight capacity	300 lb.

### HEAVY-DUTY

Width	21" – 30"
Depth	21" – 30"
A Height at adductor	4"
B Height at abductor	3.5"
C Height at leg trough	2.5"
Weight capacity	600 lb.

### CUSTOM MODIFICATIONS



Leg length discrepancies, width/depth reductions  
Foam, inner and outer covers modified to your specifications

### CUSHION RIGIDIZER



Lightweight, slim profile solid seat insert  
Contoured to eliminate "hammock effect" of sling upholstery

### PELVIC OBLIQUITY KIT



Easily installed on Invacare Matrix PS cushions  
Accommodates .5" or 1" obliquities