


## POSITIONING STRAPS

### 2 Point Padded Lap Belt:



- ▶ Center-pull push button, padded lap belt
- ▶ Two points of attachment to the wheelchair
- ▶ Suitable when the belt is supporting users with low tone or weakness
- ▶ HCPCS code EO978




- ▶ A=2" (50mm), B=3" (76mm), C=11" (28cm), D=60" (152cm)

### 4 Point Padded Lap Belt:



- ▶ Center-pull push button, padded lap belt
- ▶ Four points of attachment to the wheelchair
- ▶ Primary attachments position the belt and secondary attachments anchor it in position
- ▶ Holds pelvis in place more aggressively for situations when the person has a lot of movement
- ▶ HCPCS code EO978



- ▶ A=2" (50mm), B=3" (76mm), C=11" (28cm), D=60" (152cm)

### Air Mesh Calf Strap:



- ▶ Available in: 3½" (9cm) High
  - Medium (Fits chairs 15" -18" wide)
  - Large (Fits chairs 18" -22" wide)
- ▶ HCPCS code EO995