

INVACARE®

# Matrx® Vi Cushion

The Invacare® Matrx® Vi Cushion is designed to provide skin protection and positioning for users at moderate to high risk of skin breakdown. The unique three part foam construction targets skin protection to where it is needed most, and features improved stability, support and comfort.

**HCPCS Code**  
**E2607/E2608**



## Invacare® Matrx® Vi Cushion

Model no. MA(VVDD)-Vi

### Features and Functional Benefits

#### Outer cover:

- Moisture resistant and breathable
- Includes slip resistant base, hook and loop fasteners, lifting strap
- Available in black and Kid\*ab\*ra print

#### Inner cover:

- Moisture resistant
- Easily wipes clean
- Zippered cover allows access for cushion modifications

#### Top layer:

- Soft (ILD) High Resilient (HR) foam for optimal comfort and support
- Pre-ischial contour provides stability to pelvis
- Abductor/adductor contour promotes midline leg positioning
- Made with Ultra-Fresh™ for antibacterial and odor protection

#### Visco foam insert:

- Molded visco-elastic foam ischial relief insert
- Anatomically targeted to ischial/sacral area
- Provides immersion to sensitive bony prominences

#### Bottom layer:

- Firm (ILD) High Resilient (HR) foam for support and durability
- Trochanteric shelf provides stability to hips
- Waterfall front edge
- Rear radius for sacral support



### STANDARD

Width	14", 16", 18", 20", 22"
Depth	15" – 20"
A. Height at adductor	3.75"
B. Height at abductor	3.25"
C. Height at leg trough	2.5"
Weight capacity	300 lb.

### HEAVY-DUTY

Width	18" – 30"
Depth	18" – 26"
A. Height at adductor	3.75"
B. Height at abductor	3.25"
C. Height at leg trough	2.5"
Weight capacity	600 lb.



Kid\*ab\*ra fabric standard on pediatric sizes

### PEDIATRIC

Width	12", 14"
Depth	12", 14"
A. Height at adductor	2.875"
B. Height at abductor	2.75"
C. Height at leg trough	2"

Kid\*ab\*ra fabric standard  
black fabric optional

### DROP BASE



Includes built-in drop base, drop hooks, DuraClip fasteners

Allows adjustment of seat angle

### PELVIC OBLIQUITY KIT



Easily installed on all Invacare Matrx Vi/ Kid\*ab\*ra-Vi cushions/bases

Accommodates .5" or 1" obliquities

### CUSTOM MODIFICATIONS



Leg length discrepancies, width/depth reductions

Foam, inner and outer covers modified to your specifications