


INVACARE®

# Matrix® Elite Deep Back



MODEL NO.	A		B		C		D		E		F		TRUNK WIDTH WITH CUSHION		CUSHION DEPTH		OVERALL WIDTH WITH CUSHION		OVERALL HEIGHT WITH CUSHION		OVERALL DEPTH WITH CUSHION		BACK WEIGHT (W/O HDWR)	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs	kgs
<b>EDB1412</b>	13	33	12	31	6½	16½	9¼	23½	8	20	2½	6½	11½	29	6	15½	14	36	14	36	8	20	3.2	1.4
<b>EDB1416</b>	13	33	16	41	6½	16½	9¼	23½	8½	21½	3	7½	11½	29	6	15½	14	36	17	43	8	20	3.4	1.5
<b>EDB1418</b>	13	33	18	46	6½	16½	9¼	23½	8½	21½	3	7½	11½	29	6	15½	14	36	19	48	8	20	3.6	1.6
<b>EDB1420</b>	13	33	20	51	6½	16½	9¼	23½	8½	21½	3	7½	11½	29	6	15½	14	36	21	53	8	20	3.8	1.7
<b>EDB1512</b>	14	36	12	31	6½	16½	10¼	26½	8	20	2½	6½	12½	32	6	15½	15	38	14	36	8	20	3.3	1.5
<b>EDB1516</b>	14	36	16	41	6½	16½	10¼	26½	8½	21½	3	7½	12½	32	6	15½	15	38	17	43	8	20	3.5	1.6
<b>EDB1518</b>	14	36	18	46	6½	16½	10¼	26½	8½	21½	3	7½	12½	32	6	15½	15	38	19	48	8	20	3.7	1.7
<b>EDB1520</b>	14	36	20	51	6½	16½	10¼	26½	8½	21½	3	7½	12½	32	6	15½	15	38	21	53	8	20	3.9	1.8
<b>EDB1612</b>	15	38	12	31	6½	16½	11¼	28½	8	20	2½	6½	13½	34½	6	15½	16	41	14	36	8	20	3.4	1.5
<b>EDB1616</b>	15	38	16	41	6½	16½	11¼	28½	8½	21½	3	7½	13½	34½	6	15½	16	41	17	43	8	20	3.6	1.6
<b>EDB1618</b>	15	38	18	46	6½	16½	11¼	28½	8½	21½	3	7½	13½	34½	6	15½	16	41	19	48	8	20	3.8	1.7
<b>EDB1620</b>	15	38	20	51	6½	16½	11¼	28½	8½	21½	3	7½	13½	34½	6	15½	16	41	21	53	8	20	4.0	1.8

INVACARE®

# Matrix® Elite Deep

Back

MODEL NO.	A				B				C				D				E				F				TRUNK WIDTH		CUSHION		OVERALL WIDTH		OVERALL HEIGHT		OVERALL DEPTH		BACK WEIGHT	
	WITH CUSHION		DEPTH		WITH CUSHION		WITH CUSHION		WITH CUSHION		WITH CUSHION		WITH CUSHION		WITH CUSHION		WITH CUSHION		WITH CUSHION		WITH CUSHION		(W/O HDWR)													
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs	kgs										
<b>EDB1712</b>	16	41	12	31	6½	16½	12¼	31½	8	20	2½	6½	14½	37	6	15½	17	43	14	36	8	20	3.5	1.6												
<b>EDB1716</b>	16	41	16	41	6½	16½	12¼	31½	8½	21½	3	7½	14½	37	6	15½	17	43	17	43	8	20	3.7	1.7												
<b>EDB1718</b>	16	41	18	46	6½	16½	12¼	31½	8½	21½	3	7½	14½	37	6	15½	17	43	19	48	8	20	3.9	1.8												
<b>EDB1720</b>	16	41	20	51	6½	16½	12¼	31½	8½	21½	3	7½	14½	37	6	15½	17	43	21	53	8	20	4.1	1.9												
<b>EDB1812</b>	17	43	12	31	6½	16½	13¼	33½	8	20	2½	6½	15½	39½	6	15½	18	46	14	36	8	20	3.6	1.6												
<b>EDB1816</b>	17	43	16	41	6½	16½	13¼	33½	8½	21½	3	7½	15½	39½	6	15½	18	46	17	43	8	20	3.8	1.7												
<b>EDB1818</b>	17	43	18	46	6½	16½	13¼	33½	8½	21½	3	7½	15½	39½	6	15½	18	46	19	48	8	20	4.1	1.9												
<b>EDB1820</b>	17	43	20	51	6½	16½	13¼	33½	8½	21½	3	7½	15½	39½	6	15½	18	46	21	53	8	20	4.3	2.0												
<b>EDB1912</b>	18	46	12	31	6½	16½	14¼	36½	8	20	2½	6½	16½	42	6	15½	19	48	14	36	8	20	3.7	1.7												
<b>EDB1916</b>	18	46	16	41	6½	16½	14¼	36½	8½	21½	3	7½	16½	42	6	15½	19	48	17	43	8	20	4.0	1.8												
<b>EDB1918</b>	18	46	18	46	6½	16½	14¼	36½	8½	21½	3	7½	16½	42	6	15½	19	48	19	48	8	20	4.2	1.9												
<b>EDB1920</b>	18	46	20	51	6½	16½	14¼	36½	8½	21½	3	7½	16½	42	6	15½	19	48	21	53	8	20	4.4	2.0												
<b>EDB2012</b>	19	48	12	31	6½	16½	15¼	38½	8	20	2½	6½	17½	44½	6	15½	20	51	14	36	8	20	3.9	1.8												
<b>EDB2016</b>	19	48	16	41	6½	16½	15¼	38½	8½	21½	3	7½	17½	44½	6	15½	20	51	17	43	8	20	4.1	1.9												
<b>EDB2018</b>	19	48	18	46	6½	16½	15¼	38½	8½	21½	3	7½	17½	44½	6	15½	20	51	19	48	8	20	4.3	2.0												
<b>EDB2020</b>	19	48	20	51	6½	16½	15¼	38½	8½	21½	3	7½	17½	44½	6	15½	20	51	21	53	8	20	4.5	2.1												