




INVACARE®

# Matrix® Elite Full Back



MODEL NO	A		B		C		D		E		F		TRUNK WIDTH WITH CUSHION		CUSHION DEPTH		OVERALL WIDTH WITH CUSHION		OVERALL HEIGHT WITH CUSHION		OVERALL DEPTH WITH CUSHION		BACK WEIGHT (W/O HDWR)		
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs	kgs	
													in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs
PBEFB1620	14½	37	19½	49½	2¾	7	11¼	28½	N/A				15	38	2	5	15	38	21	53	4	10	5.4	2.4	
PBEFB1720	15½	39½	19½	49½	2¾	7	12¼	31					16	41	2	5	16	41	21	53	4	10	5.6	2.5	
PBEFB1820	16½	42	19½	49½	2¾	7	13¼	33½					17	43	2	5	17	43	21	53	4	10	5.8	2.6	
PBEFB1920	17½	44½	19½	49½	2¾	7	14¼	36					18	46	2	5	18	46	21	53	4	10	6.0	2.7	
PBEFB2020	18½	47	19½	49½	2¾	7	15¼	38½					19	48	2	5	19	48	21	53	4	10	6.2	2.8	