



INVACARE®
Matrix® Elite
 Back



MODEL NO.	A		B		C		D		E		F		TRUNK WIDTH WITH CUSHION		CUSHION DEPTH		OVERALL WIDTH WITH CUSHION		OVERALL HEIGHT WITH CUSHION		OVERALL DEPTH WITH CUSHION		BACK WEIGHT (W/O HDWR)	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs	kgs
	PBE1212	9¾	25	12	31	3	7½	6¼	16	5¾	14½	2	5	8½	21½	3	7½	11	28	14	36	5	12½	2.0
PBE1214	9¾	25	13½	34½	3	7½	6¼	16	8	20	1¼	3	8½	21½	3	7½	11	28	15	38	5	12½	2.2	1.0
PBE1216	9¾	25	16	41	3	7½	6¼	16	6½	16½	5¾	14½	8½	21½	3	7½	11	28	17	43	5	12½	2.3	1.0
PBE1218	9¾	25	17	43	3	7½	6¼	16	8¼	21	1½	4	8½	21½	3	7½	11	28	18	46	5	12½	2.5	1.1
PBE1220	9¾	25	19	48	3	7½	6¼	16	8¼	21	1½	4	8½	21½	3	7½	11	28	20	51	5	12½	2.7	1.2
PBE1312	10¼	26½	11¾	30	3½	9	8¼	21	7	18	3	7½	9½	24	3	7½	12	31	14	36	5	12½	2.2	1.0
PBE1314	10¼	26½	13½	34½	3½	9	8¼	21	7	18	3	7½	9½	24	3	7½	12	31	15	38	5	12½	2.3	1.0
PBE1316	10¼	26½	16	41	3½	9	8¼	21	8	20	3¼	8½	9½	24	3	7½	12	31	17	43	5	12½	2.5	1.1
PBE1318	10¼	26½	17	43	3	7½	8¼	21	8	20	1½	4	9½	24	3	7½	12	31	18	46	5	12½	2.7	1.2
PBE1320	10¼	26½	19	48	3	7½	8¼	21	8	20	1½	4	9½	24	3	7½	12	31	20	51	5	12½	2.9	1.3
PBE1412	12½	32	12	31	3½	9	9¼	23½	5¼	13½	4½	11½	11½	29	3	7½	14	36	13	33	5	12½	2.4	1.1
PBE1414	12½	32	13½	34½	3½	9	9¼	23½	7½	19	2¼	6	11½	29	3	7½	14	36	15	38	5	12½	2.6	1.2
PBE1416	12½	32	16	41	3½	9	9¼	23½	9¼	23½	2¼	7	11½	29	3	7½	14	36	17	43	5	12½	2.8	1.3
PBE1418	12½	32	18	46	3½	9	9¼	23½	6	15½	5¼	13½	11½	29	3	7½	14	36	19	48	5	12½	3.0	1.4
PBE1420	12½	32	20	51	3½	9	9¼	23½	7¼	18½	4¾	12	11½	29	3	7½	14	36	21	53	5	12½	3.2	1.5

INVACARE®
Matrix® Elite
 Back

MODEL NO.	A		B		C		D		E		F		TRUNK WIDTH		CUSHION		OVERALL WIDTH		OVERALL HEIGHT		OVERALL DEPTH		BACK WEIGHT	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	WITH CUSHION		DEPTH		WITH CUSHION		WITH CUSHION		WITH CUSHION		(W/O HDWR)	
													in	cm	in	cm	in	cm	in	cm	lbs	kgs		
PBE1512	13½	34½	12	31	3½	9	10¼	26½	6	15½	3½	9	12½	32	3	7½	15	38	13	33	5	12½	2.6	1.2
PBE1516	13½	34½	16	41	3½	9	10¼	26½	6	15½	5¾	14½	12½	32	3	7½	15	38	17	43	5	12½	3.0	1.4
PBE1518	13½	34½	18	46	3½	9	10¼	26½	6	15½	5¼	13½	12½	32	3	7½	15	38	19	48	5	12½	3.2	1.5
PBE1520	13½	34½	20	51	3½	9	10¼	26½	7¼	18½	4¾	12	12½	32	3	7½	15	38	21	53	5	12½	3.4	1.5
PBE1612	14½	37	12	31	3½	9	11¼	28½	6	15½	3½	9	13½	34½	3	7½	16	41	13	33	5	12½	2.8	1.3
PBE1616	14½	37	16	41	3½	9	11¼	28½	6	15½	5¾	14½	13½	34½	3	7½	16	41	17	43	5	12½	3.2	1.5
PBE1618	14½	37	18	46	3½	9	11¼	28½	6	15½	5¼	13½	13½	34½	3	7½	16	41	19	48	5	12½	3.4	1.5
PBE1620	14½	37	20	51	3½	9	11¼	28½	7¼	18½	4¾	12	13½	34½	3	7½	16	41	21	53	5	12½	3.6	1.6
PBE1712	15½	39½	12	31	3½	9	12¼	31½	6	15½	3½	9	14½	37	3	7½	17	43	13	33	5	12½	3.0	1.4
PBE1716	15½	39½	16	41	3½	9	12¼	31½	6	15½	5¾	14½	14½	37	3	7½	17	43	17	43	5	12½	3.3	1.5
PBE1718	15½	39½	18	46	3½	9	12¼	31½	6	15½	5¼	13½	14½	37	3	7½	17	43	19	48	5	12½	3.6	1.6
PBE1720	15½	39½	20	51	3½	9	12¼	31½	7¼	18½	4¾	12	14½	37	3	7½	17	43	21	53	5	12½	3.8	1.7
PBE1812	16½	42	12	31	3½	9	13¼	33½	6	15½	3½	9	15½	39½	3	7½	18	46	13	33	5	12½	3.2	1.5
PBE1816	16½	42	16	41	3½	9	13¼	33½	6	15½	5¾	14½	15½	39½	3	7½	18	46	17	43	5	12½	3.4	1.5
PBE1818	16½	42	18	46	3½	9	13¼	33½	6¾	16	4¾	12	15½	39½	3	7½	18	46	19	48	5	12½	3.7	1.7
PBE1820	16½	42	20	51	3½	9	13¼	33½	7¼	18½	4¾	12	15½	39½	3	7½	18	46	21	53	5	12½	4.0	1.8
PBE1912	17½	44½	12	31	3½	9	14¼	36½	6	15½	3½	9	16½	42	3	7½	19	48	13	33	5	12½	3.3	1.5
PBE1916	17½	44½	16	41	3½	9	14¼	36½	6	15½	5¾	13½	16½	42	3	7½	19	48	17	43	5	12½	3.6	1.6
PBE1918	17½	44½	18	46	3½	9	14¼	36½	7¼	18½	4¾	12	16½	42	3	7½	19	48	19	48	5	12½	3.9	1.8
PBE1920	17½	44½	20	51	3½	9	14¼	36½	7¼	18½	4¾	12	16½	42	3	7½	19	48	21	53	5	12½	4.2	1.9
PBE2012	18½	47	12	31	3½	9	15¼	38½	6	15½	3½	9	17½	44½	3	7½	20	51	13	33	5	12½	3.5	1.6
PBE2016	18½	47	16	41	3½	9	15¼	38½	6	15½	5¾	14½	17½	44½	3	7½	20	51	17	43	5	12½	3.8	1.7
PBE2018	18½	47	18	46	3½	9	15¼	38½	7¼	18½	4¾	12	17½	44½	3	7½	20	51	19	48	5	12½	4.1	1.9
PBE2020	18½	47	20	51	3½	9	15¼	38½	7¼	18½	4¾	12	17½	44½	3	7½	20	51	21	53	5	12½	4.4	2.0